



Larynx Needle

Trans-oral needle when treating of vocal fold injections with steroids, gel or local anesthesia in larynx. Ideal for example, treating vocal cord paresis. Single use needle with a standard luer lock hub that fits all luer syringes.

Developed together with Swedish phoneticians.

Art.No.	Description	Ward/Trsp.box
606261230	0,5 x 230 mm, 25G	10/150

Ward box = minimum order qty

CE-marked Class IIa
Material:
Hub: PP, Needle: Medical grade stainless steel
Manufacturer:
Medioplast AB / www.medioplast.com