



**Indication:**  
 Can be used as a fixation dressing, under plaster and for dermatological use. Particularly suitable over joints, amputations and severe bandaged areas. Only to be used on undamaged skin.

## Latex Free Tubular bandage, double stretch

MEDIPLAST® offers DANAFAST Double Stretch, a tubular elastic bandage that stretches in length and width and thus easy to apply. Can be cut in desired length and can be applied in single or double layer, no additional fixation is required. Good air permeability. Colour coded with a woven line which makes choice of size easy.

**DANAFAST Double Stretch:**  
 Viscose, Polyamide, Elastane  
 Length: 10 meters unstretched

Art.No	Description	Ward box	Trsp.box
6050761	RED = >3,5 cm, small arms and legs	1	35
6050762	GREEN = >5,5 cm, feet, hands, normal arms, small legs	1	28
6050763	BLUE = >8 cm, normal legs, large arms, thighs	1	20
6050764	YELLOW = >12 cm, small trunks, large thighs	1	14
6050765	BEIGE = >18 cm, adult trunks	1	7
6050766	PURPLE = >20 cm, large adult trunk	1	7

Ward box = min order qty

**DANAFast Double Stretch:**

Viscose, Polyamide, Elastane

Length: 1 meter unstretched

## Single packed

Art.No	Description	Ward box	Trsp.box
6050771	RED = >3,5 cm, small arms and legs	72	576
6050772	GREEN = >5,5 cm, feet, hands, normal arms, small legs	72	576
6050773	BLUE = >8 cm, normal legs, large arms, thighs	72	432
6050774	YELLOW = >12 cm, small trunks, large thighs	36	288
6050775	BEIGE = >18 cm, adult trunks	36	144
6050776	PURPLE = >20 cm, large adult trunk	36	144

Ward box = min order qty