



Best Quality made by a German Manufacturer

danumed® FeedBag

Empty bag for enteral nutrition and fluid balancing

The danumed® Feedbag is an empty feeding bag with a large filling volume of 1500 ml. The ENPlus port on the outlet tube allows connection to the danumed® Gravity feeding sets and to all other commercially available gravity sets for enteral feeding.

SAFETY & COMFORT

- **High quality material**
- **ENPlus port** for connection to enteral feeding systems
- **Large bag volume** of 1500 ml
- **Inscribable bag material** for clear patient identification
- **Scaling** in 100ml-steps
- With hose clamp for **quick application stop**
- Large topfiller for **convenient filling** of the bag



danumed® FeedBag

The high-quality, economical version made of latex-free, DEHP- and BPA-free bag material.

REF	Product description	Pieces/sales unit
DA 1051 1000	danumed® FeedBag	30 pieces

The danumed® FeedBag is individually sterile packed. 30 pieces in a shipping box.

The danumed® FeedBags are
• **Latex-free** • **DEHP-free** • **BPA-free**

Forhandler:
MEDI PLAST
Mediplast AS
Bjørnstjerne Bjørnsonsgate 110
3044 Drammen
T 32 88 11 00
mediplast@mediplast.no
www.mediplast.com

danumed® FeedBag Premium

The first-class version made of environmentally friendly EVA material and with premium topfiller.

REF	Product description	Pieces/sales unit
DA 1051 1001	danumed® FeedBag Premium	30 pieces

The danumed® FeedBag Premium is individually sterile packed. 30 pieces in a shipping box.

The danumed® FeedBags Premium are
• **Latex-free** • **PVC-free** • **DEHP-free** • **BPA-free**