



Flushing sets

FLUSHING SET MEDIPLAST

For flushing drip at various scopic examinations. With spike, soft drip chamber and roller clamp, universal connector in silicone for connection to the instrument. Single use, sterile packed in inner PE pouch + outer peel pack

| Art.no. | Description | Ward/trsp.box |
|----------|--|---------------|
| 60Y-280T | 1-way, phthalate free PVC tube (OD 6,8 x ID 4,8 mm) Tube length: 280 - Total length: 310 cm | 10/20 |
| 60Y-190T | Tube length: 190 - Total length: 210 cm | 10/20 |

UROLOGICAL FLUSHING SET, 2-way for pre and postoperative use,

Y-adapter with two spikes to flushing bags, roller clamp, and large chamber which allows a constant flow and minimum risk for air bubbles.

Integrated nipple for connection to the catheter. Universal connector in silicone or enclosed Luer Connector, for connecting to the cystoscope.

Single use, sterile packed in inner PE pouch + outer peel pack

| Art.no. | Description | Ward/trsp.box |
|----------|--|---------------|
| 60Y-147T | 2-way Y-connector, 2 spikes for flushing drip bags, phthalate free PVC tube (OD 6,8 x ID 4,8 mm) Tube Tube length: 150 - Total length: 210 cm | 10/20 |
| 60Y-200T | Tube length: 200 - Total length: 260 cm | 10/20 |

Material: phthalate-free PVC, softener TOTM / ABS / Silicone / PP / PE
Trsp.box = min order qty

