

Interventional
imaging

▶ SEAL ONE®

All-In-One Radial Compression Device

SEAL ONE



- Precise adjustment and control of compression level
- Targeting the radial artery
- Enhanced visibility of the puncture site
- Secured decompression

▶ SEAL ONE®



○ FEATURES

▶ Precise adjustment and control of compression level⁽¹⁾

Compression knob with scale

- Limited risk of radial artery occlusion
- Optimization of use time

▶ Targeting the radial artery

Specific design of the compression pad

- No compression of the ulnar artery
- Ensured venous return

▶ Enhanced visibility of the puncture site

Transparent device

- Better visual control

▶ Two-step decompression

Secured button

- Limited risk of unintended decompression

▶ Accessory free

Compression / Decompression knob

- Simplified handling

▶ Display of positioning time

Time indicator

- Easier monitoring of decompression protocol

Reference	Description	Packaging
0259NA	Radial compression device	30 unit-box

(1) Each half degree of graduation increases the compression by 0.35mm and reciprocally for the decompression.

This medical device is CE marked by LNE-GMED (CE 0459).

PEROUSE
MEDICAL

135 Route Neuve – 69540 Irigny – France
Tél: 33 (0)4 72 39 74 14 - Fax: 33 (0)4 78 51 89 67

www.perousemedical.com

Återförsäljare:
MEDI PLAST

Medioplast AB
Box 1004, 212 10 Malmö
T 020-78 80 35
mediplast.info@mediplast.com
www.mediplast.com

Forhandler:
MEDI PLAST

Medioplast AS
Tollbugata 115, 3041 Drammen
T 32 88 11 00
mediplast@mediplast.no
www.mediplast.com

