



Indication: For fixation around joints and trunks which are difficult to bandage. Can also be used for fixation of primary dressings and tubings. Only to be used on undamaged skin.

Fixation bandage

MEDIPLAST offers DANAMULL fixation bandages. The bandages are blood and salves repellent and have good aeration.

Double weaving with woven edges.
 Elasticity: 100%
 Viscose, Polyamide
 Length: stretched

| Art.No | Description | Ward box | Trsp.box |
|---------|------------------|----------|----------|
| 6050601 | 4 cm x 4 meters | 20 | 1440 |
| 6050602 | 6 cm x 4 meters | 20 | 960 |
| 6050603 | 8 cm x 4 meters | 20 | 720 |
| 6050604 | 10 cm x 4 meters | 20 | 600 |
| 6050605 | 12 cm x 4 meters | 20 | 480 |

Manufactured in Denmark

Ward box = min order qty