

Spiral Line



- Intended for transport and transfer of liquid drugs and drug solutions in connection with pressure or gravity infusion
- The flexible arrangement ensures an increased mobility for both patient and staff and resumes its original shape after extraction
- The tube is tested for 4 bar
- Compatible with other standardized luer-lock systems and can easily be connected to subsequent systems
- For more detailed information, read the instructions before use

1 x 2 mm (ID x OD)

| Art. no. | Description | | Ward/Trsp.box |
|----------|--------------------------------------|----------|---------------|
| 60A1226 | Priming volume 1,2 ml, length 150 cm | PVC-free | 50/100* |
| 60A1229 | Priming volume 2,4 ml, length 300 cm | PVC-free | 25/50* |
| 60A1230 | Priming volume 3,2 ml, length 400 cm | PVC-free | 25/50* |
| 60A1225 | Priming volume 4,8 ml, length 600 cm | PVC-free | 25/50* |

| | |
|--|--|
| Material: Tube: PE Connector fLL/mLL: PE Protection cap: ABS, PE Latex-free | Sterilization: EO Single use, sterile 1-pcs in peel pack * = minimum ordering quantity |
|--|--|