



Indication: For fixation around joints and trunks which are difficult to bandage. Can also be used for fixation of primary dressings and tubings. Only to be used on undamaged skin.

## Fixation bandage

MEDIPLAST offers DANALASTIC universal bandage to be used for light support or fixation. The bandage is blood and salves repellent.

Elasticity: 70%  
 Cotton, Polyamide  
 Length: stretched

Art.No 5 meters	Art.No 4 meters	Description	Ward box	Trsp.box
6050300	6050340	4 cm	10	720
6050301	6050341	6 cm	10	480
6050302	6050342	8 cm	10	320
6050303	6050343	10 cm	10	280
6050304	6050344	12 cm	10	240

Manufactured in Denmark

Ward box = min order qty