



Esmark Bandages, Latex free

MEDIPLAST® offers **ESMARK BANDAGES**. A soft and smooth synthetic rubber bandage. The bandages offer maximum and consistent pressure to extremities for exsanguination prior to surgery. The product comes in three different widths and two different lengths. Sterile product intended for single use.

ESMARK BANDAGES: Synthetic Polyisoprene

Art.No	Description	Ward box	Trsp.box
60702010	8 x 274 cm	10	150
60702011	10 x 274 cm	10	150
60702012	15 x 274 cm	10	70
60702013	8 x 365 cm	10	150
60702014	10 x 365 cm	10	150
60702015	15 x 365 cm	6	42

ESMARK BANDAGES – POWDER FREE: Synthetic Polyisoprene

Art.No	Description	Ward box	Trsp.box
60702016	8 x 274 cm	10	150
60702017	10 x 274 cm	10	150
60702018	15 x 274 cm	10	70
60702019	8 x 365 cm	10	150
60702020	10 x 365 cm	10	150
60702021	15 x 365 cm	6	42

Ward box = min order qty