

# Infusion Set

- Different length of tubing between 180-320 cm
- Vented or non-vented
- PVC-free or PVC DEHP-free material
- Roller clamp for maximum precision to ensure consistent and safe infusion rate
- Spike without silicone to ensure safe connection to infusion bag/bottle
- Liquid filter and air filter
- Flowstop with hydrophobic filter which simplifies priming and reduces the risk of leakage and contamination
- Self-ejecting rotating luer-lock
- Available with 3-way stopcock, back-check valve and/or needle-free connector
- Air-Guard prevents air passing into the tube when the infusion has run out
- Flow-Regulator for accurate flow
- MRI compatible



## Infusion Set, non-vented

Art. no.	Description		Ward/Trsp.box
60MR60500152	Infusion Set 180 cm, non-vented	PVC-free	10/100*
60MR60500152W	Infusion Set 180 cm, non-vented	PVC-free	10/100*
60MR60500150T	Infusion Set 180 cm, non-vented	DEHP-free PVC	10/100*
60MR60500182	Infusion Set 220 cm, non-vented	PVC-free	10/100*
60MR60500160	Infusion Set 320 cm, non-vented	PVC-free	10/100*
60MR60500151T	Infusion Set 320 cm, non-vented	DEHP-free PVC	10/100*

## Infusion Set, vented

Art. no.	Description		Ward/Trsp.box
60MR60500156	Infusion Set 180 cm, vented	PVC-free	10/100*
60MR60500153	Infusion Set 180 cm, vented	PVC-free	10/100*
60MR60500154T	Infusion Set 180 cm, vented	DEHP-free PVC	10/100*
60SRL60500154	Infusion Set 180 cm, vented with back-check valve	PVC-free	10/100*
60SRL60500153K	Infusion Set 180 cm, vented, y-connection with needle-free connector	PVC-free	10/100*
60SRL60500153KV	Infusion Set 180 + 12 cm, vented with needle-free connector and back-check valve	PVC-free	10/100*
60SRL60500156C	Infusion Set 180 cm, vented with closed male luer adapter and priming cap	PVC-free	50/50*
60MR60500186	Infusion Set 220 cm, vented	PVC-free	10/100*
60SRL60500186	Infusion Set 250 cm, vented	PVC-free	10/100*
60MR60500161	Infusion Set 320 cm, vented	PVC-free	10/100*
60MR60500162	Infusion Set 320 cm, vented	PVC-free	10/100*
60MR60500155T	Infusion Set 320 cm, vented	DEHP-free PVC	10/100*

## Infusion Set, non-vented, with 3-way stopcock\*\* and extension

Art. no.	Description		Ward/Trsp.box
60MR60500157	Infusion Set 180 cm, non-vented, 3-way stopcock + 15 cm extension	PVC-free	10/100*
60MR60500141T	Infusion Set 180 cm, non-vented, 3-way stopcock + 15 cm extension	DEHP-free PVC	10/100*
60MR60500183	Infusion Set 220 cm, non-vented, 3-way stopcock + 15 cm extension	PVC-free	10/100*

## Infusion Set, vented, with 3-way stopcock\*\* and extension

Art. no.	Description		Ward/Trsp.box
60MR60500138	Infusion Set 180 cm, vented, 3-way stopcock + 15 cm extension	PVC-free	10/100*
60SRL60500139	Infusion Set 180 cm, vented, 3-way stopcock + 15 cm extension	PVC-free	10/100*
60MR60500137	Infusion Set 180 cm, vented, 3-way stopcock + 15 cm extension	PVC-free	10/100*
60MR60500140T	Infusion Set 180 cm, vented, 3-way stopcock + 15 cm extension	DEHP-free PVC	10/100*
60MR60500146	Infusion Set 180 cm, vented with back-check valve, 3-way stopcock + 15 cm extension	PVC-free	10/100*
60MR60500148W	Infusion Set 180 cm, vented with back-check valve, 3-way stopcock + 15 cm extension	PVC-free	10/100*
60SRL60500148P	Infusion Set 180 cm, vented with back-check valve, 3-way stopcock with needle-free connector + 15 cm extension	PVC-free	10/100*
60MR60500140SG	Infusion Set 180 cm, vented, 2 x 3-way stopcock + 15 cm extension	PVC-free	10/50*
60MR60500140SU	Infusion Set 180 cm, vented, 2 x 3-way stopcock + 15 cm extension	PVC-free	10/50*
60MR60500187	Infusion Set 220 cm, vented, 3-way stopcock + 15 cm extension	PVC-free	10/100*
60MR60500145	Infusion Set 320 cm, vented, 3-way stopcock + 15 cm extension	PVC-free	10/100*
60MR60500147	Infusion Set 320 cm, vented with back-check valve, 3-way stopcock + 15 cm extension	PVC-free	10/100*

## Infusion Set with Flow-Regulator, vented

Art. no.	Description		Ward/Trsp.box
60MR60500156F	Infusion Set 180 cm, vented with Flow-Regulator	PVC-free	10/100*
60MR60500156FV	Infusion Set 180 cm, vented with Flow-Regulator and back-check valve	PVC-free	10/100*
60SRL60500156F	Infusion Set 180 cm, vented with Flow-Regulator and Air-Guard	PVC-free	10/100*
60MR60500161F	Infusion Set 320 cm, vented with Flow-Regulator	PVC-free	10/100*

## Infusion Set with Air-Guard, vented

Art. no.	Description		Ward/Trsp.box
60SRL60500153	Infusion Set 180 cm, vented with Air-Guard	PVC-free	10/100*
60SRL60500155AG	Infusion Set 180 cm, vented with Air-Guard and back-check valve	PVC-free	10/100*
60SRL60500156F	Infusion Set 180 cm, vented with Air-Guard and Flow Regulator	PVC-free	10/100*
60SRL60500186AG	Infusion Set 250 cm, vented with Air-Guard	PVC-free	10/100*
60SRL60500162AG	Infusion Set 320 cm, vented with Air-Guard	PVC-free	10/100*

Material:  
Tube: PVC DEHP-FREE/TPE  
Drip chamber: SB/TPE  
Spike: ABS  
Air filter: PP 0,2 or 0,8 µm  
Liquid filter: PA 15 µm  
Roller clamp: ABS/PS  
Pinch clamp: PP

Priming-stop/Flowstop: ABS/PE  
Back-check valve: SAN/MABS/Silicone  
Flow-Regulator: PC / POM  
Connector/Luer-lock: ABS  
3-way stopcock: PC/HDPE  
Needle-free connector: PC/Silicone  
Closed male luer adapter: PSU/PE  
Hydrophobic cap: PP/Acrylic

Sterilization: EO  
Infusion Set 20 drp aqua dest = 1 ml  
Single use, sterile 1 pcs in peel pack  
\* = minimum ordering quantity  
Priming volume: 0,07 ml/ cm Ø 2,9 mm  
\*\*blood, drug and lipid resistant 72 h.