




Product catalogue



Revision: July 2018



Återförsäljare:

MEDI PLAST

Mediplast AB

Box 1004, 212 10 Malmö

T 020-78 80 35

mediplast.info@mediplast.com

www.mediplast.com

Medline Surgical Gowns

Performance, Variety and Quality



Performance with Comfort

We have designed our surgical gowns to enhance the performance of your surgical team. All Medline surgical gowns are generously cut in the chest and sleeves for greater comfort and freedom of movement. Each gown features a hook & loop neck closure for an easily adjustable, gap-free custom fit. Medline's wide range of sizes and lengths satisfies the needs of every person in the O.R. with high-tech breathable fabrics keeping the wearer comfortable and dry, even during long procedures.

Abundant Choices

You want options in style, fabric and reinforcements to satisfy any preference in the OR. Medline's four gown families each feature their own distinct fabric and styles. With Medline, clinicians have the widest available range of choice in size, length, reinforcement type, and sleeve construction.

Dependable Protection

Medline's gown fabrics offer incredible protection. Low-linting and latex-free, each gown meets or exceeds performance recommendations of the EN13795* surgical gown standard. The Essential and Advanced materials' high hydrostatic head gives wearers the protection they need and the comfort they deserve. The unique Ultimate gown offers all-over impervious protection with exceptional breathability and coolness.

Colour Coding in accordance to the EN13795 standard


All Medline Surgical Gowns feature intuitive colour coding in accordance to the EN13795* standard. Product labels, neck bindings and transfer tabs display Purple for Standard Performance or Blue for High Performance. Droplets tie back to the US standard AAMI levels 1 to 4**. These colours, as well as the droplets on the gown itself, make it easy to choose the ideal gown for the individual requirements of any surgical procedure.

STANDARD PERFORMANCE SURGICAL GOWN


Non-reinforced or fabric-reinforced for **low to medium risk** of exposure

PURPLE and BLUE HELP YOU CHOOSE THE RIGHT GOWN – EVERY TIME

Colour-coding of the gowns, package labels and transfer cards indicate the level of protection from bodily fluids. Available in a variety of materials: Eclipse®, Essential and Advanced. Now it's easy to choose the right gown for every procedure.



Package Label
Standard Performance Gown



Package Label
High Performance Gown

EN13795 Compliant
Medline colour-coded surgical gowns are 100% compliant with EN13795 directives. Developed by the European Committee for Standardisation and its Technical Committee, EN13795 provides guidance on the characteristics of surgical gowns and drapes for protection from transmission of microorganisms during invasive procedures.

HIGH PERFORMANCE SURGICAL GOWN

Poly-reinforced for **high risk** of exposure

To optimize gown usage in the Operating Room, Medline has developed the OPS™ allowing you to “optimize” your gown usage according to your individual clinical needs.

Optimized Protection System for surgical gowns: OPS™

Based on extensive clinical feedback, we have identified Comfort, Protection and Cost as the most important criteria to consider when selecting a surgical gown.

- A comfortable surgical gown is desired when performing longer procedures.
- An increased protection level is necessary as fluid challenges increase during surgery.
- For short and/or low fluid intensive procedures, customers are willing to sacrifice some degree of comfort for the right cost.

Every hospital has its unique combination of clinical requirements based on the procedures it performs. For this reason, Medline has created the Optimized Protection System™ - providing the ideal portfolio to make choosing the right gown simpler and more cost effective.

OPS™ Essential


Adequate protection at lowest cost

- Light weight material, sufficient for shorter and/or low fluid procedures such as ear and mouth surgery, small orthopaedic and general surgeries etc
- Comfortable - Softened SMMS technology offers exceptional softness and breathability
- When cost is of highest concern

OPS™ Advanced

No compromise between comfort and protection

- Breathable and comfortable – quickly transfers heat away to keep the wearer cool and comfortable
- Excellent protection- excellent barrier performance with high Hydrostatic head
- Optional raglan cut- provides roominess and facilitates easier movement
- Ideal for the majority of procedures in all surgical areas



OPS™ Ultimate

The intelligent gown with breathable membrane technology

- EN13795* High performance across the entire gown
- Fully impervious to blood and blood borne pathogens (AAMI** PB70, PB71 tested)
- Breathable and cool - reacts to increasing body temperature by increasing water vapour transmission
- Ideal for all procedures where breathable all-over protection is needed such as long Cardio, Orthopaedic, Transplant surgeries, etc.

Eclipse®

A soft and comfortable gown

- Outstanding comfort and breathability - textile like appearance
- Sufficient protection



Right Protection ... Right Comfort ... Right Cost

*The EN 13795 is the European harmonized standard specifying requirements for single-use and reusable covering (i. e. surgical gowns, surgical drapes) used as medical devices for patients, clinical staff and equipment and intended to prevent the transmission of infective agents between patients and clinical staff during invasive surgical procedures.

**AAMI: Association for the Advancement of Medical Instrumentation. The AAMI PB70 is the US reference standard establishing a system of classification for protective apparel and drapes based on their liquid barrier performance and specifies related labeling requirements and standardized test methods.



Medline Surgical Drapes

Designed, Tested and Trusted By OR Professionals

Medline® Surgical Drapes were designed based on the feedback from OR Professionals across the entirety of Europe. We continue to use such feedback as the basis for introducing new products, as well as improving features and designs of existing products. Medline® surgical drapes offer the best features available to the market such as the innovative biodegradable EcoDrape™ fabric, and the clever Panther Paw™ hook & loop tube holder. Additionally, all Medline® Surgical Drapes are 100% Latex-Free and fully compliant with the EN13795* standard.

Choice

Medline® is the market's only surgical drapes and packs provider to offer a wide variety of fabrics. Medline® surgical drapes are available in the classic Eclipse® spunlace material, as well as in other cutting edge fabrics such as absorbent bi-laminate, light and drapeable SMS, or top-of-the-line tri-laminate options.

Absorbency

Medline® specialty surgical drapes are all manufactured with oversized absorbent areas. This ensures the surgical field is kept as clean as possible, and reduces post-operative clean-up time. The bi-laminate and tri-laminate fabrics provide absorbency across the whole surface of the drape. On repellent spunlace and SMS fabrics, an oversized heavy-duty absorbent reinforcement is applied to the critical areas.

Drapeability and Strength

Medline® surgical drape fabrics are designed to provide high drapeability for facilitated efficient patient draping. In addition to high drapeability, all materials are made from robust non-woven fabrics and films to ensure sterility of the surgical field during complex procedures.

Colour

Medline® specialty surgical drapes feature an intuitive labeling system complete with detailed product images and specifications. Product labels are printed in various colours based on the material being used. This aids the surgical team in selecting the correct product best fitted for each specific need.

Medline Optimized Protection System: OPS™


Protection, cost and standardization are the most important criteria when selecting a surgical drape. Medline OPS™ offer three different types of drapes to suit your individual clinical needs and “optimize” your drape usage.

- **OPS™ Essential** – For short and/or low-fluid procedures. Our most economical, yet effective drape.
- **OPS™ Advanced** – For medium fluid procedures, with absorbent reinforcement around critical areas and a soft contact with patient skin.
- **OPS™ Ultimate** – When increased protection is necessary to accommodate fluid increases during longer and more complex procedures. Combines strength, absorption and comfort.

When you use Medline’s OPS™, you can rely on:

- Appropriate protection for the type and length of the procedure
- Cost efficiency by using the appropriate drape based on the procedure
- Greater standardisation by using the same drape for the same procedure every time

Choose From Three Levels of Fluid Protection to Suit Your Procedure



Adequate protection for short procedure

- ▶ Drapes are made of bi-laminate fabric
- ▶ Meet the High Performance (HP)* level
- ▶ Fully impervious to blood and bloodborne pathogens**
- ▶ Sufficient for short and low-fluid procedures

The right choice for longer procedure with more fluids

- ▶ Drapes are made of SMS fabric
- ▶ Meet the High Performance (HP)* level for the reinforced drapes
- ▶ Fully impervious to blood and bloodborne pathogens** on reinforced area
- ▶ Specialty pouches and reinforcements are available

Maximum protection for longer, more complex procedures

- ▶ Drapes are made of tri-laminate material
- ▶ Meet the High Performance (HP)* level
- ▶ Fully impervious to blood and bloodborne pathogens**
- ▶ Recommended for long and high-fluid procedures
- ▶ Specialty pouches and reinforcements are available

Fluid protection levels according to standards.

Each OPS™ drape is colour-coded and the product label identifies Standard Performance (SP) level or High Performance (HP) level per the EN13795 standard*. The droplets correspond to the 4 AAMI classification levels**.

Protection	AAMI**	EN13795*
6 droplets	4	HP
5 droplets	3	SP/HP
3 droplets	2	SP

Table of contents

Drapes & Packs

Angiography	10
Cardiosurgery	17
ENT Surgery	26
General Surgery	40
Lap Drapes and Packs	41
Lap/Gyn Drapes and Packs	46
Thyroid Drapes and Packs	52
Universal Drapes and Packs	54
Gynaecological Surgery	69
Baby Blanket	70
Obstetrical Drapes and Packs	71
Peri-Gynae Packs	80
Speciality	87
NeuroSurgery	90
Ophtalmic Surgery	95
Orthopedic Surgery	101
Arthroscopy Drapes and Packs	102
Extremity Drapes and Packs	109
Hand Drapes and Packs	127
Hip Drapes and Packs	131
Shoulder Drapes and Packs	136
Split Drapes and Packs	142
Paediatric Surgery	161
Other Drapes and Packs	161
Split Packs	164
Universal Drapes and Packs	165
Supplementary Program	167
Adhesive drape sheets	167
Cotton Towel	185
Drape sheets	188
Fenestrated Drapes	197
Fenestrated Drapes w. Adh.	199
Leggings	206
Perineal Cover	208

Table of contents

Pouches	209
Protective Sleeves	213
Ring Cover	215
Stockinettes	217
Tapes/Velcro ´s	220
Tourniquet cover	223
Towels	224
Tube pocket	225
Under-Buttocks Drape	226
Urological Surgery	229
T.U.R. Drapes and Packs	229
Urology packs	232

Equipment Covers

Back Table Cover	236
Banded Bags	241
C-Arm Cover	246
Grip Cover	253
Mayo Stand Cover	254
Microscope Cover	257
Other Covers	262
Probe Covers	268
Video Camera Cover	270

Invisishield

Clear plastic drapes	274
Incise drapes	280
Instrument and fluid collection pouches	284

Drapes & Packs

Drape sheet 112x145 cm with adhesive

Item no.	29225CEB	Description	Dimensions	Quantity
eClass-Nr.		spunlace drape 112x145 cm, adhesive		1
Material	Eclipse - Spunl			
PU	40			
SU	0			



Drape Sheet, Three Quarter

Item no.	ES9349CE	Description	Dimensions	Quantity
eClass-Nr.		Three quarter drape	145x193cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	24			
SU	0			



Drape with adhesive


Item no.	ES933CE	Description	Dimensions	Quantity
eClass-Nr.		Drape with adhesive	300x221 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	24			
SU	0			



Drapes & Packs

Ophthalmology Pack

Item no.	DYJPEOPHPSM3	Description	Dimensions	Quantity
eClass-Nr.		Eye drape, fen. w/ incise film 11x13 cm, 2 pouches	165x178 cm	1
Material	ADV - SMMMS	Back table cover, 50µm	112 x 110 cm	1
PU	15	SMS wrap		1
SU	0			




POUCH, INSTRUMENT, PE, ADH, 16, 5X75CM

Item no.	96880CE	Description	Dimensions	Quantity
eClass-Nr.				
Material	PE			
PU	230			
SU	555			



Underbuttocks Drape with Pouch and Drainage Port

Item no.	E81684CS	Description	Dimensions	Quantity
eClass-Nr.		under buttocks drape, 51x23 cm reinforcement, graduated fluid col..	102x113 cm	1
Material	PE, Cellulose			
PU	80			
SU	0			




Drapes & Packs

Angiography

Angiographie drape

Item no.	DYJPEFEADSM	Description	Dimensions	Quantity
eClass-Nr.		Wrap	100x100 cm	1
Material	ADV - SMMMS	Femoral angiography drape, 2 fen. diam. 15 cm w/ incise film	216x325 cm	1
PU	5			
SU	0			




Angiography drape

Item no.	TBAG202CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-04	Angiography drape w/ adh, reinf and fen	221x290 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	12			
SU	0			



Angiography drape

Item no.	SPT-ESAG354	Description	Dimensions	Quantity
eClass-Nr.		angiography drape	220x350cm	1
Material	ESS - 2 layers			
PU	10			
SU	0			




Drapes & Packs

Angiography


Angiography drape with 2 clear pannels with pouch

Item no.	TBAG252CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-04	Angiography drape, 2 clear pannels with pouch	221x290 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	12			
SU	0			




Angiography drape with one clear panel

Item no.	TBAG292CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-04	Angiography drape, 1 clear panel	226x290 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	12			
SU	0			



Angiography drape, femoral, two clear panels

Item no.	TBAG352CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-04	Angiography drape, 2 fenestrations and two clear PE panels	221x343 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	12			
SU	0			



Drapes & Packs

Angiography

Angiography Pack

Item no.	29464FRB	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-04	CSR wrap	85x85 cm	1
Material	Ultimate - 3 la	Absorbent Hand towels	20x30 cm	2
PU	7	OR gown, Eclipse, XL	130 cm	1
SU	0	Hand towels	30x20 cm	1
		Angiography drape, 1 clear panel	226x290 cm	1
		Back table cover	140x190 cm	1

No picture available

Angiography pack

Item no.	29463FRC	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-04	Back table cover	140x190 cm	1
Material	Ultimate - 3 la	CSR wrap	95x95 cm	1
PU	8	Absorbent Hand towels	20x30 cm	2
SU	0	Drape angio femoral, w/2fen (15cm), INCM NS	221 x 290cm	1

No picture available

Angiography pack

Item no.	TBAG392CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-04	Angiography drape, w/ 2 fenestrations, 1 clear panel	241x343 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	8			
SU	0			




Drapes & Packs

Angiography


Angiography pack

Item no.	CP29462CE	Description	Dimensions	Quantity
eClass-Nr.		Hand towels	30x39 cm	4
Material	Ultimate - 3 la	Angiography drape, 1 clear panel	226x290 cm	1
PU	8	Back table cover	140x190 cm	1
SU	0			




Angiography pack

Item no.	CP29046CE	Description	Dimensions	Quantity
eClass-Nr.		Absorbent Towel	39x58 cm	4
Material	Ultimate - 3 la	Banded bag, L-50cm, D-100cm, circular	50x100 cm	1
PU	6	Angiography drape, 2 clear pannels with pouch	221x290 cm	1
SU	0	Back table cover	140x190 cm	1



Angiography pack

Item no.	CPU29051CE	Description	Dimensions	Quantity
eClass-Nr.		Hand towels	30x20 cm	2
Material	Ultimate - 3 la	Back table cover, 50µm	120x140 cm	1
PU	6	OR gown, Essential, XL	130 cm	2
SU	0	Banded bag with sewn elastic, bordeaux	A=Ø56 cm, B=40	2
		Banded bag with sewn elastic, brown	95x85 cm	1
		Angiography drape, w/ 2 fenestrations, 1 clear panel	241x343 cm	1




Drapes & Packs

Angiography


Angiography pack

Item no.	CPU29053CE	Description	Dimensions	Quantity
eClass-Nr.		Adhesive tape	9x50 cm	1
Material	Ultimate - 3 la	Angiography drape, 2 fenestrations and two clear PE panels	221x343 cm	1
PU	0	Half drape	112x145 cm	2
SU	0	Back table cover, 50µm	120x140 cm	1
		Banded bag with sewn elastic, brown	95x85 cm	1




Brachial angiography drape

Item no.	TB29458CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-04	Brachial angiography drape w/adh and fen	97x112 cm	1
Material	Ultimate - 3 la	CSR wrap	60x60 cm	1
PU	50			
SU	0			



Brachial angiography drape

Item no.	DYNJPE4106	Description	Dimensions	Quantity
eClass-Nr.		Wrap	60x60 cm	1
Material	Eclipse - Spunl	Brachial angiography drape, 7 cm diam. adh. fen. w/ incise film	97x112 cm	1
PU	40			
SU	0			




Drapes & Packs

Angiography


Central Line Drape

Item no.	DYNJPE4121SM	Description	Dimensions	Quantity
eClass-Nr.		Central Line Drape, 11.5 cm diam. adh. fen., 46x46 cm clear windo..	178x253 cm	1
Material	ADV - SMMMS	SMS wrap		1
PU	24			
SU	0			




Femoral angiography drape

Item no.	DYNJPE4113	Description	Dimensions	Quantity
eClass-Nr.		Femoral angiography drape	201x279 cm	1
Material	Eclipse - Spunl	SMS wrap		1
PU	12			
SU	0			



Femoral angiography drape

Item no.	TBAG3521CE	Description	Dimensions	Quantity
eClass-Nr.		Aniography drape, 2 fenestrations and two clear PE panels	196x343 cm	1
Material	Ultimate - 3 la	SMS wrap		1
PU	12			
SU	0			




Drapes & Packs

Angiography


Interventional cardiology pack

Item no.	TB29560CE	Description	Dimensions	Quantity
eClass-Nr.		Side sheet, adhesive	75x90 cm	1
Material	Ultimate - 3 la	Angiography drape	226/243x352 cm	1
PU	6	Hand towels	30x20 cm	2
SU	0	CSR wrap	86x86 cm	1
		Back table cover, reinforced, 50µm	140x190 cm	1



PICC line drape

Item no.	DYNJPE4122SM	Description	Dimensions	Quantity
eClass-Nr.		SMS wrap		1
Material	ADV - SMMMS	PICC Line Drape, 2 fenestration diam.11.5 cm	163/284.5x264 c	1
PU	12			
SU	0			




Drapes & Packs

Cardiosurgery


Bootie pack

Item no.	DYJPEBOPSM	Description	Dimensions	Quantity
eClass-Nr.		Wrap	61x61 cm	1
Material	ADV - SMMMS	Booties with 5x10 cm adhesive tape	32x43 cm	2
PU	30			
SU	0			




Cardiac collar drape

Item no.	TB27007CE	Description	Dimensions	Quantity
eClass-Nr.		Cardiac collar drape		1
Material	Ultimate - 3 la	SMS wrap		1
PU	60			
SU	0			



Cardiac Pack

Item no.	29C285FRA	Description	Dimensions	Quantity
eClass-Nr.		Adhesive tape, PE	8x40 cm	6
Material	Eclipse - Spunl	Adhesive drape	75x91 cm	4
PU	3	Perineal sheet, adhesive	16x37 cm	1
SU	0	Drape, w/adh, reinf heavy	196x340 cm	1
		Side drape, reinforced, adhesive	98x110 cm	2
		Drape sheet, reinf, adhesive	290x240 cm	1
		Drape Sheet	145x240 cm	1
		CSR wrap	112x112 cm	1
		Back table cover, folded, 90 µm	203x240 cm	1
		Back table cover, 75µm	160x240 cm	1




Drapes & Packs

Cardiosurgery

Cardiac surgery pack


Item no. 29288CEA
eClass-Nr. 34-35-01-03
Material ECL - Spunlac
PU 2
SU 0



Description	Dimensions	Quantity
Drape, adhesive	35x16 cm	1
Side drape w/adh, reinf, light	98x110 cm	8
Back table cover, folded, 50µm	140x190 cm	1
Head Drape, w/ adh, velcro	240x196 cm	1
Drape, w/adh, reinf heavy	196x340 cm	1
Drape, adhesive	180x98 cm	2
Adhesive tape	9x50 cm	6
Pouch, Instmt, double, PE	33x41 cm	1
Back table cover, 90µm	160x240 cm	1
CSR wrap	112x112 cm	1

Cardiovascular drape


Item no. DYJPECADSM
eClass-Nr.
Material ADV - SMMMS
PU 10
SU 0



Description	Dimensions	Quantity
Wrap	100x100 cm	1
Cardiovascular drape, fen. 28x38 cm with full incise film	196/300x349 cm	1

Cardiovascular pack

Item no. DYJPECAPSM1
eClass-Nr.
Material ADV - SMMMS
PU 4
SU 0




Description	Dimensions	Quantity
Peri-groin drape, adh. split 45x92 cm, elastic perineal cover	196x311 cm	1
Booties with 5x10 cm adhesive tape	32x43 cm	2
Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
Back table cover,	152x190 cm	1
Cardiovascular drape, fen. 28x38 cm with full incise film	196/300x349 cm	1

Drapes & Packs

Cardiosurgery


Cardiovascular pack

Item no.	ES29562CE	Description	Dimensions	Quantity
eClass-Nr.		Perineal sheet, adhesive	16x37 cm	1
Material	ESS - 2 layers	Adhesive tape	9x50 cm	2
PU	7	Drape adhesive	75x75 cm	1
SU	0	Instrument pouch in clear plastic with adhesive	33x41cm	1
		Absorbent Hand towels	20x30 cm	6
		Back table cover	152x229 cm	1
		SMS wrap, large	90x90 cm	1
		Mayo stand cover, ultra absorbent reinforcement 75x100 cm	80x145 cm	1
		Cardiovascular drape, 1 fen. with incise film, 38x65cm	200/290x338 cm	1



Cardiovascular pack

Item no.	ES29561CE	Description	Dimensions	Quantity
eClass-Nr.		Instrument pouch in clear plastic with adhesive	33x41cm	1
Material	ESS - 2 layers	Peri-groin drape, adhesive split, perineal drape, reinforcement	200x260 cm	1
PU	4	Adhesive tape	9x50 cm	3
SU	0	Drape Booties w/ adh	29x33 cm	2
		Drape Booties w/ adh	29x33 cm	2
		Drape adhesive	75x75 cm	2
		Back table cover, folded, 50µm	140x190 cm	1
		Back table cover	152x229 cm	1
		SMS wrap, large	90x90 cm	1
		Absorbent Hand towels	20x30 cm	6
		cardiac T drape, 32x40 cm fenestration with full incise film, 50.5x1..	200/290x338cm	1
		Mayo stand cover, ultra absorbent reinforcement 75x100 cm	80x145 cm	1



Cardiovascular pack,

Item no.	29209FRC	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-03	Perineal sheet, adhesive	16x37 cm	1
Material	ECL - Spunlac	Adhesive tape, PE	8x40 cm	4
PU	3	Drape, reinforced, adhesive	240x196 cm	1
SU	0	Instrument pouch, PE, adhesive	30,5x29 cm	1
		Back table cover, 90µm	160x240 cm	1
		CV split drape, 236/254x312cm	236/254x312cm	1
		Hand towels	30x39 cm	2
		Back table cover, 75µm, folded	140x190cm	1
		Drape	98x98 cm	2
		CSR wrap	95x95 cm	1


No picture available

Drapes & Packs

Cardiosurgery


Cardiovascular split drape

Item no.	TB29154	Description	Dimensions	Quantity
eClass-Nr.		Drape CVmU w/ adh	254x312 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	8			
SU	0			




CPB pack

Item no.	CP29225CE	Description	Dimensions	Quantity
eClass-Nr.		Drape	98x98 cm	1
Material	Ultimate - 3 la	Foot sheet w/ adh	183 x 191 cm	1
PU	3	Instrument pouch, PE, adhesive	30,5x29 cm	1
SU	0	Perineal sheet, adhesive	16x37 cm	1
		Hand towels	30x39 cm	3
		Drape sheet, adhesive	98x98 cm	1
		Back table cover, folded, 90 µm	203x240 cm	1
		Adhesive tape, PE	8x40 cm	4
		Drape, adhesive	49x65 cm	2
		Polyrayon panel	49x50 cm	2
		CV split drape, 236/254x312cm	236/254x312cm	1
		Back table cover, 90µm	160x240 cm	1
		CSR wrap	60x60 cm	1



Pacemaker drape

Item no.	DYNJPE4107	Description	Dimensions	Quantity
eClass-Nr.		Wrap	76x76 cm	1
Material	ECL - Spunlac	Pacemaker drape, 2 rectangular fen. 15x18 cm w/ fen. incise film, c..	208x292 cm	1
PU	20			
SU	0			



Drapes & Packs

Cardiosurgery


Pacemaker Drape

Item no.	ES10020CE	Description	Dimensions	Quantity
eClass-Nr.		Pacemaker Drape, 15x15 cm adhesive fenestration	200x280 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	24			
SU	555			




Pacemaker Drape

Item no.	ES9409CE	Description	Dimensions	Quantity
eClass-Nr.		Pacemaker Drape	262x360cm	1
Material	ESS - 2 layers	SMS wrap, large	90x90 cm	1
PU	14			
SU	0			



Pacemaker pack

Item no.	TB30020CE	Description	Dimensions	Quantity
eClass-Nr.		Pacemaker drape, adhesive fenestration	220x240 cm	1
Material	Ultimate - 3 la	Adhesive tape, PE	8x40 cm	2
PU	10	Mayo stand cover, reinforced	80x142 cm	1
SU	0	Back table cover, reinforced, 50µm	140x190 cm	1
		Hand towels	30x20 cm	2




Drapes & Packs

Cardiosurgery


Pacemaker pack

Item no.	TB9408CE	Description	Dimensions	Quantity
eClass-Nr.		Pacemaker drape, 15x15 fenestration with incise film	210x316 cm	1
Material	Ultimate - 3 la	Back table cover	140x190 cm	1
PU	8	CSR wrap	86x86 cm	1
SU	0			




Peri-Groin drape

Item no.	TB29500CE	Description	Dimensions	Quantity
eClass-Nr.		Peri-Groin drape, split, perineal drape, reinforcement	196x240 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	12			
SU	0			



Peri-groin drape

Item no.	DYJPEPERDSM	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	ADV - SMMMS	Peri-groin drape, adh. split 45x92 cm, elastic perineal cover	196x311 cm	1
PU	6			
SU	0			




Drapes & Packs

Cardiosurgery

Peri-groin drape

Item no.	DYJPEPERDEC	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	Eclipse - Spunl	Peri-groin drape, adh. split 45x92 cm, adh. perineal cover w/ elastic	201x236 cm	1
PU	6			
SU	0			




Peri-groin Drape

Item no.	SPTPEPERDSM	Description	Dimensions	Quantity
eClass-Nr.		Peri-groin drape, adh. split 45x92 cm, elastic perineal cover	196x311 cm	1
Material	ADV - SMMMS			
PU	20			
SU	0			



TAVI drape

Item no.	TB9159CE	Description	Dimensions	Quantity
eClass-Nr.		SMS wrap		1
Material	Ultimate - 3 la	TAVI drape, 3 fenestration w/full incise film, 2 clear panels, reinf	240/290x448 cm	1
PU	0			
SU	555			




Drapes & Packs

Cardiosurgery


TAVI pack

Item no.	TB29159CE	Description	Dimensions	Quantity
eClass-Nr.		Back table cover, reinforced, 50µm	140x190 cm	1
Material	Ultimate - 3 la	Drape adhesive	75x75 cm	4
PU	5	TAVI drape, 3 fenestration w/full incise film, 2 clear panels, reinf	240/290x448 cm	1
SU	555	SMS wrap		1




Thoracotomy drape,

Item no.	TB29420CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-03	Drape with adhesive fenestration 36x28cm, armboard covers	254/196x310 cm	1
Material	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU	15			
SU	0			



Varicose Vein Pack

Item no.	ES10120CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-03	Hand towels	30x20 cm	2
Material	ESS - 2 layers	CSR wrap	86x86 cm	1
PU	8	Adhesive tapes	9x50 cm	1
SU	0	Drape	200x150 cm	1
		Split sheet, adhesive	150x230 cm	1
		Head sheet, adhesive	260x160 cm	1
		Back table cover	140x190 cm	1




Drapes & Packs

Cardiosurgery

Varicose vein pack,

Item no. TB30120CE
eClass-Nr. 34-35-01-03
Material Ultimate - 3 la
PU 8
SU 0


Description	Dimensions	Quantity
Absorbent Hand towels	20x30 cm	1
Velcro strip	2.5x13 cm	1
Perineal sheet, adhesive	15 x 37 cm	1
U-drape, split 18x60cm, adhesive	180x240 cm	1
Head sheet, adhesive	180x240 cm	1
Back table cover	140x190 cm	1
CSR wrap	86x86 cm	1



Varicose veins pack

Item no. ES29120CE
eClass-Nr.
Material ESS - 2 layers
PU 8
SU 0

Description	Dimensions	Quantity
Absorbent Hand towels	20x30 cm	1
Perineal sheet, adhesive	16x37 cm	1
Hook and loop tube holder	14x2.5 cm	1
Back table cover, reinforced, 50µm	140x190 cm	1
Head drape with adhesive	240X180 CM	1
Split drape, adhesive, reinforcement	240x183 cm	1
SMS wrap		1




Drapes & Packs

ENT Surgery


Basic face pack

Item no.	TB291060CE	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
Material	Ultimate - 3 la	Adhesive tapes	9x50 cm	1
PU	6	Back table cover	140x190 cm	1
SU	0	ENT drape 196x300cm, 12.5x15 cm adh. fen.	196x300 cm	1
		OR gown, Eclipse, XLL	150 cm	1



Basic face pack

Item no.	TB30460CE	Description	Dimensions	Quantity
eClass-Nr.		Absorbent Hand towels	20x30 cm	4
Material	Ultimate - 3 la	Adhesive tapes	9x50 cm	1
PU	10	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
SU	555	ENT drape 196x300cm, 12.5x15 cm adh. fen.	196x300 cm	1
		Back table cover	140x190 cm	1



Dermatology pack

Item no.	ES30249CE	Description	Dimensions	Quantity
eClass-Nr.		SMS wrap		1
Material	ESS - 2 layers	Dermatology drape, 12 cm diam. adh. fen.	100x150 cm	1
PU	36	Back table cover, reinforced, 50µm	140x190 cm	1
SU	0			



Drapes & Packs

ENT Surgery

Drape with fenestration

Item no.	TB30668CE	Description	Dimensions	Quantity
eClass-Nr.		Drape with 10cm dia fenestration, with soft adhesive tape	150x98 cm	1
Material	Ultimate - 3 la	CSR wrap	60x60 cm	1
PU	42			
SU	0			




Drape, fenestration diameter 10cm, adhesive

Item no.	TB30669CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Drape sheet w/ adhesive fen. diam.10cm	98x150 cm	1
Material	Ultimate - 3 la	CSR wrap	60x60 cm	1
PU	42			
SU	0			



EENT drape

Item no.	DYNJPE7012SM	Description	Dimensions	Quantity
eClass-Nr.		Wrap	80x80 cm	1
Material	ADV - SMMMS	EENT drape, instrument pad 41x41 cm w/ adh. fen. 5x6 cm	203x229 cm	1
PU	10			
SU	0			




Drapes & Packs

ENT Surgery


EENT pack

Item no.	DYJPEENTPSM	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
Material	ADV - SMMMS	Head turban drape with adhesive	102x69 cm	1
PU	8	EENT split drape, adh. split 9x61 cm	196x279 cm	1
SU	0	Back table cover,	152x190 cm	1




EENT pack

Item no.	DYJPEENTPEC	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, SMS reinf. 55x88 cm	80x142 cm	1
Material	Eclipse - Spunl	Head turban drape with adhesive	102x70 cm	1
PU	8	Hand towel	34x36 cm	2
SU	0	EENT split drape, adh. split 9x61 cm	198x279 cm	1
		Back table cover,	152x190 cm	1



EENT pack III

Item no.	DYNJPE7020	Description	Dimensions	Quantity
eClass-Nr.		Head turban drape with adhesive	102x70 cm	1
Material	Eclipse - Spunl	Suture bag, adh.	16x31 cm	1
PU	9	Mayo stand cover, SMS reinf. 46x61 cm	61x135 cm	1
SU	0	EENT split drape, adh. split 9x61 cm	198x279 cm	1
		Back table cover	125x229 cm	1




Drapes & Packs

ENT Surgery


ENT drape

Item no.	29670CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-05	Ent drape, reinforced, fenestration diameter 6cm, adhesive	150x196 cm	1
Material	ECL - Spunlac	CSR wrap	76x76 cm	1
PU	34			
SU	0			



ENT drape

Item no.	TB30456CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-05	Fenestrated drape, fenestration diameter 12cm, adhesive	220x240 cm	1
Material	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU	18			
SU	0			



ENT Drape

Item no.	ES30454CE	Description	Dimensions	Quantity
eClass-Nr.	ESS - 2 layers	ENT drape, 11 cm circular adhesive fenestration	175x250 cm	1
Material		SMS wrap		1
PU	32			
SU	555			



Drapes & Packs

ENT Surgery


ENT drape

Item no.	TB30665CE	Description	Dimensions	Quantity
eClass-Nr.		ENT drape, 12cm adhesive fenestration	158x220cm	1
Material	Ultimate - 3 la	SMS wrap		1
PU	22			
SU	555			




ENT Drape

Item no.	ES30458CE	Description	Dimensions	Quantity
eClass-Nr.		ENT drape, 7 cm circular adhesive fenestration	200x240cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	24			
SU	555			



ENT Drape

Item no.	ES10453CE	Description	Dimensions	Quantity
eClass-Nr.		Bi-laminate drape, 220x240 cm with 12 cm circular adhesive fenestration	220x240 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	24			
SU	555			




Drapes & Packs

ENT Surgery


ENT pack

Item no.	TB9192CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-05	Mayo stand cover, reinforced	80x142 cm	1
Material	Ultimate - 3 la	Head bar drape, adhesive	86x107 cm	1
PU	10	Split sheet, reinforced, split 6x64cm	196x279 cm	1
SU	0	Back table cover	140x190 cm	1




ENT Pack

Item no.	29232CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-06	Ent drape, reinforced, fenestration diameter 6cm, adhesive	150x196 cm	1
Material	ECL - Spunlac	Velcro strip	2.5x13 cm	1
PU	16	Absorbent Hand towels	20x30 cm	2
SU	0	Back table cover	140x190 cm	1
		CSR wrap	76x76 cm	1



ENT pack

Item no.	TB39450CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-05	Adhesive w/ round hole 5 cm	13 x 16,5 cm	1
Material	Ultimate - 3 la	Drape EENT/Opht., fen 12 x 7 cm oval ecc adh	270x196 cm	1
PU	12	Back table cover	140x190 cm	1
SU	0	CSR wrap	76x76 cm	1




Drapes & Packs

ENT Surgery

ENT pack

Item no.	CP29547CE	Description	Dimensions	Quantity
eClass-Nr.		Split Sheet	196x279 cm	1
Material	ECL - Spunlac	Head/Bar Drape, adhesive	82.5x107 cm	1
PU	12	CSR wrap	86x86 cm	1
SU	0			




ENT Pack

Item no.	CP29193CE	Description	Dimensions	Quantity
eClass-Nr.		Absorbent Towel	39x58 cm	3
Material	Ultimate - 3 la	OR gown, Essential, XL	130 cm	2
PU	4	OR gown, Essential, L	120 cm	1
SU	0	Head bar drape, adhesive	86x107 cm	1
		Split sheet, reinforced, split 6x64cm	196x279 cm	1
		Back table cover	140x190 cm	1



ENT Pack

Item no.	CP29023CE	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
Material	Ultimate - 3 la	Head bar drape, adhesive	86x107 cm	1
PU	17	Velcro strip	2.5x13 cm	1
SU	0	Three quarter drape	145x193 cm	1




Drapes & Packs

ENT Surgery


ENT Pack

Item no.	CP29036CE	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
Material	Ultimate - 3 la	Velcro strip	2.5x13 cm	1
PU	9	Side drape	75x98 cm	1
SU	0	Split sheet, reinforced, split 6x94cm, adhesive	196 x 275 cm	1
		Back table cover	140x190 cm	1




ENT Pack

Item no.	CP30167CE	Description	Dimensions	Quantity
eClass-Nr.		Adhesive tape, PE	8x40 cm	2
Material	ECL - Spunlac	Hand towels	30x39 cm	2
PU	0	Drape, w/ adhesive	150x196 cm	1
SU	0	Side drape, reinforced, adhesive	98x110 cm	3
		Back table cover	140x190 cm	1
		CSR wrap	86x86 cm	1



ENT pack

Item no.	9192CEA	Description	Dimensions	Quantity
eClass-Nr.		Head/Bar Drape, adhesive	82.5x107 cm	1
Material	Eclipse - Spunl	Split Sheet	196x279 cm	1
PU	8	Back table cover	140x190 cm	1
SU	0	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1




Drapes & Packs

ENT Surgery


ENT Pack

Item no.	TB9193CE	Description	Dimensions	Quantity
eClass-Nr.		Head bar drape, adhesive	86x107 cm	1
Material	Ultimate - 3 la	Split sheet, reinforced, split 6x64cm	196x279 cm	1
PU	9	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
SU	0	CSR wrap	86x86 cm	1
		Back table cover	140x190 cm	1




ENT Pack

Item no.	9193CEA	Description	Dimensions	Quantity
eClass-Nr.		Adhesive tape	9x50 cm	1
Material	Eclipse - Spunl	Drape	145x193 cm	1
PU	15	2 drape sheets	112x147 cm + 11	1
SU	0	CSR wrap	76x76 cm	1



ENT Pack

Item no.	ES10913CE	Description	Dimensions	Quantity
eClass-Nr.		Back table cover	140x190 cm	1
Material	ESS - 2 layers	Absorbent Hand towels	20x30 cm	2
PU	20	SMS wrap		1
SU	555	ENT drape, adhesive fenestration 13x17 cm	150x175 cm	1




Drapes & Packs

ENT Surgery


ENT pack

Item no.	ES10912CEB	Description	Dimensions	Quantity
eClass-Nr.		Absorbent Hand towels	20x30 cm	2
Material	ADV - SMMMS	Back table cover, reinforced, 50µm	140x190 cm	1
PU	18	CSR wrap	60x60 cm	1
SU	0	ENT drape, oval fenestration 8x10cm, adhesive	150x200 cm	1




ENT Pack

Item no.	ES9193CE	Description	Dimensions	Quantity
eClass-Nr.		Absorbent Hand towels	20x30 cm	2
Material	ESS - 2 layers	Head turban drape, adhesive	85x110 cm	1
PU	10	Split sheet, reinforced, split 6x64cm	198x280 cm	1
SU	555	SMS wrap		1
		Back table cover, reinforced, 50µm	140x190 cm	1
		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1



ENT Pack

Item no.	ES30460CE	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	2
Material	ESS - 2 layers	Adhesive tape, PE	8x40 cm	2
PU	10	ENT drape, adhesive fenestration 12x16 cm	228x300 cm	1
SU	0	Back table cover	140x190 cm	1




Drapes & Packs

ENT Surgery

Head bar drape

Item no.	TB9556CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-05	Head bar drape, adhesive	86x107 cm	1
Material	Ultimate - 3 la	CSR wrap	60x60 cm	1
PU	20			
SU	0			




Head bar drape

Item no.	DYNJPE7008	Description	Dimensions	Quantity
eClass-Nr.		Wrap	61x61 cm	1
Material	ECL - Spunlac	Head bar drape, w/ cuff 15 cm, adh.	102x69 cm	1
PU	20			
SU	0			



Head bar drape

Item no.	9556CEA	Description	Dimensions	Quantity
eClass-Nr.		Head/Bar Drape, adhesive	82.5x107 cm	1
Material	Eclipse - Spunl	CSR wrap	60x60 cm	1
PU	20			
SU	0			



Drapes & Packs

ENT Surgery


Jaw drape

Item no.	TB29004CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-05	Jaw drape, 15x11cm fenestration	95x145 cm	1
Material	Ultimate - 3 la			
PU	36			
SU	0			




Minor procedure drape

Item no.	ES10459CE	Description	Dimensions	Quantity
eClass-Nr.		Bi-laminate drape, 220x240 cm with 12 cm circular adhesive fenestration	220x240 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	24			
SU	555			



Nose surgery pack

Item no.	TB30C172CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-05	Adhesive tape, PE	8x40 cm	1
Material	Ultimate - 3 la	Absorbent Hand towels	20x30 cm	1
PU	12	Ent drape, fenestration diameter 12cm, adhesive	220x240 cm	1
SU	0	Back table cover	140x190 cm	1




Drapes & Packs

ENT Surgery


Pack small interventions

Item no.	TB30013CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Adhesive tape	9x50 cm	1
Material	Ultimate - 3 la	Absorbent Hand towels	20x30 cm	2
PU	12	Fenestrated drape, fenestration diameter 12cm, adhesive	220x240 cm	1
SU	0	Back table cover	140x190 cm	1
		CSR wrap	76x76 cm	1




Plastic surgery pack

Item no.	CP39250CE	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Hand towels</td> <td>30x39 cm</td> <td>2</td>		Hand towels	30x39 cm	2
Material	Eclipse - Spunl	CSR wrap	76x76 cm	1
PU	14	sponge swabs, gauze, woven, 12 ply, 17 threads	10x10 cm	16
SU	0	DRAPE EENT, FEN 13X17CM	145 x 145cm	1
		Back table cover, 50µm	120x140 cm	1



Set for small procedures

Item no.	ES10113CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-03-02	Absorbent Hand towels	20x30 cm	1
Material	ESS - 2 layers	CSR wrap	60x60 cm	1
PU	16	Fenestrated drape, adhesive oval fenestration (14*7cm)	240x200 cm	1
SU	0	Back table cover	140x190 cm	1




Drapes & Packs

ENT Surgery


Stomatology pack

Item no.	CP30002CE	Description	Dimensions	Quantity
eClass-Nr.		Drape, 220x240cm, 15x5 cm adhesive oval fenestration	220x240 cm	1
Material	Ultimate - 3 la	Back table cover	140x190 cm	1
PU	14	CSR wrap	86x86 cm	1
SU	0			




Surgery Pack, Sterile wrap

Item no.	TB30248CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-05	Drape sheet w/ adhesive fen. diam.10cm	98x150 cm	1
Material	Ultimate - 3 la	Back table cover	140x190 cm	1
PU	30	CSR wrap	76x76 cm	1
SU	0			



V shaped Head turban drape

Item no.	CP3954CE	Description	Dimensions	Quantity
eClass-Nr.		Adh extendable, turquoise	9x50 cm	1
Material	ECL - Spunlac	EENT drape sheet, 145x170cm		1
PU	28	CSR wrap	76x76 cm	1
SU	0			



Drapes & Packs

General Surgery

Absorbent drape with adhesive

Item no.	ABS5035CE	Description	Dimensions	Quantity
eClass-Nr.		absorbent drape, with adhesive	50x35cm	1
Material	Polyester and			
PU	88			
SU	0			



Leggings

Item no.	ES8420CE	Description	Dimensions	Quantity
eClass-Nr.		Leggings	75x120 cm	2
Material	ESS - 2 layers	SMS wrap		1
PU	36			
SU	0			




Drapes & Packs

General Surgery

► Lap Drapes and Packs


Abdominal pack

Item no.	DYJPEABPEC	Description	Dimensions	Quantity
eClass-Nr.		Laparoscopic abdominal drape, adh. fen. 25x28 cm, 2 pouches 25x76 cm	198/260x310 cm	1
Material	Eclipse - Spunl	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
PU	6	Back table cover,	152x190 cm	1
SU	0	Hand towel	34x36 cm	4




Endoscopy Drape, fen 11" x 14"

Item no.	TB9458CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Endoscopy Drape, 28x36 cm adhesive fenestration	196/254x314 cm	1
Material	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU	10			
SU	0			



Laparoscopic cholecystectomy drape

Item no.	DYNJPE3102	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	Eclipse - Spunl	Laparoscopic cholecystectomy drape, adh. fen. 25x28 cm, 2 trough..	198/259x310 cm	1
PU	8			
SU	0			




Drapes & Packs

General Surgery

► Lap Drapes and Packs


Laparoscopy drape,

Item no.	TB9429CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Drape, leggings, reinf, fen 25x28cm, adh	196/269x304 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	7			
SU	0			




Laparoscopy Pack

Item no.	TB9432CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-07	Absorbent Towel	39x58 cm	2
Material	Ultimate - 3 la	Drape, leggings, reinf, fen 25x28cm, adh	196/269x304 cm	1
PU	7	CSR wrap	86x86 cm	1
SU	0	Back table cover	140x190 cm	1



Laparoscopy Pack

Item no.	DYJPELAPPSM	Description	Dimensions	Quantity
eClass-Nr.	ADV - SMMM5	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
Material	8	Laparoscopic cholecystectomy drape, adh. fen. 30x33 cm	195/254x292 cm	1
PU	0	Back table cover,	152x190 cm	1
SU	0			



Drapes & Packs

General Surgery

► Lap Drapes and Packs


Laparoscopy pack

Item no.	ES10499CE	Description	Dimensions	Quantity
eClass-Nr.		Absorbent Hand towels	20x30 cm	4
Material	ESS - 2 layers	Video camera cover with pre-cut tip, with cardboard tab for applicat..	13x240 cm	1
PU	4	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
SU	555	Back table cover, reinforced, 50µm	140x190 cm	1
		Laparoscopy drape, abdominal fenestration 30x35cm	255x301 cm	1




Laparoscopy pack,

Item no.	TB9427CE	Description	Dimensions	Quantity
eClass-Nr. <td>34-35-01-07</td> <td>Drape, leggings, reinf, fen 25x28cm, adh</td> <td>196/269x304 cm</td> <td>1</td>	34-35-01-07	Drape, leggings, reinf, fen 25x28cm, adh	196/269x304 cm	1
Material <td>Ultimate - 3 la</td> <td>Mayo stand cover, absorbent reinf. 55x88 cm</td> <td>80x142 cm</td> <td>1</td>	Ultimate - 3 la	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
PU <td>6</td> <td>Back table cover</td> <td>140x190 cm</td> <td>1</td>	6	Back table cover	140x190 cm	1
SU <td>0</td> <td></td> <td></td> <td></td>	0			



Laparotomy drape

Item no.	TB29410CE	Description	Dimensions	Quantity
eClass-Nr. <td>34-35-01-01</td> <td>Drape with armboard covers, reinforced, fenestration 12x32cm, adhesive</td> <td>196/254x310 cm</td> <td>1</td>	34-35-01-01	Drape with armboard covers, reinforced, fenestration 12x32cm, adhesive	196/254x310 cm	1
Material <td>Ultimate - 3 la</td> <td>CSR wrap</td> <td>76x76 cm</td> <td>1</td>	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU <td>12</td> <td></td> <td></td> <td></td>	12			
SU <td>0</td> <td></td> <td></td> <td></td>	0			




Drapes & Packs

General Surgery

► Lap Drapes and Packs


Laparotomy drape

Item no.	DYNJPE3002	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	Eclipse - Spunl	Laparotomy drape, fen. w/ incise film 23x36 cm, instrument pad	198/259x307 cm	1
PU	14			
SU	0			




Laparotomy pack

Item no.	TB29138CE	Description	Dimensions	Quantity
eClass-Nr.		Drape with armboard covers, reinforced, fenestration 12x32cm, adhesive	196/254x310 cm	1
Material	Ultimate - 3 la	Suture bag, adh.	16x31 cm	1
PU	8	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
SU	0	Back table cover, reinforced, 50µm	140x190 cm	1
		Hand towels	30x20 cm	1



Laparotomy Sheet with Incise Film

Item no.	29511CEA	Description	Dimensions	Quantity
eClass-Nr.		Laparotomy drape, adh. fen. 23x37 cm	196/254x314 cm	1
Material	Eclipse - Spunl	CSR wrap	86x86 cm	1
PU	10			
SU	0			




Drapes & Packs

General Surgery

► Lap Drapes and Packs

Major abdominal drape

Item no.	DYNJPE3103	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	ECL - Spunlac	Laparoscopic abdominal drape, adh. fen. 25x28 cm, 2 pouches 25x76 cm	198/260x310 cm	1
PU	8			
SU	0			




Drapes & Packs

General Surgery

► Lap/Gyn Drapes and Packs


Abdominal drape

Item no.	DYJPEABDEC	Description	Dimensions	Quantity
eClass-Nr.		Wrap	100x100 cm	1
Material	ECL - Spunlac	Abdominal drape, abdominal fen. 25x30 cm, perineal fen. 12x24 cm..	240/250x247 cm	1
PU	6			
SU	0			




Coelioscopy pack

Item no.	29421CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-07	Coelioscopy drape, fenestration 35x27cm, adhesive	196x240 cm	1
Material	Eclipse - Spunl	Velcro strip	2.5x13 cm	2
PU	10	Hand towels	30x20 cm	2
SU	0	Leggings	140x79cm	2
		Back table cover	140x190 cm	1
		CSR wrap	86x86 cm	1



Laparoscopic assisted vaginal hysterectomy drape,

Item no.	TB29474CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-08	drape, reinf pouch adh abdo fen 18x23x23cm, fen with flap 15x8cm..	254/284x277 cm	1
Material	Ultimate - 3 la	CSR wrap	112x112 cm	1
PU	5			
SU	0			




Drapes & Packs

General Surgery

► Lap/Gyn Drapes and Packs


Laparoscopy drape

Item no.	DYNJPE9130	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	ECL - Spunlac	Leggings w/ cuff 15 cm	74x109 cm	2
PU	10	Laparoscopy pack, adh. triangular abdo. fen. 22x20x20 cm, adh. re..	89/259x246 cm	1
SU	0			




Laparoscopy Drape with leggings, 2 adhesive fenestrations

Item no.	TB29496CE	Description	Dimensions	Quantity
eClass-Nr. <td>34-35-01-01</td> <td>Absorbent Towel</td> <td>39x58 cm</td> <td>2</td>	34-35-01-01	Absorbent Towel	39x58 cm	2
Material <td>Ultimate - 3 la</td> <td>Laparoscopic drape, XL leggings, 34x38cm and 11x14cm adh. fen...</td> <td>255x310 cm</td> <td>1</td>	Ultimate - 3 la	Laparoscopic drape, XL leggings, 34x38cm and 11x14cm adh. fen...	255x310 cm	1
PU <td>5</td> <td>Back table cover</td> <td>140x190 cm</td> <td>1</td>	5	Back table cover	140x190 cm	1
SU <td>0</td> <td>CSR wrap</td> <td>112x112 cm</td> <td>1</td>	0	CSR wrap	112x112 cm	1



Laparoscopy Drape with leggings, 2 adhesive fenestrations

Item no.	TB27052CEA	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Under buttocks drape with cuff 24cm, reinforced</td> <td>96x126 cm</td> <td>1</td>		Under buttocks drape with cuff 24cm, reinforced	96x126 cm	1
Material <td>Ultimate - 3 la</td> <td>Drape adhesive</td> <td>75x75 cm</td> <td>1</td>	Ultimate - 3 la	Drape adhesive	75x75 cm	1
PU <td>4</td> <td>Adhesive tape</td> <td>9x50 cm</td> <td>1</td>	4	Adhesive tape	9x50 cm	1
SU <td>555</td> <td>CSR wrap</td> <td>86x86 cm</td> <td>1</td>	555	CSR wrap	86x86 cm	1
		Back table cover, folded, 50µm	140x190 cm	1
		Laparoscopic drape, XL leggings, 34x38cm and 11x14cm adh. fen...	255x310 cm	1
		Back table cover	140x190 cm	1




Drapes & Packs

General Surgery

► Lap/Gyn Drapes and Packs


Laparoscopy pack

Item no.	TB29142CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-07	Leggings	140x79 cm	2
Material	Ultimate - 3 la	Hand towels	30x20 cm	1
PU	9	Back table cover, reinforced, 50µm	125x230 cm	1
SU	0	T drape, reinf, abdo trapezoid fen. adh, perineal fen 8x15cm	248/98 x 246 cm	1




Laparoscopy Pack

Item no.	CP29047CEA	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>instrument pouch, 3 compartments, adhesive</td> <td>25x48 cm</td> <td>1</td>		instrument pouch, 3 compartments, adhesive	25x48 cm	1
Material	Ultimate - 3 la	back table cover, 90 µm	203x240 cm	1
PU	4	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
SU	0	Absorbent Towel	39x58 cm	4
		Video camera cover telescopic folded	13 x 235 cm	1
		Velcro strip	2.5x13 cm	1
		Laparoscopic drape, XL leggings, 34x38cm and 11x14cm adh. fen...	255x310 cm	1
		CSR wrap	112x112 cm	1



Laparoscopy pack

Item no.	ES10496CE	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Absorbent Towel</td> <td>39x58 cm</td> <td>2</td>		Absorbent Towel	39x58 cm	2
Material	ESS - 2 layers	Adhesive tapes	9x50 cm	1
PU	7	Instrument pouch adhesive, 2 compartments	61x42.5 cm	1
SU	0	SMS wrap		1
		Back table cover	140x190 cm	1
		Laparoscopy drape	255x306 cm	1




Drapes & Packs

General Surgery

► Lap/Gyn Drapes and Packs

Laparoscopy pack


Item no. ES29496CE
eClass-Nr.
Material ESS - 2 layers
PU 6
SU 555



Description	Dimensions	Quantity
Absorbent Hand towels	20x30 cm	2
Adhesive tapes	9x50 cm	1
Side sheet, adhesive	75x90 cm	1
SMS wrap		1
Back table cover, reinforced, 50µm	140x190 cm	1
Laparoscopy drape	255x301 cm	1

Laparoscopy pack I


Item no. DYNJPE9100SM
eClass-Nr.
Material ADV - SMMMS
PU 10
SU 0



Description	Dimensions	Quantity
Adhesive drape	38x66 cm	1
Leggings w/ cuff 15 cm	74x109 cm	2
Suture bag, adh.	16x31 cm	1
Back table cover	125x229 cm	1
Hand towel	51x43 cm	1
Laparoscopy drape, adh. abdo. fen. 19x20cm, peri. fen. 8x15cm	89/259x246 cm	1

LAVH drape

Item no. DYNJPE9103
eClass-Nr.
Material Eclipse - Spunl
PU 8
SU 0



Description	Dimensions	Quantity
Wrap	100x100 cm	1
LAVH drape, adh. trapezoidal fen. 17/23x23 cm, adh. perineal fen. ...	254/277x274 cm	1


Drapes & Packs

General Surgery

► Lap/Gyn Drapes and Packs


LAVH pack

Item no.	Material	Description	Dimensions	Quantity
Item no.	CPU29063CE			
eClass-Nr.	Ultimate - 3 la	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
Material	4	instrument pouch, 3 compartments, adhesive	25x48 cm	1
PU	0	Absorbent Towel	39x58 cm	4
SU		Adhesive tape	9x50 cm	2
		Video camera cover telescopic folded	13 x 235 cm	1
		Under buttocks drape with cuff 24cm, reinforced	96x126 cm	1
		Velcro strip	2.5x13 cm	1
		drape, reinf pouch adh abdo fen 18x23x23cm, fen with flap 15x8cm..	254/284x277	1
		Back table cover, 75µm	160x240 cm	1
		CSR wrap	60x60 cm	1



LAVH pack with gowns

Item no.	Material	Description	Dimensions	Quantity
Item no.	CP29121CE			
eClass-Nr.	Ultimate - 3 la	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
Material	2	Adhesive tape	9x50 cm	2
PU	0	Video camera cover telescopic folded	13 x 235 cm	1
SU		Under buttocks drape with cuff 24cm, reinforced	96x126 cm	1
		Velcro strip	2.5x13 cm	1
		instrument pouch, 3 compartments, adhesive	25x48 cm	1
		Absorbent Towel	39x58 cm	4
		OR gown, Eclipse, XL	130 cm	2
		drape, reinf pouch adh abdo fen 18x23x23cm, fen with flap 15x8cm..	254/284x277	1
		Back table cover, 75µm	160x240 cm	1
		CSR wrap	60x60 cm	1
		Gown Eclipse, L size	120 cm	1



Rectum Pack

Item no.	Material	Description	Dimensions	Quantity
Item no.	TB27050CE			
eClass-Nr.	34-35-01-01	Absorbent Hand towels	20x30 cm	1
Material	Ultimate - 3 la	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
PU	4	Under buttocks drape with cuff 24cm, reinforced	96x126 cm	1
SU	0	Back table cover	140x190 cm	1
		Laparoscopic drape, XL leggings, 34x38cm and 11x14cm adh. fen...	255x310 cm	1
		Adh extendable, turquoise	9x50 cm	3




Drapes & Packs

General Surgery

► Lap/Gyn Drapes and Packs


Rectum pack

Item no.	DYJPEREPSM	Description	Dimensions	Quantity
eClass-Nr.		Rectum drape, adh. abdo. fen. 29x29cm, adh. rectum fen. 16x7cm	258/250x337 cm	1
Material	ADV - SMMMS	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
PU	6			
SU	0	Back table cover,	152x190 cm	1



Robotics drape

Item no.	DYNJPE9106SM	Description	Dimensions	Quantity
eClass-Nr.		Robotics drape, adb. fen. 33x38cm, per. fen. 8x15cm	254/277x310 cm	1
Material	ADV - SMMMS	SMS wrap		1
PU	5			
SU	0			




Drapes & Packs

General Surgery

► Thyroid Drapes and Packs


Thyroid drape

Item no.	TB29522CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Thyroid drape 196/269x309cm, w/diamond shaped fen. 13x13cm	196/269X309 cm	1
Material	Ultimate - 3 la	Back table cover	140x190 cm	1
PU	10			
SU	0			



Thyroid drape

Item no.	DYNJPE7002	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	Eclipse - Spunl	Thyroid drape, adh. fen. 13x18 cm	196x307 cm	1
PU	10			
SU	0			



Thyroid drape,

Item no.	TB29422CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Thyroid drape, reinforced, diamond fenestration 13x13cm, adhesive	196x311 cm	1
Material	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU	12			
SU	0			




Drapes & Packs

General Surgery

► Thyroid Drapes and Packs


Thyroid pack

Item no.	DYJPETHYPEC2	Description	Dimensions	Quantity
eClass-Nr.		Foot split sheet, adh. split 20x25 cm	180x240 cm	1
Material	ECL - Spunlac	Split drape, adh. split 20x72 cm	150x240 cm	1
PU	5	Back table cover,	152x190 cm	1
SU	0	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1



Thyroid pack

Item no.	TB277041CE	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
Material	Ultimate - 3 la	Foot split sheet, adhesive	183x240 cm	1
PU	6	Head split sheet, adhesive	240x145 cm	1
SU	0	Back table cover, reinforced, 50µm	140x190 cm	1




Drapes & Packs

General Surgery

► Universal Drapes and Packs


2 Side drapes, reinforced - 98x112cm

Item no.	TB29400CE	Description	Dimensions	Quantity
eClass-Nr.		Side drape, reinforced, 98 x 112 cm	98 x 112 cm	2
Material	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU	24			
SU	0			




Basic drape pack

Item no.	ES30015CE	Description	Dimensions	Quantity
eClass-Nr.		Absorbent Hand towels	20x30 cm	2
Material	ESS - 2 layers	side drape, adh	90x100 cm	4
PU	22	SMS wrap		1
SU	555	Back table cover	140x190 cm	1



Basic pack

Item no.	TB30014CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Side drape	75x98 cm	4
Material	Ultimate - 3 la	Absorbent Hand towels	20x30 cm	2
PU	16	Back table cover	140x190 cm	1
SU	0	CSR wrap	76x76 cm	1




Drapes & Packs

General Surgery

► Universal Drapes and Packs


Basic pack

Item no.	ES15307CE	Description	Dimensions	Quantity
eClass-Nr.		drape sheet	75x90cm	1
Material	ESS - 2 layers	Drape adhesive	75x75 cm	4
PU	16			
SU	0			




Minor surgery pack

Item no.	ES15303CE	Description	Dimensions	Quantity
eClass-Nr.		drape sheet	75x90cm	1
Material	ESS - 2 layers	Drape, adhesive	45x50 cm	4
PU	40			
SU	0			



Pack for small Procedures

Item no.	TB30198CE	Description	Dimensions	Quantity
eClass-Nr.		Back table cover, folded, 75µm	140x120 cm	1
Material	Ultimate - 3 la	Back table cover, 50µm	120x140 cm	1
PU	12	Drape	47x49 cm	1
SU	0	Drape, w/fen, w/adh	240x196 cm	1




Drapes & Packs

General Surgery

► Universal Drapes and Packs

Pack for small procedures with gowns

Item no.	CP29151CE	Description	Dimensions	Quantity
eClass-Nr.		Hand towels	30x20 cm	2
Material	ECL - Spunlac	OR gown, Essential, XL	130 cm	2
PU	12	Small drape, with adhesive fenestration 5cm	50x60cm	1
SU	0	Back table cover	140x190 cm	1




Pubic drape

Item no.	ES3007CE	Description	Dimensions	Quantity
eClass-Nr.		Perineal sheet, adhesive	16x37 cm	1
Material	ESS - 2 layers			
PU	200			
SU	0			



Sterile barrier pack

Item no.	TB29547CE	Description	Dimensions	Quantity
eClass-Nr.		Tri-laminate panel, 200x187 cm, 4 loops	200x187 cm	1
Material	Ultimate - 3 la	Tri-laminate panel, 200x187 cm, adhesive, 4 loops	200x187 cm	1
PU	10	95x145 cm trilaminate panel with adhesive	95x145 cm	1
SU	0	CSR wrap	76x76 cm	1




Drapes & Packs

General Surgery

► Universal Drapes and Packs


Universal pack

Item no.	TB29105CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Velcro strip	2.5x13 cm	2
Material	Ultimate - 3 la	Absorbent Hand towels	20x30 cm	4
PU	6	Side sheets, reinforced, adhesive	98x120 cm	2
SU	0	Foot sheet, reinforced, adhesive	196x210 cm	1
		Head sheet, reinforced, adhesive	160x260 cm	1
		Back table cover	140x190 cm	1
		CSR wrap	112x112 cm	1




Universal pack

Item no.	TB30105CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Side drape	75x98 cm	2
Material	Ultimate - 3 la	Velcro strip	2.5x13 cm	1
PU	10	Absorbent Hand towels	20x30 cm	2
SU	0	Foot sheet, adhesive	196x200 cm	1
		Head drape, adhesive	260x160 cm	1
		CSR wrap	86x86 cm	1



Universal pack

Item no.	TB30114CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Side drape	75x98 cm	2
Material	Ultimate - 3 la	Velcro strip	2.5x13 cm	1
PU	8	Absorbent Hand towels	20x30 cm	4
SU	0	Mayo stand cover, reinforced	80x142 cm	1
		Foot sheet, adhesive	196x200 cm	1
		Head drape, adhesive	260x160 cm	1
		Back table cover	140x190 cm	1




Drapes & Packs

General Surgery

► Universal Drapes and Packs


Universal pack

Item no.	TB30110CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Side drape	75x98 cm	2
Material	Ultimate - 3 la	Velcro strip	2.5x13 cm	1
PU	8	Absorbent Hand towels	20x30 cm	4
SU	0	Foot sheet, adhesive	196x200 cm	1
		Head drape, adhesive	260x160 cm	1
		Back table cover	140x190 cm	1
		CSR wrap	86x86 cm	1




Universal Pack

Item no.	29114CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Head sheet, reinforced, adhesive	183x254 cm	1
Material	ECL - Spunlac	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
PU	6	Side Sheet	98x112 cm	2
SU	0	Foot sheet w/ adh	183 x 191 cm	1
		Back table cover	140x190 cm	1



Universal Pack

Item no.	29906CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Head sheet, reinforced, adhesive	183x254 cm	1
Material	ECL - Spunlac	Side Sheet	98x112 cm	1
PU	4	Drape, w/adh, w/reinf	300x183 cm	1
SU	0	Adhesive tape, PE	8x40 cm	2
		Back table cover, 90µm	160x240 cm	1
		Back table cover, folded, 50µm	140x190 cm	1
		Hand towels	30x20 cm	2
		Instrument pouch in clear plastic with adhesive	33x41cm	1
		CSR wrap	112x112 cm	1




Drapes & Packs

General Surgery

► Universal Drapes and Packs


Universal Pack

Item no.	29205CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Drape	98x150 cm	1
Material	ECL - Spunlac	Drape adhesive	75x75 cm	4
PU	6	Adhesive tape, PE	8x40 cm	2
SU	0	Back table cover, folded, 50µm	140x120 cm	1
		Absorbent Towel	39x58 cm	2
		Side drape, reinforced, adhesive	110x98 cm	2
		Head sheet w/ reinf	196x240 cm	2
		Back table cover	140x190 cm	1
		CSR wrap	112x112 cm	1




Universal Pack

Item no.	CP29136CE	Description	Dimensions	Quantity
eClass-Nr. <td>Ultimate - 3 Ia</td> <td>Pouch, Instmt, double, PE</td> <td>33x41 cm</td> <td>1</td>	Ultimate - 3 Ia	Pouch, Instmt, double, PE	33x41 cm	1
Material		Hand towels	30x20 cm	4
PU	4	OR gown, Advanced, L	130 cm	1
SU	0	OR gown, Advanced, XL	140 cm	2
		Cotton drape absorbent, blue, Xray	43x66 cm	4
		Velcro strip	2.5x13 cm	1
		Side drape	75x98 cm	2
		Foot sheet, adhesive	196x200 cm	1
		Head drape, adhesive	260x160 cm	1
		Back table cover	140x190 cm	1



Universal pack

Item no.	DYJPEUNPSMG	Description	Dimensions	Quantity
eClass-Nr. <td>34-35-01-01</td> <td>Foot drape, adh.</td> <td>183x175 cm</td> <td>1</td>	34-35-01-01	Foot drape, adh.	183x175 cm	1
Material	ADV - SMMMS	Side drapes, adh.	75x93 cm	2
PU	4	Surgical gown, XL		1
SU	0	Surgical gown, L		1
		Hand towel	34x36 cm	4
		Back table cover,	152x190 cm	1
		Head drape, adh.	260x157 cm	1
		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1




Drapes & Packs

General Surgery

► Universal Drapes and Packs


Universal pack

Item no.	DYJPEUNPSM2	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Head drape, adh.	152x250 cm	1
Material	ADV - SMMMS	Side drapes, adh.	75x93 cm	2
PU	5	Foot drape, adh.	183x175 cm	1
SU	0	Back table cover,	152x190 cm	1
		Hand towel	34x36 cm	4
		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1




Universal pack

Item no.	SM29105CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Side drape, adh.	102x112 cm	2
Material	ADV - SMMMS	Foot drape, adh.	196x196 cm	1
PU	5	Wrap	91x91 cm	1
SU	0	Back table cover,	152x190 cm	1
		Hand towel	34x36 cm	4
		Head drape, adh.	152x260 cm	1



Universal pack

Item no.	DYJPEUNPSM1	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Side drape, adh.	112x102 cm	2
Material	ADV - SMMMS	Foot drape, adh.	183x175 cm	1
PU	5	Head drape, adh.	260x157 cm	1
SU	0	Back table cover,	152x190 cm	1
		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1




Drapes & Packs

General Surgery

► Universal Drapes and Packs


Universal pack

Item no.	DYJPEUNMPM1	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Side drape, adh.	112x102 cm	2
Material	ADV - SMMMS	Foot drape, adh.	183x175 cm	1
PU	5	Head drape, adh.	300x183 cm	1
SU	0	Back table cover,	152x190 cm	1
		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1




Universal Pack

Item no.	CP29114CEA	Description	Dimensions	Quantity
eClass-Nr.		Pouch, Instmt, double, PE	33x41 cm	1
Material	ULT- 3 laggig	Drape sheet, adhesive	196x145 cm	1
PU	4	Absorbent Towel	39x58 cm	4
SU	0	Foot drape, adhesive, reinforced	190x183 cm	1
		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
		Velcro strip	2.5x13 cm	1
		Adhesive tape	9x50 cm	2
		Side sheets, reinforced, adhesive	98x120 cm	2
		Head sheet, reinforced, adhesive	254x183 cm	1
		back table cover, 90 µm	203x240 cm	1
		CSR wrap	60x60 cm	1



Universal pack

Item no.	DYJPEUNMPEC1	Description	Dimensions	Quantity
eClass-Nr.	ECL - Spunlac	Side drape, adh.	89x112 cm	2
Material		Foot drape, adh.	196x178 cm	1
PU	5	Head drape, adh.	300x178 cm	1
SU	0	Back table cover,	152x190 cm	1
		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1



Drapes & Packs

General Surgery

► Universal Drapes and Packs


Universal Pack

Item no.	CP29115CE	Description	Dimensions	Quantity
eClass-Nr.		Side Sheet	98x112 cm	2
Material	Eclipse - Spunl	Head sheet, reinforced, adhesive	183x254 cm	1
PU	3	Instrument pouch with adhesive, 2 compartments	43,5x49 cm	1
SU	0	Drape	75x90cm	1
		Foot sheet w/ adh	183 x 191 cm	1
		OR gown, Eclipse, XLL, fabric reinforced	150 cm	2
		OR gown, XL, 130 cm, fabric reinforced	130 cm	1
		CSR wrap	119x179 cm	2
		Absorbent Towel	39x58 cm	6
		Mayo stand cover, ultra absorbent reinforcement 75x100 cm	80x145 cm	1
		Adhesive tape	9x50 cm	2




Universal pack

Item no.	29823CEA	Description	Dimensions	Quantity
eClass-Nr.		Instrument pouch, PE, adhesive	30,5x29 cm	1
Material	ECL - Spunlac	Side Drape reinforced	98x110 cm	2
PU	3	Drape	98x98 cm	1
SU	0	spunlace foot drape, 183x190 cm, 37x67 cm reinforcement, 2 cord..	183 x 190 cm	1
		OR gown, Eclipse, XLL, fabric reinforced	150 cm	2
		OR gown, Eclipse, XL	130 cm	1
		Head drape	251x183 cm	1
		Mayo stand cover, ultra absorbent reinforcement	80x145 cm	1
		Absorbent Hand towels	20x30 cm	1
		Back table cover	140x190 cm	1



Universal pack

Item no.	ES29114CE	Description	Dimensions	Quantity
eClass-Nr.		bi-laminate side drape	110x95 cm	2
Material	ESS - 2 layers	bi-laminate head drape	260x170 cm	1
PU	8	foot drape	200x200 cm	1
SU	0	Hook and loop tube holder	14x2,5 cm	1
		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
		Back table cover,	152x190 cm	1
		Adhesive tape	9x50 cm	1
		Absorbent Hand towels	20x30 cm	4




Drapes & Packs

General Surgery

► Universal Drapes and Packs


Universal pack

Item no.	ES29105CE	Description	Dimensions	Quantity
eClass-Nr.		Hook and loop tube holder	14x2.5 cm	1
Material	ESS - 2 layers	Adhesive tape	9x50 cm	1
PU	8	Absorbent Hand towels	20x30 cm	2
SU	0	bi-laminate side drape	110x95 cm	2
		bi-laminate head drape	260x170 cm	1
		foot drape	200x200 cm	1
		SMS wrap		1
		Back table cover, reinforced, 50µm	140x190 cm	1




Universal pack

Item no.	ES29104CE	Description	Dimensions	Quantity
eClass-Nr.		Hook and loop tube holder	14x2.5 cm	1
Material	ESS - 2 layers	foot drape	175x190 cm	1
PU	8	bi-laminate side drape	90x75 cm	2
SU	555	Adhesive tape	9x50 cm	1
		bi-laminate head drape	240x150 cm	1
		Absorbent Hand towels	20x30 cm	4
		SMS wrap		1
		Back table cover, reinforced, 50µm	140x190 cm	1



Universal Pack Large

Item no.	ES10116CE	Description	Dimensions	Quantity
eClass-Nr.		Head drape, adhesive	160x260 cm	1
Material	ESS - 2 layers	bi-laminate side drape	110x100 cm	2
PU	10	Foot sheet, adhesive	190x190 cm	1
SU	555	Adhesive tape	9x50 cm	1
		Back table cover	140x190 cm	1
		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
		Absorbent Hand towels	20x30 cm	2




Drapes & Packs

General Surgery

► Universal Drapes and Packs


Universal pack large

Item no.	SM29107CE	Description	Dimensions	Quantity
eClass-Nr.		head drape	180x300 cm	1
Material	ADV - SMMMS	side drape	100x110 cm	2
PU	8	foot drape	240x190 cm	1
SU	0	Adhesive tape	9x50 cm	2
		SMS wrap		1
		Back table cover	140x190 cm	1
		Absorbent Hand towels	20x30 cm	2




Universal pack w/ mayo stand cover

Item no.	ES10114CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
Material	ESS - 2 layers	Hand towels	30x20 cm	2
PU	10	Adhesive tapes	9x50 cm	1
SU	0	Side sheet, adhesive	75x90 cm	2
		Foot drape, adhesive	175x175 cm	1
		Head sheet, adhesive	260x160 cm	1
		Back table cover	140x190 cm	1



Universal pack w/o mayo stand cover

Item no.	ES10110CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Foot drape, adhesive	175x175 cm	1
Material	ESS - 2 layers	Absorbent Hand towels	20x30 cm	2
PU	10	CSR wrap	86x86 cm	1
SU	0	Adhesive tapes	9x50 cm	1
		Side sheet, adhesive	75x90 cm	2
		Head sheet, adhesive	260x160 cm	1
		Back table cover	140x190 cm	1




Drapes & Packs

General Surgery

► Universal Drapes and Packs

Universal pack w/o mayo stand cover

Item no.	ES10111CE	Description	Dimensions	Quantity
eClass-Nr.		Side sheet, adhesive	98x75 cm	2
Material	ESS - 2 layers	Foot drape, adhesive	190x200 cm	1
PU	10	Hook and loop tube holder	14x2.5 cm	1
SU	555	Head drape, adhesive	160x260 cm	1
		SMS wrap		1
		Back table cover, reinforced, 50µm	140x190 cm	1
		Absorbent Hand towels	20x30 cm	4




Universal pack,

Item no.	TB29106CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Hook and loop tube holder	14x2.5 cm	1
Material	Ultimate - 3 la	Absorbent Hand towels	20x30 cm	4
PU	6	Side sheets, reinforced, adhesive	98x100 cm	2
SU	0	Foot sheet, reinforced, adhesive	196x200 cm	1
		Head sheet, reinforced, adhesive	160x260 cm	1
		Back table cover	140x190 cm	1

No picture available

Universal pack,

Item no.	TB29114CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Velcro strip	2.5x13 cm	1
Material	Ultimate - 3 la	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
PU	6	Absorbent Hand towels	20x30 cm	4
SU	0	Side sheets, reinforced, adhesive	98x100 cm	2
		Foot sheet, reinforced, adhesive	196x200 cm	1
		Head sheet, reinforced, adhesive	160x260 cm	1
		Back table cover	140x190 cm	1




Drapes & Packs

General Surgery

► Universal Drapes and Packs


Universal pack,

Item no.	ES10115CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Side sheet, adhesive	75x90 cm	2
Material	ESS - 2 layers	Foot drape, adhesive	175x175 cm	1
PU	12	Head sheet, adhesive	240x150 cm	1
SU	0	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
		Adhesive tapes	9x50 cm	1
		Absorbent Hand towels	20x30 cm	4
		Back table cover	140x190 cm	1




Universal pack,

Item no.	ES29115CE	Description	Dimensions	Quantity
eClass-Nr.		bi-laminate head drape	240x150 cm	1
Material	ESS - 2 layers	bi-laminate side drape	90x75 cm	2
PU	8	foot drape	175x190 cm	1
SU	555	Back table cover,	152x190 cm	1
		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
		Adhesive tapes	9x50 cm	1
		Absorbent Hand towels	20x30 cm	4



Universal pack,

Item no.	TB29007CE	Description	Dimensions	Quantity
eClass-Nr.		head drape, large	300 x 183cm	1
Material	Ultimate - 3 la	Adhesive tape	9x50 cm	1
PU	4	Side sheets, reinforced, adhesive	98x100 cm	3
SU	555	Foot drape, adhesive, reinforced	190x183 cm	1
		Instrument pouch in clear plastic with adhesive	33x51cm	1
		Back Table Cover, reinforced 50µ	110x100 cm	1
		Back table cover, folded, 50µm	140x190 cm	1
		Back table cover, reinforced, 50µm	140x190 cm	1
		Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1
		Half drape	112x145 cm	1




Drapes & Packs

General Surgery

► Universal Drapes and Packs


Universal Pack, reinforced

Item no.	CP29120CE	Description	Dimensions	Quantity
eClass-Nr.		Adhesive tapes	9x50 cm	2
Material	ECL - Spunlac	Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1
PU	2	Drape, w/adh, reinf heavy	196x340 cm	1
SU	0	spunlace foot drape, 183x190 cm, 37x67 cm reinforcement, 2 cord..	183 x 190 cm	1
		OR gown, Eclipse, XLL, fabric reinforced	150 cm	3
		OR gown, Essential, XL	130 cm	1
		Hand towels	30x39 cm	3
		Back table cover, 75µm	160x240 cm	1
		Side Sheet	98x112 cm	2




Universal Pack, reinforced

Item no.	ES29106CE	Description	Dimensions	Quantity
eClass-Nr.		head drape with armboard covers	269x182 cm	1
Material	ESS - 2 layers	bi-laminate side drape	110x95 cm	2
PU	8	foot drape	200x200 cm	1
SU	0	Back table cover, reinforced, 50µm	140x190 cm	1
		Absorbent Hand towels	20x30 cm	2
		Adhesive tape	9x50 cm	1



Universal Pack, w/ gown

Item no.	TB29115CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	head drape, large	300 x 183cm	1
Material	ULT- 3 laggig	Mayo stand cover, reinforced	80x142 cm	1
PU	3	Foot drape, adhesive, reinforced	190x183 cm	1
SU	0	Side drape, reinforced, 98 x 112 cm	98 x 112 cm	2
		OR gown, Eclipse, XLL	150 cm	2
		OR gown, Eclipse, XL	130 cm	1
		Hand towels	30x39 cm	6
		Back table cover	140x190 cm	1




Drapes & Packs

General Surgery

► Universal Drapes and Packs

Universal set pack, Maxi

Item no.	TB29005CE	Description	Dimensions	Quantity
eClass-Nr.	Ultimate - 3 la	Side drape, reinforced, 98 x 112 cm	98 x 112 cm	2
Material	4	head drape, large	300 x 183cm	1
PU	0	Suture bag	25x33 cm	1
SU	0	Foot drape, adhesive, reinforced	190x183 cm	1
		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
		Back table cover, 75µm	160x240 cm	1
		Adhesive tapes	9x50 cm	1
		Absorbent Towel	39x58 cm	2




Drapes & Packs

Gynaecological Surgery

Under Buttocks Drape

Item no.	SPT-4576	Description	Dimensions	Quantity
eClass-Nr.		Under buttocks drape, 96x20 cm reinf, 2 pouches	96x120 cm	1
Material	PE			
PU	32			
SU	0			




Drapes & Packs

Gynaecological Surgery

► Baby Blanket


Baby blanket

Item no.	1001CEB	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-08	Drape blanket Baby, White	100x105 cm	1
Material	ECL - Spunlac	CSR wrap	60x60 cm	1
PU	50			
SU	0			




Baby blanket

Item no.	SPTE-1001	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-08	Drape blanket Baby, White	100x105 cm	1
Material	ECL - Spunlac			
PU	40			
SU	0			



Baby Blanket

Item no.	RG966B	Description	Dimensions	Quantity
eClass-Nr.		Paper Wrap	75 x 75 cm	1
Material	Spunlace	Baby Blanket, Spunlace	100 x 100 cm	1
PU	36			
SU	0			



Drapes & Packs

Gynaecological Surgery

► Obstetrical Drapes and Packs


Cæsarean Birth Drape with 270° Pouch

Item no.	TB29452CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-08	Cesarean birth drape, w/adh, w/FCP	196/249x300 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	8			
SU	0			




Cæsarean Birth Drape with 360° Pouch

Item no.	TB29448CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-08	Cesarean birth drape, w/adh, w/FCP	196/249x314 cm	1
Material	Ultimate - 3 la	CSR wrap	112x112 cm	1
PU	7			
SU	0			



Cæsarean Birth Drape with 360° Pouch

Item no.	ES10852CE	Description	Dimensions	Quantity
eClass-Nr.	ESS - 2 layers	SMS wrap		1
Material				
PU	8			
SU	0			




Drapes & Packs

Gynaecological Surgery

► Obstetrical Drapes and Packs


Cæsarean Birth Drape with Incise Film

Item no.	TB29446CE	Description	Dimensions	Quantity
eClass-Nr.		drape, armboard covers, reinf pouch, pear-shape fen 30x36cm soli..	196/249x314	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	10			
SU	0			




Cæsarean birth pack

Item no.	CP29045CE	Description	Dimensions	Quantity
eClass-Nr.		Absorbent Towel	39x58 cm	4
Material	ULT- 3 laggig	Drape blanket Baby, White	100x105 cm	1
PU	5	Back table cover, 75µm, folded	140x190cm	1
SU	0	Cesarean birth drape, w/adh, w/FCP	196/249x314 cm	1
		Back table cover	140x190 cm	1
		CSR wrap	60x60 cm	1



Cæsarean birth pack

Item no.	CP29071CE	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Absorbent Towel</td> <td>39x58 cm</td> <td>4</td>		Absorbent Towel	39x58 cm	4
Material	ECL - Spunlac	OR gown, Eclipse, L	120 cm	1
PU	4	OR gown, Eclipse, XL	130 cm	2
SU	0	Drape blanket Baby, White	100x105 cm	1
		drape, armboard covers, reinf pouch, pear-shape fen 30x36cm soli..	196/249x314	1
		Back table cover, 90µm	160x240 cm	1




Drapes & Packs

Gynaecological Surgery

► Obstetrical Drapes and Packs


Cæsarean birth pack

Item no.	ES29445CE	Description	Dimensions	Quantity
eClass-Nr.		Absorbent Hand towels	20x30 cm	4
Material	ESS - 2 layers	Instrument pouch in clear plastic with adhesive	33x51cm	1
PU	5	Drape blanket baby, white	100x105 cm	1
SU	555	C-Section drape with 360° fluid collection pouch	196/259x314 cm	1
		Back table cover	140x190 cm	1
		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1




Cesarean birth pack,

Item no.	TB29453CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-08	Cesarean birth drape, w/adh, w/FCP	196/249x300 cm	1
Material	Ultimate - 3 la	Drape blanket Baby, White	100x105 cm	1
PU	6	Mayo stand cover, reinforced	80x142 cm	1
SU	0	Absorbent Hand towels	20x30 cm	2
		Back table cover	140x190 cm	1



Cesarean pack birth, attached leggings,

Item no.	TB29454CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-08	Drape blanket Baby, White	100x105 cm	1
Material	Ultimate - 3 la	Mayo stand cover, reinforced	80x142 cm	1
PU	5	drape, armboard covers, leggings, reinf, pouch, shaped fen 30x36cm, adh	260/196x297 cm	50
SU	0	Back table cover	140x190 cm	1




Drapes & Packs

Gynaecological Surgery

► Obstetrical Drapes and Packs


Cesarean pack,

Item no.	TB29194CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-08	Adhesive tape	9x50 cm	1
Material	Ultimate - 3 la	Absorbent Hand towels	20x30 cm	2
PU	6	Back table cover, reinforced, 50µ, 140x120cm	140x120 cm	1
SU	0	C-section drape, armboard covers, reinf, pouch, shaped adh fen 30x36cm	249/196 x 300 cm	1
		Baby blanket	100 x 105 cm	1
		Back table cover	140x190 cm	1
		CSR wrap	86x86 cm	1




C-section drape

Item no.	DYJPECSEPEC2	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Back table cover,</td> <td>152x190 cm</td> <td>1</td>		Back table cover,	152x190 cm	1
Material	Eclipse - Spunl	Mayo stand cover, SMS reinf. 55x88 cm	80x142 cm	1
PU	6	C-section drape, 33x38 cm incise film w/ pear shaped fen., pouch	115/259x307 cm	1
SU	0	Baby blanket	76x86 cm	1



C-Section drape

Item no.	DYNJPE6102SM	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Wrap</td> <td>91x91 cm</td> <td>1</td>		Wrap	91x91 cm	1
Material	ADV - SMMMS	C-section drape, incise film w/ pear shaped fen. 18x18 cm	196/259x307 cm	1
PU	7			
SU	0			




Drapes & Packs

Gynaecological Surgery

► Obstetrical Drapes and Packs


C-Section Drape with fenestrated incise and pouch

Item no.	DYNJPE6102	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	ECL - Spunlac	C-section drape, incise film w/ pear shaped fen. 18x18 cm	198/259x307 cm	1
PU	7			
SU	0			




C-section drape with incise and pouch

Item no.	DYNJPE6101	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	ECL - Spunlac	C-section drape, fen. w/ incise film 15x15 cm, pouch	198/259x307 cm	1
PU	7			
SU	0			



C-Section drape with suction port

Item no.	DYNJPE6105SM	Description	Dimensions	Quantity
eClass-Nr.		SMS wrap		1
Material	ADV - SMMMS	C-Section drape with suction port	196/259x307 cm	1
PU	7			
SU	0			




Drapes & Packs

Gynaecological Surgery

► Obstetrical Drapes and Packs


C-section pack

Item no.	DYJPECSEPSM1	Description	Dimensions	Quantity
eClass-Nr.		Back table cover,	152x190 cm	1
Material	ADV - SMMMS	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
PU	8	C-section drape, incise film w/ pear shaped fen. 18x18 cm	196/259x307 cm	1
SU	0	Baby blanket	76x86 cm	1




C-section pack

Item no.	DYJPECSEPEC1	Description	Dimensions	Quantity
eClass-Nr. <td> <td>Hand towel</td> <td>34x36 cm</td> <td>4</td> </td>	<td>Hand towel</td> <td>34x36 cm</td> <td>4</td>	Hand towel	34x36 cm	4
Material	ECL - Spunlac	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
PU	8	Back table cover,	152x190 cm	1
SU	0	C-section drape, incise film w/ pear shaped fen. 18x18 cm	198/259x307 cm	1
		Baby blanket	76x86 cm	1



Delivery pack

Item no.	CP29162CE	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Absorbent Towel</td> <td>39x58 cm</td> <td>1</td>		Absorbent Towel	39x58 cm	1
Material	ESS - 2 layers	OR gown, Eclipse, XL	130 cm	1
PU	5	Drape blanket Baby, White	100x105 cm	1
SU	0	Under bottock drape with pouch	96,5x118 cm	1
		Leggings	140x79 cm	2
		Small drape, 75x75cm, without adhesive	75x75 cm	1
		Side sheet, adhesive	75x90 cm	1
		Back table cover	140x190 cm	1




Drapes & Packs

Gynaecological Surgery

► Obstetrical Drapes and Packs


Delivery pack

Item no.	ES10311CE	Description	Dimensions	Quantity
eClass-Nr.		Leggings	75x120 cm	2
Material	ESS - 2 layers	Adhesive tape	9x50 cm	1
PU	12	Drape blanket baby, white	100x105 cm	1
SU	555	Drape adhesive	75x75 cm	1
		SMS wrap		1
		UBD, 85x114cm	85x114 cm	1
		Absorbent Hand towels	20x30 cm	2
		Back table cover	140x190 cm	1




Delivery pack

Item no.	CP29162CEA	Description	Dimensions	Quantity
eClass-Nr. <td> <td>Absorbent Towel</td> <td>39x58 cm</td> <td>1</td> </td>	<td>Absorbent Towel</td> <td>39x58 cm</td> <td>1</td>	Absorbent Towel	39x58 cm	1
Material	ESS - 2 layers	OR gown, Eclipse, XL	130 cm	1
PU	7	Drape blanket Baby, White	100x105 cm	1
SU	0	Under buttock drape with pouch	96,5x118 cm	1
		Leggings	140x79 cm	2
		Small drape, 75x75cm, without adhesive	75x75 cm	1
		Side sheet, adhesive	75x90 cm	1
		Back table cover	140x190 cm	1



Delivery pack

Item no.	ES10310CEB	Description	Dimensions	Quantity
eClass-Nr. <td> <td>Absorbent Hand towels</td> <td>20x30 cm</td> <td>2</td> </td>	<td>Absorbent Hand towels</td> <td>20x30 cm</td> <td>2</td>	Absorbent Hand towels	20x30 cm	2
Material	ESS - 2 layers	Leggings	140x79cm	2
PU	7	Drape blanket Baby, White	100x105 cm	1
SU	0	Back table cover, 50µm	120x140 cm	1
		UBD, 85x114cm	85x114 cm	1
		OR gown, Essential, XL	130 cm	1
		CSR wrap	60x60 cm	1
		Drape	75x90cm	1




Drapes & Packs

Gynaecological Surgery

► Obstetrical Drapes and Packs

OB pack

Item no.	DYJPEOBPSM	Description	Dimensions	Quantity
eClass-Nr.		Adhesive drape	38x66 cm	1
Material	ADV - SMMMS	Leggings w/ cuff 15 cm	74x109 cm	2
PU	10	Under buttocks drape, 1 fluid collection pouch	102x112 cm	1
SU	0	Abdominal drape, adh.	135x97 cm	1
		Baby blanket	76x86 cm	1
		Back table cover, reinforced, 50µm	140x190 cm	1
		Hand towel	34x36 cm	4



Drapes & Packs

Gynaecological Surgery

► Obstetrical Drapes and Packs

Birth ball cover



Item no.	SPT-ESBB01
eClass-Nr.	
Material	ESS - 2 layers
Dimensions	
PU	96
SU	0

Drapes & Packs

Gynaecological Surgery

► Peri-Gynae Packs


Abdominal/Perineal pack

Item no.	CP29176CE	Description	Dimensions	Quantity
eClass-Nr.		Pouch, Instmt, double, PE	33x41 cm	1
Material	ECL - Spunlac	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
PU	4	Hand towels	30x20 cm	4
SU	0	Velcro strip	2.5x13 cm	1
		Side sheet, adhesive	75x90 cm	2
		Adhesive tape	9x50 cm	2
		Under buttocks drape with cuff 24cm, reinforced	96x126 cm	1
		LAVH drape, fenestration 24x29/7x15 cm with leggings, adhesive	250x337 cm	1
		Back table cover, 75µm	160x240 cm	1
		CSR wrap	60x60 cm	1




Drape Lithotomy

Item no.	29014CEA	Description	Dimensions	Quantity
eClass-Nr. <td>34-35-01-08</td> <td>Lithotomy drape, 160x300cm, with adh, with fen, with FCP</td> <td>160/267x300 cm</td> <td>1</td>	34-35-01-08	Lithotomy drape, 160x300cm, with adh, with fen, with FCP	160/267x300 cm	1
Material <td>Eclipse - Spunl</td> <td>Absorbent Hand towels</td> <td>20x30 cm</td> <td>1</td>	Eclipse - Spunl	Absorbent Hand towels	20x30 cm	1
PU <td>7</td> <td>Back table cover</td> <td>140x190 cm</td> <td>0</td>	7	Back table cover	140x190 cm	0
SU <td>0</td> <td>CSR wrap</td> <td>86x86 cm</td> <td>1</td>	0	CSR wrap	86x86 cm	1



Gynaecology Pack

Item no.	TB9175CEA	Description	Dimensions	Quantity
eClass-Nr. <td>34-35-01-08</td> <td>Under buttocks drape with cuff 24cm, reinforced</td> <td>96x126 cm</td> <td>1</td>	34-35-01-08	Under buttocks drape with cuff 24cm, reinforced	96x126 cm	1
Material <td>Ultimate - 3 la</td> <td>Drape, adhesive</td> <td>45x50 cm</td> <td>1</td>	Ultimate - 3 la	Drape, adhesive	45x50 cm	1
PU <td>10</td> <td>Leggings</td> <td>140x79 cm</td> <td>2</td>	10	Leggings	140x79 cm	2
SU <td>0</td> <td>Back table cover, reinforced, 50µm</td> <td>125x230 cm</td> <td>1</td>	0	Back table cover, reinforced, 50µm	125x230 cm	1
		Hand towels	30x20 cm	1
		Drape, adhesive	98x145 cm	1



Drapes & Packs


Gynaecological Surgery

► Peri-Gynae Packs

Gynae-uro pack,

Item no. TB29960CEA
eClass-Nr. 34-35-01-08
Material Ultimate - 3 la
PU 6
SU 0


Description	Dimensions	Quantity
Lithotomy drape with attached leggings, reinforced, fenestration 6x16cm	89/229x312 cm	1
Hand towels	30x20 cm	2
Under buttocks drape with cuff 24cm, reinforced	96x126 cm	1
Drape, adhesive	49x65 cm	1
Back table cover	140x190 cm	1



Gyn-cystoscopy pack

Item no. ES10360CEA
eClass-Nr. 34-35-01-08
Material ESS - 2 layers
PU 17
SU 0


Description	Dimensions	Quantity
Absorbent Hand towels	20x30 cm	2
CSR wrap	86x86 cm	1
Leggings	140x79cm	2
Cystoscopy sheet, fenestration with adhesive 8*18cm	120x75 cm	1
Back table cover	140x190 cm	1



Gyn-cystoscopy pack

Item no. ES10362CE
eClass-Nr. ESS - 2 layers
Material 14
PU 555
SU 0

Description	Dimensions	Quantity
Adhesive tape	9x50 cm	1
Leggings	75x120 cm	2
Drape adhesive	75x75 cm	2
SMS wrap		1
Back table cover, reinforced, 50µm	140x190 cm	1
Absorbent Hand towels	20x30 cm	2




Drapes & Packs

Gynaecological Surgery

► Peri-Gynae Packs


Gynecology/Urology Pack

Item no.	SM8007CE	Description	Dimensions	Quantity
eClass-Nr.		Hand towels	30x20 cm	1
Material	ADV - SMMMS	gynecological/urology drape, 80/255x289cm	80/255x289cm	1
PU	10	SMS wrap	60x60 cm	1
SU	555	Back table cover, reinforced, 50µm	140x190 cm	1




Gynecology/Urology Pack

Item no.	ES10361CE	Description	Dimensions	Quantity
eClass-Nr.		Adhesive tape	9x50 cm	1
Material	ESS - 2 layers	Gynecology drape, 12x25 cm fenestration with adhesive, leggings	280x175/242cm	1
PU	8	Gynecology drape, 9x12cm heart shaped adhesive fenestration, wi..	35/40x114cm	1
SU	555	Back table cover, reinforced, 50µm	140x190 cm	1
		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
		Absorbent Hand towels	20x30 cm	2



Gynecology/Urology Pack

Item no.	ES10805CE	Description	Dimensions	Quantity
eClass-Nr.		Adhesive tape	9x50 cm	1
Material	ESS - 2 layers	Drape Lithotomy fen 10 x 17 cm ecc w/ adh leg	280/175 x 242 cm	1
PU	10	drape sheet	75x90cm	1
SU	0	SMS wrap		1
		Back table cover, reinforced, 50µm	140x190 cm	1
		Absorbent Hand towels	20x30 cm	2




Drapes & Packs

Gynaecological Surgery

► Peri-Gynae Packs


Lithotomy drape

Item no.	DYNJPE9004	Description	Dimensions	Quantity
eClass-Nr.		Wrap	100x100 cm	1
Material	ECL - Spunlac	Lithotomy drape, adh. abdominal fen. 15x30cm, lower fen. 8x15cm	254x333 cm	1
PU	8			
SU	0			




Lithotomy drape

Item no.	DYNJPE9001	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	Eclipse - Spunl	Lithotomy drape, adh. perineal fen. 8x15 cm, pouch	74/163x279 cm	1
PU	8			
SU	0			



Lithotomy pack II

Item no.	DYNJPE9010SM	Description	Dimensions	Quantity
eClass-Nr.		Adhesive drape	38x66 cm	1
Material	ADV - SMMMS	Leggings w/ cuff 15 cm	74x109 cm	2
PU	12	Blue PE under buttocks drape, wetlaid reinf.	89x119 cm	1
SU	0	Back table cover	125x229 cm	1
		Hand towel	51x43 cm	2
		Abdominal drape, adh.	135x97 cm	1




Drapes & Packs

Gynaecological Surgery

► Peri-Gynae Packs


Lithotomy pack II

Item no.	DYNJPE9000T	Description	Dimensions	Quantity
eClass-Nr.		Lithotomy drape, adh. fen. 8x15 cm, leggings	71/259x203 cm	1
Material	Eclipse - Spunl	Blue PE under buttocks drape, wetlaid reinf.	89x119 cm	1
PU	10	Back table cover	125x229 cm	1
SU	0	Hand towel	51x43 cm	1




Lithotomy pack II

Item no.	DYNJPE9010	Description	Dimensions	Quantity
eClass-Nr.		Adhesive drape	38x66 cm	1
Material	Eclipse - Spunl	Leggings w/ cuff 15 cm	74x109 cm	2
PU	12	Blue PE under buttocks drape, wetlaid reinf.	89x119 cm	1
SU	0	Back table cover	125x229 cm	1
		Hand towel	51x43 cm	2
		Abdominal drape, adh.	135x97 cm	1



Peri-Gyn / TVT pack

Item no.	ES10806CE	Description	Dimensions	Quantity
eClass-Nr.		Peri-Gyn drape, w/leggings, w/pouch	220/300x279cm	1
Material	ESS - 2 layers	Velcro strip	2.5x13 cm	2
PU	6	Under buttocks drape with cuff 24cm, reinforced	96x126 cm	1
SU	0	Adh extendable, turquoise	9x50 cm	1
		CSR wrap	60x60 cm	1
		Back table cover	140x190 cm	1
		Absorbent Hand towels	20x30 cm	2



Drapes & Packs


Gynaecological Surgery

► Peri-Gynae Packs

Peri-gynae pack,

Item no. TB29170CEA
eClass-Nr. 34-35-01-08
Material Ultimate - 3 la
PU 8
SU 0


Description	Dimensions	Quantity
Lithotomy drape with attached leggings, reinforced, fenestration 6x16cm	89/229x312 cm	1
Under buttocks drape with cuff 24cm, reinforced	96x126 cm	1
Drape, adhesive	45x50 cm	1
SMS wrap		1
Back table cover, reinforced, 50µm	140x190 cm	1
Hand towels	30x20 cm	1



Peri-gynae pack,

Item no. TB29176CEA
eClass-Nr. 34-35-01-08
Material Ultimate - 3 la
PU 14
SU 0


Description	Dimensions	Quantity
Lithotomy T drape	137x192 cm	1
Drape, adhesive	45x50 cm	1
Leggings	140x79cm	2
Back table cover, reinforced, 50µm	125x230 cm	1
Hand towels	30x20 cm	1



Perineal pack

Item no. CP29129CE
eClass-Nr. ECL - Spunlac
Material 6
PU 0
SU 0

Description	Dimensions	Quantity
Hand towels	30x39 cm	1
Small drape	45x50 cm	1
Lithotomy drape, 160x300cm, with adh, with fen, with FCP	160/267x300 cm	1
Back table cover	140x190 cm	1
CSR wrap	86x86 cm	1



Drapes & Packs

Gynaecological Surgery

► Peri-Gynae Packs


TVT drape, pack

Item no.	DG29152	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-08	Absorbent Towel	39x58 cm	2
Material	ECL - Spunlac	Incise film	13 x 15cm	1
PU	6	Back table cover	140x190 cm	1
SU	0	CSR wrap	86x86 cm	1
		TVT drape	160/ 300 x 90cm	1



Universal pack, Lithotomy position

Item no.	CP29148CE	Description	Dimensions	Quantity
eClass-Nr.		Back table cover, folded, 50µm	140x190 cm	1
Material	ECL - Spunlac	Leggings	140x79cm	2
PU	4	Under buttocks drape with cuff 24cm, reinforced	96x126 cm	1
SU	0	Side drape, reinforced, adhesive	98x110 cm	3
		Head Drape, w/ adh, velcro	240x196 cm	1
		Instrument pouch with adhesive, 2 compartments	43,5x49 cm	1
		Back table cover	140x190 cm	1
		CSR wrap	112x112 cm	1



Drapes & Packs

Gynaecological Surgery

► Speciality


Delivery pack

Item no.	DG29C293A	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-08	Under buttocks drape with cuff 24cm, reinforced	96x126 cm	1
Material	ECL - Spunlac	Drape, reinforced, adhesive	240x196 cm	1
PU	8	Drape Sheet, 98x110cm, w/adh MED 5x95cm 2 sides	98x110 cm	1
SU	0	Back table cover, reinforced, 50µm	140x190 cm	1
		Adhesive tape, PE	8x40 cm	3
		Leggings	140x79 cm	2




Gynaecology pack

Item no.	TB27051CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-08	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
Material	Ultimate - 3 la	Absorbent Hand towels	20x30 cm	1
PU	10	Velcro strip	2.5x13 cm	1
SU	0	Under buttocks drape with cuff 24cm, reinforced	96x126 cm	1
		Adhesive drape	75x91 cm	1
		Gynecology drape, 10x11 cm adhesive fenestration	260/278x260 cm	1
		Back table cover	140x190 cm	1



Gynaecology pack

Item no.	DYJPEGYNPEC	Description	Dimensions	Quantity
eClass-Nr.	Eclipse - Spunl	Blue PE under buttocks drape, wetlaid reinf.	89x119 cm	1
Material	7	Gynecology drape, adh., adh. fen. 10x12 cm	260x311 cm	1
PU	0	Mayo stand cover, SMS reinf. 55x88 cm	80x142 cm	1
SU		Back table cover,	152x190 cm	1
		Hand towel	34x36 cm	1




Drapes & Packs

Gynaecological Surgery

► Speciality


Gynecology pack

Item no.	DYJPEGYNPSM	Description	Dimensions	Quantity
eClass-Nr.		Blue PE under buttocks drape, wetlaid reinf.	89x119 cm	1
Material	ADV - SMMMS	Gynecology drape, adh. fen. 10x12 cm	260x304 cm	1
PU	7	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
SU	0	Back table cover,	152x190 cm	1
		Hand towel	34x36 cm	1




Pfannenstiel pack

Item no.	DYJPEPLPSM1	Description	Dimensions	Quantity
eClass-Nr.		Pfannenstiel drape, fen. w/ full incise film	193/240x326 cm	1
Material	ADV - SMMMS	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
PU	8	Back table cover,	152x190 cm	1
SU	0			



Pfannenstiel pack

Item no.	DYJPEPLPEC1	Description	Dimensions	Quantity
eClass-Nr.		Pfannenstiel drape, fen. w/ fen. incise film	193/240x325 cm	1
Material	Eclipse - Spunl	Mayo stand cover, SMS reinf. 55x88 cm	80x142 cm	1
PU	8	Back table cover,	152x190 cm	1
SU	0			




Drapes & Packs

Gynaecological Surgery

► Speciality

Pfannenstiel pack I,

Item no.	TB27043CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-08	Pfannenstiel drape, 27.5x21 cm adhesive fenestration	250/198X317 cm	1
Material	Ultimate - 3 la	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
PU	8	Back table cover	140x190 cm	1
SU	0			




Drapes & Packs

NeuroSurgery


Craniotomy drape

Item no.	DYNJPE10040	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	Eclipse - Spunl	Craniotomy drape, 30x20 cm adh. fen. w/ full incise film, pouch	188/297x343 cm	1
PU	9			
SU	0			




Craniotomy pack

Item no.	CP29544CE	Description	Dimensions	Quantity
eClass-Nr.		Hand towels	30x39 cm	2
Material	ECL - Spunlac	Adhesive tape, PE	8x40 cm	1
PU	4	Drape, adhesive	49x65 cm	4
SU	0	Back table cover, 75µm, folded	140x190cm	1
		Drape sheet	196x249 cm	1
		Neurological Drape, with incise film and pouch	241x254 cm	1
		Back table cover, 75µm	140x190 cm	1
		CSR wrap	112x112 cm	1



Craniotomy pack

Item no.	CP29008CE	Description	Dimensions	Quantity
eClass-Nr.		Neurological Drape, with incise film and pouch	241x254 cm	1
Material	Eclipse - Spunl	Drape sheet, adhesive	45x49 cm	6
PU	4	Drape, reinforced, adhesive	240x196 cm	1
SU	0	Adhesive tape, PE	8x40 cm	4
		CSR wrap	112x112 cm	1
		Back table cover, folded, 50µm	140x120 cm	1
		Back table cover, reinforced, 50µm	140x190 cm	1




Drapes & Packs

NeuroSurgery

Craniotomy pack

Item no.	ES10548CE	Description	Dimensions	Quantity
eClass-Nr.		Craniotomy drape	236x406cm	1
Material	ESS - 2 layers	Drape, adhesive	45x50 cm	4
PU	6	Adhesive tape, PE	8x40 cm	2
SU	0	Back table cover, reinforced, 50µm	140x190 cm	1
		Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	2




Craniotomy pack,

Item no.	29244FRD	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Neurological Drape, with incise film and pouch	241x254 cm	1
Material	ECL - Spunlac	Drape sheet	196x249 cm	1
PU	6	Adhesive tape, PE	8x40 cm	1
SU	0	Back table cover, 75µm, folded	140x190cm	1
		Hand towels	30x39 cm	2
		Drape , oval neck diameter 15cm adhesive	48 x 63 cm	2
		Back table cover, reinforced, 50µm	140x190 cm	1
		CSR wrap	85x85 cm	1

No picture available

Neurological drape,

Item no.	TB9448CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	drape with pouch, reinf, oval fenestration with incise film 28x18cm	254x241 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	12			
SU	0			




Drapes & Packs

NeuroSurgery


Neurological pack

Item no.	CP9448CEA	Description	Dimensions	Quantity
eClass-Nr.		Pouch, Instmt, double, PE	33x41 cm	1
Material	Spunlace	Craniotomy drape, 30x20 cm adh. fen. w/ full incise film, pouch	188/297x343 cm	1
PU	4	Head drape, reinforcement, adhesive, 1 cord holding tab	254x183 cm	1
SU	0	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
		back table cover, 90 µm	203x240 cm	1
		Hand towels	30x20 cm	4
		Velcro strip	2.5x13 cm	1
		Adhesive tape	9x50 cm	3
		Drape sheet, adhesive	45x49 cm	4
		CSR wrap	60x60 cm	1




Neurological pack

Item no.	DYJPENEUPSM	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
Material	ADV - SMMMS	Craniotomy drape, fen. w/ full incise film 28x20 cm, adh.	264x252 cm	1
PU	5	Back table cover,	152x190 cm	1
SU	0	Hand towel	34x36 cm	1



Neurological pack

Item no.	DYJPENEUEPC	Description	Dimensions	Quantity
eClass-Nr.		Craniotomy drape, fen. w/ full incise film 28x20 cm, adh.	264x254 cm	1
Material	ECL - Spunlac	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
PU	5	Back table cover,	152x190 cm	1
SU	0			




Drapes & Packs

NeuroSurgery


Spine drape

Item no.	DYNJPE10011	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	ECL - Spunlac	Spine drape, fen. 10x41 cm w/ full incise film	196/265x307 cm	1
PU	12			
SU	0			




Spine pack

Item no.	CP29125CE	Description	Dimensions	Quantity
eClass-Nr.		Head drape, reinforced, adhesive	270x183 cm	2
Material	Eclipse - Spunl	Side drape, reinforced, adhesive	98x110 cm	2
PU	6	Adhesive tape, PE	8x40 cm	2
SU	0	CSR wrap	112x112 cm	1
		Back table cover	140x190 cm	1



Spine pack

Item no.	CPU29242CE	Description	Dimensions	Quantity
eClass-Nr.		Adhesive tape, PE	8x40 cm	2
Material	Ultimate - 3 la	not found...	49x65 cm	4
PU	8	Laparotomy drape, adh. fen. 23x37 cm	196/254x314 cm	1
SU	0	CSR wrap	102x102 cm	1
		Back table cover	140x190 cm	1




Drapes & Packs

NeuroSurgery

Spine pack

Item no.	CPU29033CE	Description	Dimensions	Quantity
eClass-Nr.		spine drape	220x320cm	1
Material	Ultimate - 3 la	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
PU	10	Back table cover, 75µm	140x190 cm	1
SU	0			



Universal spine drape,

Item no.	TB29418CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	drape with armboard covers, reinf, fen 10x30cm adh 2 pouches, 4 velcro	285/251x351 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	7			
SU	0			




Drapes & Packs

Ophtalmic Surgery


Eye drape

Item no.	DYNJPE7004	Description	Dimensions	Quantity
eClass-Nr.		Eye drape, full incise film 11x13 cm, 2 pouches	165x178 cm	1
Material	wetlaid	Wrap	76x76 cm	1
PU	10			
SU	0			



Eye drape

Item no.	DYNJPE7005	Description	Dimensions	Quantity
eClass-Nr.		Eye drape, fen. w/ incise film 11x13 cm, 2 pouches	165x178 cm	1
Material	wetlaid	Wrap	76x76 cm	1
PU	10			
SU	0			



Ophthalmic drape

Item no.	TB9562CE	Description	Dimensions	Quantity
eClass-Nr.		Ophthalmic drape with pouch	196x279 cm	1
Material	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU	16			
SU	0			




Drapes & Packs

Ophthalmic Surgery

Ophthalmic drape

Item no.	DYNJPE7011	Description	Dimensions	Quantity
eClass-Nr.		Wrap	66x66 cm	1
Material	Eclipse - Spunl	Ophthalmic drape, incise film w/4x6 cm fen., pouch	196x254 cm	1
PU	20			
SU	0			



Ophthalmic drape,

Item no.	S8441CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-05	Ophthalmic drape, oval fenestration 8x5cm, adhesive	113x158 cm	1
Material	ADV - SMMMS	CSR wrap	112x112 cm	1
PU	60			
SU	0			



Ophthalmic drape, apertured incise, w/pouch

Item no.	S9563CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-05	Ophthalmic Drape, incise film fenestrated, pouch	157x183 cm	1
Material	ADV - SMMMS	CSR wrap	76x76 cm	1
PU	24			
SU	0			




Drapes & Packs

Ophtalmic Surgery


Ophthalmic incise drape with fluid pouch

Item no.	IDE1062	Description	Dimensions	Quantity
eClass-Nr.		Ophthalmic incise drape with fluid pouch	122x129 cm	1
Material	SMMMS			
PU	10			
SU	0			




Ophthalmic pack

Item no.	CP29024CE	Description	Dimensions	Quantity
eClass-Nr.		Wrap	72x72 cm	1
Material	ADV - SMMMS	drape sheet	75x90cm	3
PU	16	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
SU	0	Ophthalmic Drape, incise film fenestrated, pouch	157x183 cm	1



Ophthalmic pack

Item no.	DYJPEOPHPSM2	Description	Dimensions	Quantity
eClass-Nr.		Eye drape, fen. w/ incise film 11x13 cm, 2 pouches	165x178 cm	1
Material	ADV - SMMMS	Back table cover	140x190 cm	1
PU	12	SMS wrap		1
SU	0			




Drapes & Packs

Ophthalmic Surgery

Ophthalmic pack

Item no.	DYJPEOPHPSM1	Description	Dimensions	Quantity
eClass-Nr.		Eye drape, full incise film 11x13 cm, 2 pouches	165x178 cm	1
Material	ADV - SMMMS	Back table cover	140x190 cm	1
PU	12	SMS wrap		1
SU	0			




Ophthalmology Pack

Item no.	CP29117CE	Description	Dimensions	Quantity
eClass-Nr.		Absorbent Towel	39x58 cm	3
Material	ADV - SMMMS	OR gown, Essential, L	120 cm	1
PU	5	OR gown, Essential, XL	130 cm	2
SU	0	Ophthalmic drape with pouch	196x279 cm	1
		Back table cover	140x190 cm	1



Ophthalmology pack

Item no.	CP95622CE	Description	Dimensions	Quantity
eClass-Nr.		Ophthalmic Drape, incise film fenestrated, pouch	157x183 cm	1
Material	ADV - SMMMS	Back table cover, reinforced, 50µm	140x190 cm	152
PU	16	CSR wrap	76x76 cm	1
SU	0			




Drapes & Packs

Ophthalmic Surgery


Ophthalmology pack,

Item no.	TB25621CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-05	Ophthalmic drape, "Alcon" style folding	196x275cm	1
Material	Ultimate - 3 la	Back table cover	140x190 cm	1
PU	16	CSR wrap	86x86 cm	1
SU	0			



Strabismus pack

Item no.	CP30458CE	Description	Dimensions	Quantity
eClass-Nr.		Ophthalmic drape, w/ adh fen ovale 15x5cm	110x240cm	1
Material	Ultimate - 3 la	Hand towels	30x20 cm	2
PU	20	Back table cover	140x190 cm	1
SU	0	CSR wrap	76x76 cm	1



Drapes & Packs

Ophtalmic Surgery

Ophthalmic drape, Apertd Inc
4 x 6 cm, w/ pouch



Item no.	S9562CE
eClass-Nr.	34-35-01-05
Material	ADV - SMMMS
Dimensions	190x279 cm
PU	20
SU	0

Ophthalmic drape



Item no.	S9565CE
eClass-Nr.	34-35-01-05
Material	ADV - SMMMS
Dimensions	196x291 cm
PU	16
SU	0

Ophthalmic drape with
pouch, full incise



Item no.	S7444CE
eClass-Nr.	34-35-01-05
Material	ADV - SMMMS
Dimensions	122x157 cm
PU	34
SU	0

Ophthalmic drape,
fenestrated incise



Item no.	S7445CE
eClass-Nr.	34-35-01-05
Material	ADV - SMMMS
Dimensions	122x157 cm
PU	34
SU	0

Drapes & Packs

Orthopedic Surgery

Split drape

Item no.	ES10425CE	Description	Dimensions	Quantity
eClass-Nr.		SMS wrap		1
Material	ESS - 2 layers	U-drape, split 65x7cm, adhesive	150 x 200 cm	1
PU	32			
SU	555			



Split drape

Item no.	ES29425CE	Description	Dimensions	Quantity
eClass-Nr.		SMS wrap		1
Material	ESS - 2 layers	U-drape, split 65x7cm, reinforced, adhesive	150 x 200 cm	1
PU	42			
SU	555			



Drapes & Packs

Orthopedic Surgery

▶ Arthroscopy Drapes and Packs


Arthroscopy drape

Item no.	DYJPEARDTB	Description	Dimensions	Quantity
eClass-Nr.		arthroscopy drape	231/290x315 cm	1
Material	Ultimate - 3 la	CSR wrap	112x112 cm	1
PU	6			
SU	0			



Arthroscopy drape

Item no.	DYNJPE8101	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	ECL - Spunlac	Arthroscopy drape, elastic fen. 6 cm diam., pouch	229x300 cm	1
PU	12			
SU	0			



Arthroscopy drape

Item no.	ES9414CE	Description	Dimensions	Quantity
eClass-Nr.		Arthroscopy drape with pouch	230x300 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	16			
SU	555			




Drapes & Packs

Orthopedic Surgery

▶ Arthroscopy Drapes and Packs

Arthroscopy Drape with Armboard Covers

Item no.	TB9423CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Arthroscopy drape	234/295x447 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	8			
SU	0			



Arthroscopy drape,

Item no.	TB9414CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Arthroscopy drape, elastic fenestration diameter 6cm	234x302 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	10			
SU	0			



Arthroscopy drape,

Item no.	TB9413CE	Description	Dimensions	Quantity
eClass-Nr.	Ultimate - 3 la	Arthroscopy drape, elastic fenestration diameter 6cm	234x302 cm	1
Material	12	CSR wrap	86x86 cm	1
PU	0			
SU	0			




Drapes & Packs

Orthopedic Surgery

▶ Arthroscopy Drapes and Packs

Arthroscopy Pack

Item no.	CP29055CE	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
Material	Ultimate - 3 la	Absorbent Towel	39x58 cm	4
PU	6	Camera cover, with elastoflex	17.5 x 240 cm	1
SU	0	Botee	55x37.5 cm	1
		Adhesive tape	9x50 cm	2
		Foot drape, adhesive	175x175 cm	1
		Pouch, Instmt, double, PE	33x41 cm	1
		Arthroscopy drape, elastic fenestration diameter 6cm	234x302 cm	1
		Back table cover	140x190 cm	1
		CSR wrap	86x86 cm	1




Arthroscopy pack

Item no.	DYJPEARPSM1	Description	Dimensions	Quantity
eClass-Nr.		Impervious stockinette	23x91 cm	1
Material	ADV - SMMMS	Arthroscopy drape, elastic fen. 6 cm, w/ pouch	221/290x310 cm	1
PU	6	Drape sheet	135x196 cm	1
SU	0	Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1
		OP tape	10x50 cm	2
		Back table cover,	152x190 cm	1
		Hand towel	34x36 cm	2



Arthroscopy pack

Item no.	DYJPEARPEC1	Description	Dimensions	Quantity
eClass-Nr.		Impervious stockinette	23x91 cm	1
Material	ECL - Spunlac	Arthroscopy drape, elastic fen. 6cm, fluid pouch	221/290x310 cm	1
PU	6	Drape sheet	135x196 cm	1
SU	0	OP tape	10x50 cm	2
		Back table cover,	152x190 cm	1
		Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1




Drapes & Packs

Orthopedic Surgery

► Arthroscopy Drapes and Packs


Arthroscopy Pack

Item no.	CP29145CE	Description	Dimensions	Quantity
eClass-Nr.		Velcro strip	2.5x13 cm	1
Material	Eclipse - Spunl	Mayo stand cover, ultra absorbent reinforcement 75x100 cm	80x145 cm	1
PU	4	Side drape	75x98 cm	1
SU	0	Extremity sheet, elastic fenestration diameter 6cm	221x326 cm	1
		Video camera cover with cardboard	17.5x240 cm	1
		OR gown, Eclipse, XLL, fabric reinforced	150 cm	1
		OR gown, XL, 130 cm, fabric reinforced	130 cm	1
		CSR wrap	112x112 cm	2
		Hand towels	30x20 cm	2
		Half drape	112x145 cm	2




Arthroscopy pack

Item no.	TB29128CE	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Velcro strip</td> <td>2.5x13 cm</td> <td>2</td>		Velcro strip	2.5x13 cm	2
Material	Ultimate - 3 la	Arthroscopy drape, elastic fenestration diameter 6cm	234x302 cm	1
PU	8	Back table cover, folded, 50µm	140x190 cm	1
SU	0	CSR wrap	86x86 cm	1
		Back table cover	140x190 cm	1



Arthroscopy Pack

Item no.	TB9416CE	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Arthroscopy drape, elastic fenestration diameter 6cm</td> <td>234x302 cm</td> <td>1</td>		Arthroscopy drape, elastic fenestration diameter 6cm	234x302 cm	1
Material	Ultimate - 3 la	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
PU	8	Back table cover	140x190 cm	1
SU	0			



Drapes & Packs

Orthopedic Surgery

▶ Arthroscopy Drapes and Packs


Arthroscopy Pack

Item no.	TB9185CE	Description	Dimensions	Quantity
eClass-Nr.		Adhesive tape	9x50 cm	4
Material	Ultimate - 3 la	Mayo stand cover, reinforced	80x142 cm	1
PU	4	Hand towels	30x39 cm	4
SU	0	Suture bag, adh.	16x31 cm	1
		Arthroscopy drape, elastic fenestration diameter 6cm	234x302 cm	1
		Impervious stockinette	23x122 cm	1
		OR gown, Eclipse, L	120 cm	1
		U-drape PE, split 15x53cm, adhesive	152x183 cm	1
		R gown, Eclipse, XXL, polyreinforced	150 cm	1
		Back table cover, reinforced, 50µm	125x230 cm	1
		CSR wrap	60x60 cm	1
		CSR wrap	112x112 cm	1




Arthroscopy Pack

Item no.	ES10128CE	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Arthroscopy drape with pouch</td> <td>230x300 cm</td> <td>1</td>		Arthroscopy drape with pouch	230x300 cm	1
Material <td>ESS - 2 layers</td> <td>Back table cover</td> <td>140x190 cm</td> <td>1</td>	ESS - 2 layers	Back table cover	140x190 cm	1
PU <td>12</td> <td>SMS wrap</td> <td></td> <td>1</td>	12	SMS wrap		1
SU <td>555</td> <td></td> <td></td> <td></td>	555			



Arthroscopy Pack

Item no.	ES10129CE	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Arthroscopy T-Drape with Pouch</td> <td>290/221x322 cm</td> <td>1</td>		Arthroscopy T-Drape with Pouch	290/221x322 cm	1
Material <td>ESS - 2 layers</td> <td>Drape sheet</td> <td>150x150 cm</td> <td>1</td>	ESS - 2 layers	Drape sheet	150x150 cm	1
PU <td>10</td> <td>SMS wrap</td> <td></td> <td>1</td>	10	SMS wrap		1
SU <td>555</td> <td>Back table cover</td> <td>140x190 cm</td> <td>1</td>	555	Back table cover	140x190 cm	1




Drapes & Packs

Orthopedic Surgery

► Arthroscopy Drapes and Packs


Arthroscopy Pack

Item no.	TB29051CE	Description	Dimensions	Quantity
eClass-Nr.		Impervious stockinette	20x97 cm	1
Material	Ultimate - 3 la	Arthroscopy drape, elastic fenestration diameter 6cm	234x302 cm	1
PU	8	Back table cover, reinforced, 50µm	140x190 cm	1
SU	0	Adhesive tape	9x50 cm	2




Arthroscopy Pack

Item no.	ES10416CE	Description	Dimensions	Quantity
eClass-Nr.		Video camera cover with pre-cut tip, with cardboard tab for applicat..	13x240 cm	1
Material	ESS - 2 layers	Arthroscopy T-Drape with Pouch	290/221x322 cm	1
PU	8	Stockinette	36.5x72 cm	1
SU	555	Adhesive tape	9x50 cm	2
		Back table cover,	152x190 cm	2
		Mayo Stand Cover, absorbent reinf. 55x88 cm	80x142 cm	1
		Absorbent Hand towels	20x30 cm	2



Arthroscopy Pack

Item no.	ES29051CE	Description	Dimensions	Quantity
eClass-Nr.		Arthroscopy drape with pouch	230x300 cm	1
Material	ESS - 2 layers	Drape sheet	150x150 cm	1
PU	10	Stockinette	36.5x72 cm	1
SU	555	Adhesive tapes	9x50 cm	2
		SMS wrap		1
		Back table cover, reinforced, 50µm	140x190 cm	1
		Absorbent Hand towels	20x30 cm	2




Drapes & Packs

Orthopedic Surgery

▶ Arthroscopy Drapes and Packs

Arthroscopy T-drape

Item no.	DYNJPE8102	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	ECL - Spunlac	Arthroscopy drape, elastic fen. 6cm, fluid pouch	221/290x310 cm	1
PU	6			
SU	0			




Arthroscopy T-Drape with pouch

Item no.	DYNJPE8102SM	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	ADV - SMMMS	Arthroscopy drape, elastic fen. 6 cm, w/ pouch	221/290x310 cm	1
PU	6			
SU	0			



Cross ligament pack

Item no.	RG29114	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Impervious stockinette	23x122 cm	1
Material	Ultimate - 3 la	Mayo stand cover, reinforced	80x142 cm	1
PU	3	Arthroscopy drape	234/295x447 cm	1
SU	0	Instrument pouch in clear plastic with adhesive	33x41cm	1
		CSR wrap	95x95 cm	1
		U-drape PE, split 15x53cm, adhesive	152x183 cm	1
		Hand towels	30x20 cm	1
		Back table cover, 90µm	160x240 cm	1
		Adhesive tape	9x50 cm	2



Drapes & Packs

Orthopedic Surgery

► Extremity Drapes and Packs


Bilateral Extremity Drape

Item no.	29417CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Bilateral Extremity Drap E, 273 x 351cm	273x351 cm	1
Material	ECL - Spunlac	CSR wrap	76x76 cm	1
PU	10			
SU	0			




Bilateral extremity drape

Item no.	TB29417WCE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Extremity drape bilateral	273x381 cm	1
Material	Ultimate - 3 la	Back table cover	140x190 cm	1
PU	8	CSR wrap	86x86 cm	1
SU	0			



Bilateral extremity drape

Item no.	DYNJPE8004	Description	Dimensions	Quantity
eClass-Nr.		Wrap	66x66 cm	1
Material	Eclipse - Spunl	Bilateral extremity drape, 2 elastic fen. 6 cm diam.	229x307 cm	1
PU	12			
SU	0			




Drapes & Packs

Orthopedic Surgery

► Extremity Drapes and Packs

Bilateral Extremity pack

Item no.	CP29550CE	Description	Dimensions	Quantity
eClass-Nr.		Back table cover, folded, 50µm	140x120 cm	1
Material	Ultimate - 3 la	Extremity drape bilateral	273x381 cm	1
PU	6	Back table cover	140x190 cm	1
SU	0	CSR wrap	86x86 cm	1




Bilateral Limb Drape

Item no.	29416CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Bilateral Limb Drape 183 x 305 cm	183x305 cm	1
Material	ECL - Spunlac	CSR wrap	76x76 cm	1
PU	10			
SU	0			



Bilateral limb drape

Item no.	TB29416WCE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Bilateral limb drape with fenestration	183x305 cm	1
Material	Ultimate - 3 la	Back table cover	140x190 cm	1
PU	7	CSR wrap	76x76 cm	1
SU	0			




Drapes & Packs

Orthopedic Surgery

► Extremity Drapes and Packs


Bilateral Limb T-drape

Item no.	DYNJPE8006SM	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	ADV - SMMMS	Bilateral limb T-drape, 2 elastic fen. 6 cm	221/272x351 cm	1
PU	8			
SU	0			




Bilateral Limb T-drape

Item no.	DYNJPE8006	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	ECL - Spunlac	Bilateral limb T-drape, 2 elastic fen. 6 cm	229/272x351 cm	1
PU	8			
SU	0			



Extrem. Pck, fen/ circ/ elas 6 cm

Item no.	29414CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Extrem. Drape, fen/ circ/ elas 6 cm	221 x 326 cm	1
Material	ECL - Spunlac	CSR wrap	76x76 cm	1
PU	12			
SU	0			



Drapes & Packs

Orthopedic Surgery

► Extremity Drapes and Packs


Extremity drape

Item no.	DYJPEEXDSM	Description	Dimensions	Quantity
eClass-Nr.		Extremity drape, elastic fen. 6 cm	221x328 cm	1
Material	ADV - SMMMS	Wrap	91x91 cm	1
PU	10			
SU	0			




Extremity drape

Item no.	DYJPEEXDEC2	Description	Dimensions	Quantity
eClass-Nr.		Extremity drape, 4 cm diam. elastic fen.	229x328 cm	1
Material	ECL - Spunlac	Wrap	91x91 cm	1
PU	10			
SU	0			



Extremity drape

Item no.	DYJPEEXDEC	Description	Dimensions	Quantity
eClass-Nr.		Extremity drape, elastic fen. 6 cm	229x328 cm	1
Material	Eclipse - Spunl	Wrap	91x91 cm	1
PU	10			
SU	0			



Drapes & Packs

Orthopedic Surgery

► Extremity Drapes and Packs

Extremity drape

Item no.	TB29414CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Extremity drape, reinforced, elastic fenestration diameter 6cm	230x325 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	11			
SU	0			



Extremity drape

Item no.	ES10414CE	Description	Dimensions	Quantity
eClass-Nr.		Standard extremity drape, fen diam 6cm	325x230 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	18			
SU	555			



Extremity drape

Item no.	ES29414CE	Description	Dimensions	Quantity
eClass-Nr.		Extremity drape, reinforced, elastic fenestration diameter 6cm	230x325 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	14			
SU	555			



Drapes & Packs

Orthopedic Surgery

► Extremity Drapes and Packs


Extremity drape,

Item no.	TB29405CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Extremity drape, reinforced, elastic fenestration diameter 6cm	269/221x325 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	9			
SU	0			




Extremity Drape, XL

Item no.	RG29413	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Extremity sheet, reinforced, elastic fenestration 27x29cm	300x303 cm	1
Material	ECL - Spunlac	CSR wrap	85x85 cm	1
PU	10			
SU	0			



Extremity pack

Item no.	TB29188CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Hand towels	30x20 cm	2
Material	Ultimate - 3 la	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
PU	7	Back table cover, reinforced, 50µm	125x230 cm	1
SU	0	CSR wrap	86x86 cm	1
		Back table cover, folded, 50µm	140x190 cm	1
		Suture bag	25x34 cm	1
		Extremity drape, reinforced, elastic fenestration diameter 6cm	230x325 cm	1



Drapes & Packs


Orthopedic Surgery

► Extremity Drapes and Packs

Extremity pack

Item no. TB30012CE
eClass-Nr. 34-35-01-02
Material Ultimate - 3 la
PU 10
SU 0


Description	Dimensions	Quantity
Hook and loop tube holder	14x2.5 cm	1
Extremity drape, elastic fenestration diameter 5cm	242x300 cm	1
Back table cover	140x190 cm	1
CSR wrap	86x86 cm	1



Extremity pack

Item no. TB9186CE
eClass-Nr. 34-35-01-02
Material Ultimate - 3 la
PU 8
SU 0


Description	Dimensions	Quantity
Extremity drape, reinforced, elastic fenestration diameter 6cm	230x325 cm	1
Back table cover, folded, 50µm	140x190 cm	1
Back table cover	140x190 cm	1
CSR wrap	86x86 cm	1



Extremity pack

Item no. ES10409CEB
eClass-Nr. 34-35-01-02
Material ESS - 2 layers
PU 14
SU 0

Description	Dimensions	Quantity
Hand towels	30x20 cm	2
CSR wrap	86x86 cm	1
Adhesive tapes	9x50 cm	1
Extremity drape, elast. Fenestration 3.5cm	300 x 230 cm	1
Back table cover	140x190 cm	0




Drapes & Packs

Orthopedic Surgery

► Extremity Drapes and Packs


Extremity pack

Item no.	29020CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Drape	90x150 cm	1
Material	ECL - Spunlac	U-drape PE, split 20x76cm, adhesive	152x213 cm	1
PU	6	Extremity sheet w/ armboard covers	221/269x325cm	1
SU	0	Back table cover	140x190 cm	1
		CSR wrap	86x86 cm	1




Extremity pack

Item no.	29173CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Extremity sheet, elastic fenestration diameter 6cm	221x326 cm	1
Material	ECL - Spunlac	Hand towels	30x20 cm	2
PU	8	U-drape PE, split 20x76cm, adhesive	152x213 cm	1
SU	0	Back table cover	140x190 cm	1



Extremity Pack

Item no.	TB29118CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Side sheet, adhesive	75x90 cm	1
Material	Ultimate - 3 la	Extremity drape, 6 cm elastic fenestration	218x323 cm	1
PU	8	Adhesive tapes	9x50 cm	1
SU	0	Back table cover, reinforced, 50µm	140x190 cm	1
		Mayo stand cover, reinforced	80x142 cm	1




Drapes & Packs

Orthopedic Surgery

► Extremity Drapes and Packs


Extremity pack

Item no.	ES10410CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Absorbent Hand towels	20x30 cm	2
Material	ESS - 2 layers	Back table cover,	152x190 cm	2
PU	8	Adhesive tape	9x50 cm	1
SU	0	Mayo Stand Cover, absorbent reinf. 55x88 cm	80x142 cm	1
		Extremity drape, elast. Fenestration 5cm	300 x 200 cm	1
		Drape adhesive	75x75 cm	1



Extremity pack

Item no.	DYJPEEXPEC	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Impervious stockinette</td> <td>23x91 cm</td> <td>1</td>		Impervious stockinette	23x91 cm	1
Material	ECL - Spunlac	Drape sheet	135x196 cm	1
PU	5	Extremity drape, elastic fen. 6 cm	229x328 cm	1
SU	0	OP tape	10x50 cm	2
		Back table cover,	152x190 cm	1
		Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1



Extremity pack

Item no.	DYJPEEXPSM	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Impervious stockinette</td> <td>23x91 cm</td> <td>1</td>		Impervious stockinette	23x91 cm	1
Material	ADV - SMMMS	Drape sheet	135x196 cm	1
PU	5	Extremity drape, elastic fen. 6 cm	221x328 cm	1
SU	0	OP tape	10x50 cm	2
		Back table cover,	152x190 cm	1
		Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1




Drapes & Packs

Orthopedic Surgery

► Extremity Drapes and Packs


Extremity pack

Item no.	DYJPEEXPSM2	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1
Material	ADV - SMMMS	Extremity drape, 1 elastic fen. 4 cm	221x328 cm	1
PU	1	OP tape	10x50 cm	1
SU	0	Back table cover,	152x190 cm	1




Extremity pack

Item no.	DYJPEEXPSM1	Description	Dimensions	Quantity
eClass-Nr.		Drape sheet	135x196 cm	1
Material	ADV - SMMMS	Extremity drape, elastic fen. 6 cm	221x328 cm	1
PU	4	Mayo stand cover, SMS reinf. 55x88 cm	80x142 cm	1
SU	0	OP tape	10x50 cm	1
		Back table cover,	152x190 cm	1



Extremity pack

Item no.	DYJPEEXPEC3	Description	Dimensions	Quantity
eClass-Nr.		Impervious stockinette	23x91 cm	1
Material	ECL - Spunlac	Drape sheet	135x196 cm	1
PU	6	Extremity T-drape, elastic fen. 6 cm	229/307x325 cm	1
SU	0	OP tape	10x50 cm	2
		Back table cover,	152x190 cm	1
		Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1




Drapes & Packs

Orthopedic Surgery

► Extremity Drapes and Packs


Extremity pack

Item no.	DYJPEEXPSM3	Description	Dimensions	Quantity
eClass-Nr.		Impervious stockinette	23x91 cm	1
Material	ADV - SMMMS	Drape sheet	135x196 cm	1
PU	6	Extremity T-drape, elastic fen. 6 cm	221/307x325 cm	1
SU	0	Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1
		OP tape	10x50 cm	2
		Back table cover,	152x190 cm	1
		Hand towel	34x36 cm	4




Extremity Pack

Item no.	CP29032CE	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Mayo stand cover, ultra absorbent reinforcement 75x100 cm</td> <td>80x145 cm</td> <td>1000</td>		Mayo stand cover, ultra absorbent reinforcement 75x100 cm	80x145 cm	1000
Material	Ultimate - 3 la	Three quarter drape	145x193 cm	500
PU	8	Extremity drape, 6 cm elastic fenestration	218x323 cm	1
SU	0	Back table cover, 50µm	120x140 cm	1
		Back table cover, 75µm	140x190 cm	1
		Hand towels	30x39 cm	300



Extremity pack

Item no.	TB29020CE	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Drape</td> <td>90x150 cm</td> <td>1</td>		Drape	90x150 cm	1
Material	Ultimate - 3 la	Extremity drape, reinforced, elastic fenestration diameter 6cm	269/221x325 cm	1
PU	6	Back table cover	140x190 cm	1
SU	0	U-drape PE, split 20x76cm, adhesive	152x213 cm	1
		CSR wrap	86x86 cm	1




Drapes & Packs

Orthopedic Surgery

► Extremity Drapes and Packs

Extremity pack

Item no.	DYJPEEXPEC1	Description	Dimensions	Quantity
eClass-Nr.		Drape sheet	135x196 cm	1
Material	Eclipse - Spunl	Extremity drape, elastic fen. 6 cm	229x328 cm	1
PU	4	Mayo stand cover, SMS reinf. 55x88 cm	80x142 cm	1
SU	0	OP tape	10x50 cm	1
		Back table cover,	152x190 cm	1




Extremity pack

Item no.	TB29184CE	Description	Dimensions	Quantity
eClass-Nr.		Drape	90x150 cm	1
Material	Ultimate - 3 la	Drape, 3 cm diam. elastic fen.	303x300 cm	1
PU	6	CSR wrap	86x86 cm	1
SU	0	Back table cover	140x190 cm	1



Extremity pack

Item no.	TB29973CE	Description	Dimensions	Quantity
eClass-Nr.		Impervious stockinette	23x122 cm	1
Material	Ultimate - 3 la	Three quarter drape	145x193 cm	1
PU	6	Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1
SU	0	Adhesive tapes	9x50 cm	3
		Extremity drape, reinforced, elastic fenestration diameter 6cm	230x325 cm	1
		Back table cover	140x190 cm	1




Drapes & Packs

Orthopedic Surgery

► Extremity Drapes and Packs

Extremity pack

Item no.	DYJPEEXPEC2	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, SMS reinf. 55x88 cm	80x142 cm	1
Material	Eclipse - Spunl	Extremity drape, 4 cm diam. elastic fen.	229x328 cm	1
PU	5	OP tape	10x50 cm	1
SU	0	Back table cover,	152x190 cm	1




Extremity pack

Item no.	DYJPEEXPEC3G	Description	Dimensions	Quantity
eClass-Nr.		Surgical gown, L, PE reinforced		1
Material	Eclipse - Spunl	Drape sheet	135x196 cm	1
PU	5	Extremity T-drape, elastic fen. 6 cm	229/307x325 cm	1
SU	0	Surgical gown, XL, PE reinforced		1
		Back table cover,	152x190 cm	1
		OP tape	10x50 cm	2
		Impervious stockinette	23x91 cm	1
		Hand towel	34x36 cm	4



Extremity pack

Item no.	ES10412CE	Description	Dimensions	Quantity
eClass-Nr.		Drape sheet	150x150 cm	1
Material	ESS - 2 layers	Extremity drape, elast. Fenestration 6cm	300 x 230 cm	1
PU	10	Adhesive tape	9x50 cm	1
SU	555	SMS wrap		1
		Back table cover	140x190 cm	1




Drapes & Packs

Orthopedic Surgery

► Extremity Drapes and Packs


Extremity pack

Item no.	ES29962CE	Description	Dimensions	Quantity
eClass-Nr.		Absorbent Hand towels	20x30 cm	2
Material	ESS - 2 layers	Adhesive tape	9x50 cm	2
PU	6	Stockinette	36.5x72 cm	1
SU	555	Reinforced extremity drape w/ armboard covers, fen diam 6cm	221/295x335 cm	1
		Back table cover,	152x190 cm	2
		Mayo Stand Cover, absorbent reinf. 55x88 cm	80x142 cm	1



Extremity pack

Item no.	ES29973CE	Description	Dimensions	Quantity
eClass-Nr.		Adhesive tape	9x50 cm	2
Material	ESS - 2 layers	Drape adhesive	75x75 cm	1
PU	8	Extremity drape, reinforced, elastic fenestration diameter 6cm	230x325 cm	1
SU	555	Back table cover,	152x190 cm	2
		Mayo Stand Cover, absorbent reinf. 55x88 cm	80x142 cm	1
		Absorbent Hand towels	20x30 cm	2



Extremity Pack

Item no.	AG29C637B	Description	Dimensions	Quantity
eClass-Nr.		Hand towels	30x20 cm	2
Material	Eclipse - Spunl	Adhesive tape, PE	8x40 cm	2
PU	4	Drape with table extension, reinforced, elastic fenestration diam 3cm	196x371 cm	1
SU	0	Stockinette	20x96 cm	1
		OR gown, Eclipse, XLL, fabric reinforced	150 cm	2
		Back table cover, reinforced, 50µm	140x190 cm	1




Drapes & Packs

Orthopedic Surgery

► Extremity Drapes and Packs


Extremity Pack

Item no.	AG29C636B	Description	Dimensions	Quantity
eClass-Nr.		Hand towels	30x20 cm	2
Material	Eclipse - Spunl	Adhesive tape, PE	8x40 cm	2
PU	5	Extremity sheet, reinforced, elastic fenestration diameter 6cm	218x321 cm	1
SU	0	Stockinette	20x96 cm	1
		OR gown, Eclipse, XLL, fabric reinforced	150 cm	2
		Back table cover, reinforced, 50µm	140x190 cm	1




Extremity pack

Item no.	TB29185CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Extremity drape, reinforced, elastic fenestration diameter 6cm	230x325 cm	1
Material	Ultimate - 3 la	Back table cover	140x190 cm	1
PU	8	CSR wrap	86x86 cm	1
SU	0			



Extremity pack 2

Item no.	ES10411CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Adhesive tape	9x50 cm	1
Material	ESS - 2 layers	Stockinette	36.5x72 cm	1
PU	8	Extremity drape, elast. Fenestration 5cm	300 x 200 cm	1
SU	0	Absorbent Hand towels	20x30 cm	2
		Back table cover	140x190 cm	1
		Back table cover, folded, 50µm	140x190 cm	1
		Mayo Stand Cover, absorbent reinf. 55x88 cm	80x142 cm	1




Drapes & Packs

Orthopedic Surgery

► Extremity Drapes and Packs


Extremity pack,

Item no.	TB30C012CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Drape	91x150 cm	1
Material	Ultimate - 3 la	Extremity drape, elastic fenestration diameter 5cm	242x300 cm	1
PU	10	Velcro strip	2.5x13 cm	1
SU	0	Back table cover	140x190 cm	1
		CSR wrap	86x86 cm	1




Extremity pack,

Item no.	TB962CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Extremity drape, reinforced, elastic fenestration diameter 6cm	269/221x325 cm	1
Material	Ultimate - 3 la	Mayo stand cover, ultra absorbent reinforcement 75x100 cm	80x145 cm	1
PU	5	Absorbent Hand towels	20x30 cm	2
SU	0	Back table cover, folded, 50µm	140x190 cm	1
		Adhesive tapes	9x50 cm	1
		Back table cover, reinforced, 50µm	140x190 cm	1



Extremity Pack, reinforced

Item no.	ES29186CE	Description	Dimensions	Quantity
eClass-Nr.	ESS - 2 layers	Drape sheet	150x150 cm	1
Material	10	Extremity drape, reinforced, elastic fenestration diameter 6cm	230x325 cm	1
PU	555	SMS wrap		1
SU		Back table cover, reinforced, 50µm	140x190 cm	1




Drapes & Packs

Orthopedic Surgery

► Extremity Drapes and Packs


Extremity T-drape

Item no.	DYNJPE8003SM	Description	Dimensions	Quantity
eClass-Nr.		Extremity T-drape, elastic fen. 6 cm	221/307x325 cm	1
Material	ADV - SMMMS	Wrap	91x91 cm	1
PU	6			
SU	0			




Extremity T-drape

Item no.	DYNJPE8003	Description	Dimensions	Quantity
eClass-Nr.		Extremity T-drape, elastic fen. 6 cm	229/307x325 cm	1
Material	ECL - Spunlac	Wrap	91x91 cm	1
PU	6			
SU	0			



Limb Pack

Item no.	TB29C024CEA	Description	Dimensions	Quantity
eClass-Nr.		Extremity drape, reinforced, elastic fenestration diameter 6cm	269/221x325 cm	1
Material	Ultimate - 3 la	Back table cover, folded, 50µm	140x190 cm	1
PU	6	Adhesive tape, PE	8x40 cm	2
SU	0	Back table cover	140x190 cm	1
		CSR wrap	86x86 cm	1




Drapes & Packs

Orthopedic Surgery

► Extremity Drapes and Packs


Lower extremity pack

Item no.	TB29041CE	Description	Dimensions	Quantity
eClass-Nr.		Adhesive tape, PE	8x40 cm	4
Material	Ultimate - 3 la	Extremity drape, reinforced, elastic fenestration diameter 6cm	230x325 cm	1
PU	6	Back table cover, reinforced, 50µm	125x230 cm	1
SU	0	Stockinette	20x96 cm	1
		CSR wrap	86x86 cm	1



Lower extremity pack, reinforced

Item no.	ES29024CE	Description	Dimensions	Quantity
eClass-Nr.		Adhesive tape	9x50 cm	2
Material	ESS - 2 layers	Drape sheet	150x150 cm	1
PU	8	Back table cover, reinforced, 50µm	140x190 cm	1
SU	555	Reinforced extremity drape w/ armboard covers, fen diam 6cm	221/295x335 cm	1
		SMS wrap		1



Drapes & Packs

Orthopedic Surgery

► Hand Drapes and Packs


Hand drape

Item no.	TB29419CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Drape with table extension, reinforced, elastic fenestration diam 3 cm	196x371 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	9			
SU	0			




Hand drape

Item no.	DYNJPE8005SM	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Wrap</td> <td>91x91 cm</td> <td>1</td>		Wrap	91x91 cm	1
Material	ADV - SMMMS	Hand drape, elastic fen. 5 cm diam.	196/371x292 cm	1
PU	18			
SU	0			



Hand drape

Item no.	SM29419CE	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Hand drape</td> <td>162.5/373x196cm</td> <td>1</td>		Hand drape	162.5/373x196cm	1
Material	ADV - SMMMS	SMS wrap		1
PU	12			
SU	0			




Drapes & Packs

Orthopedic Surgery

► Hand Drapes and Packs


Hand Drape, fen 3 cm

Item no.	29419CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Drape with table extension, reinforced, elastic fenestration diam 3cm	196x371 cm	1
Material	ECL - Spunlac	CSR wrap	86x86 cm	1
PU	10			
SU	0			




Hand pack

Item no.	TB29095CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Back table cover, folded, 50µm	140x190 cm	1
Material	Ultimate - 3 la	Drape with table extension, reinforced, elastic fenestration diam 3 cm	196x371 cm	1
PU	6	Back table cover	140x190 cm	1
SU	0	CSR wrap	86x86 cm	1



Hand pack

Item no.	TB29409CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Hand towels	30x20 cm	1
Material	Ultimate - 3 la	Mayo stand cover, reinforced	80x142 cm	1
PU	7	Drape with table extension, reinforced, elastic fenestration diam 3 cm	196x371 cm	1
SU	0	Back table cover	140x190 cm	1
		CSR wrap	112x112 cm	1




Drapes & Packs

Orthopedic Surgery

► Hand Drapes and Packs


Hand pack

Item no.	TB39418CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Drape with table extension, reinforced, elastic fenestration diam 3 cm	196x371 cm	1
Material	Ultimate - 3 la	Back table cover	140x190 cm	1
PU	7	CSR wrap	86x86 cm	1
SU	0			



Hand Pack

Item no.	RGBS504	Description	Dimensions	Quantity
eClass-Nr.		Drape with table extension, reinforced, elastic fenestration diam 3 cm	196x371 cm	1
Material	Ultimate - 3 la	Back table cover, folded, 50µm	140x190 cm	1
PU	6	Back table cover	140x190 cm	1
SU	0			



Hand pack

Item no.	DYJPEHAPECG	Description	Dimensions	Quantity
eClass-Nr.		Drape sheet	135x196 cm	1
Material	Eclipse - Spunl	Surgical gown, L		1
PU	6	Surgical gown, XL		1
SU	0	Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1
		OP tape	10x50 cm	1
		Back table cover,	152x190 cm	1
		Hand towel	34x36 cm	4




Drapes & Packs

Orthopedic Surgery

► Hand Drapes and Packs


Hand Pack

Item no.	SM39418CE	Description	Dimensions	Quantity
eClass-Nr.		Hand drape	162.5/373x196cm	1
Material	ADV - SMMMS	Back table cover	140x190 cm	1
PU	10	SMS wrap		1
SU	555			




Hand Pack

Item no.	SM29095CEA	Description	Dimensions	Quantity
eClass-Nr.		Hand drape	162.5/373x196cm	1
Material	ADV - SMMMS	Drape sheet	150x150 cm	1
PU	8	SMS wrap		1
SU	555	Back table cover, reinforced, 50µm	140x190 cm	1



Hand pack

Item no.	SM29409CEA	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
Material	ADV - SMMMS	Drape sheet	150x150 cm	1
PU	8	Hand towels	30x39 cm	2
SU	555	Back table cover,	152x190 cm	1
		Hand drape	162.5/373x196cm	1



Drapes & Packs

Orthopedic Surgery

► Hip Drapes and Packs


Hip drape

Item no.	TB29430CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Hip drape, reinforced, elastic fenestration diameter 20cm	221x326 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	10			
SU	0			




Hip drape

Item no.	DYNJPE8201SM	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Wrap</td> <td>91x91 cm</td> <td>1</td>		Wrap	91x91 cm	1
Material	ADV - SMMMS	Hip drape, adh. split 20x75 cm, elastic fen. w/ adh. 16x20 cm	221/285x352 cm	1
PU	5			
SU	0			



Hip drape

Item no.	DYNJPE8201	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Wrap</td> <td>91x91 cm</td> <td>1</td>		Wrap	91x91 cm	1
Material	ECL - Spunlac	Hip drape, elastic pear shaped fen. 16x20 cm, 2 pouches	226/284x348 cm	1
PU	5			
SU	0			



Drapes & Packs

Orthopedic Surgery

► Hip Drapes and Packs

Hip drape

Item no.	DYNJPE8202	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	ECL - Spunlac	Hip drape, 23 cm diam. elastic fen., 1 instrument pad	226x325 cm	1
PU	12			
SU	0			




Hip drape for anterior approach

Item no.	ES29431CEA	Description	Dimensions	Quantity
eClass-Nr.		Hip drape with stockinette	360x400 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	5			
SU	555			



Hip Drape with Pouches

Item no.	TB29429CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Hip drape; w/adh, w/fen, w/pouches	217/290x350 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	6			
SU	0			



Drapes & Packs

Orthopedic Surgery

► Hip Drapes and Packs

Hip Pack

Item no.	TB29187CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Impervious surgical gown poly reinforced, XL, 130 cm	130 cm	2
Material	Ultimate - 3 la	Impervious stockinette	30x122 cm	1
PU	4	Suture bag, adh.	16x31 cm	1
SU	0	Drape Hip, fen w/ kraton	325 x 221 cm	1
		Back table cover, reinforced, 50µm	125x230 cm	1
		U-drape PE, split 20x76cm, adhesive	152x213 cm	1
		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
		Hand towels	30x39 cm	4



Hip pack

Item no.	DYJPEHIPPSM	Description	Dimensions	Quantity
eClass-Nr.	ADV - SMMMS	Hip drape, adh. split 20x75 cm, elastic fen. w/ adh. 16x20 cm	221/285x352 cm	1
Material	ADV - SMMMS	Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1
PU	6	Drape sheet	135x196 cm	1
SU	0	Impervious stockinette	30x122 cm	1
		OP tape	10x50 cm	2
		U drape, adh. split 20x75 cm	152x213 cm	1
		Hand towel	34x36 cm	4
		Back table cover,	152x190 cm	1



Hip Pack

Item no.	ES29428CE	Description	Dimensions	Quantity
eClass-Nr. <td>ESS - 2 layers</td> <td>Absorbent Hand towels</td> <td>20x30 cm</td> <td>4</td>	ESS - 2 layers	Absorbent Hand towels	20x30 cm	4
Material	ESS - 2 layers	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
PU	6	Impervious stockinette	30x122 cm	1
SU	555	Adhesive tape	9x50 cm	2
		U-drape PE, split 20x76cm, adhesive	152x213 cm	1
		Hip drape, reinforcement, pouch	217/290x350 cm	1
		Back table cover, reinforced, 50µm	140x190 cm	1
		SMS wrap		1




Drapes & Packs

Orthopedic Surgery

► Hip Drapes and Packs


Luxation pouch

Item no.	CP8691CE	Description	Dimensions	Quantity
eClass-Nr.		Luxation pouch, 77x110cm, with adhesive	77x100 cm	1
Material	PE	CSR wrap	60x60 cm	1
PU	80			
SU	0			




Traumatology pack

Item no.	CP29054CEA	Description	Dimensions	Quantity
eClass-Nr.		Side sheet, adhesive	75x90 cm	6
Material	PE	Absorbent Towel	39x58 cm	4
PU	4	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
SU	0	Adh extendable, turquoise	9x50 cm	2
		Pouch, Instmt, double, PE	33x41 cm	1
		Vertical Isolation Drape, PE, 330x245cm with fen	246x330cm	1
		back table cover, 90 µm	203x240 cm	1
		CSR wrap	112x112 cm	1



Vertical Isolation Drape

Item no.	5803CE	Description	Dimensions	Quantity
eClass-Nr.		vertical isolation drape	330x246 cm	1
Material	ESS - 2 layers	CSR wrap	86x86 cm	1
PU	10			
SU	0			




Drapes & Packs

Orthopedic Surgery

► Hip Drapes and Packs


Vertical isolation drape, PE

Item no.	5801CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Vertical isolation drape, incise 78x26 cm, inst. pouch 25.5x28 cm	270x330 cm	1
Material	PE	CSR wrap	85x85 cm	1
PU	16			
SU	0			




Vertical Isolation drape, Set

Item no.	TB29119CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Hand towels	30x20 cm	1
Material	PE	Mayo stand cover, ultra absorbent reinforcement	80x145 cm	1
PU	10	Vertical isolation drape, incise 78x26 cm, inst. pouch 25.5x28 cm	270x330 cm	1
SU	0	Back table cover, reinforced, 50µm	140x190 cm	1



Vertical pack

Item no.	CP5802CE	Description	Dimensions	Quantity
eClass-Nr.		Hand towels	30x20 cm	2
Material	PE	Adhesive tape, PE	8x40 cm	1
PU	12	Vertical isolation drape, incise 78x26 cm, inst. pouch 25.5x28 cm	270x330 cm	1
SU	0	Back table cover, 75µm	140x190 cm	1
		CSR wrap	86x86 cm	1



Drapes & Packs

Orthopedic Surgery

► Shoulder Drapes and Packs

Beach chair shoulder drape,

Item no.	TB29369CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Drape with pouch, reinforced, pear-shaped fenestration, adhesive, velcro	262x411 cm	1
Material	Ultimate - 3 la	CSR wrap	112x112 cm	1
PU	6			
SU	0			




Beach chair shoulder drape,

Item no.	DYNJPE8412	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	Eclipse - Spunl	Shoulder drape, elastic and adh. fen. 9x10 cm	404x264 cm	1
PU	5			
SU	0			



Shoulder Drape Pack

Item no.	TB29365CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Shoulder split drape, pe-reinforced, split 25x64cm, adhesive	196 x 295 cm	1
Material	Ultimate - 3 la	Shoulder Drape		1
PU	9	CSR wrap	76x76 cm	1
SU	0			




Drapes & Packs

Orthopedic Surgery

► Shoulder Drapes and Packs

Shoulder drape pack

Item no.	DYNJPE8401SM	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	ADV - SMMMS	Body U drape for shoulder, adh. split 25x63 cm	196x290 cm	1
PU	10	Shoulder drape, shoulder shape adh. fen. w/ split 28x35 cm	254x168 cm	1
SU	0			




Shoulder Drape with Pouch for Sitting Position

Item no.	TB29370CE	Description	Dimensions	Quantity
eClass-Nr.		Shoulder drape, top 127 cm in clear film, 8 cm diameter elastic	411x262 cm	1
Material	Ultimate - 3 la	SMS wrap		1
PU	8			
SU	555			



Shoulder drape,

Item no.	TB29367CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Shoulder drape with pouch, pear-shaped split with incise	160x254 cm	1
Material	Ultimate - 3 la	Shoulder split drape, pe-reinforced, split 25x64cm, adhesive	196 x 295 cm	1
PU	8	CSR wrap	76x76 cm	1
SU	0			




Drapes & Packs

Orthopedic Surgery

► Shoulder Drapes and Packs


Shoulder pack

Item no.	DYJPESHPEC3	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1
Material	ECL - Spunlac	Shoulder drape, elastic and adh. fen. 9x10 cm	404x264 cm	1
PU	6	U drape, adh. split 15x61 cm	137x193 cm	1
SU	0	Back table cover,	152x190 cm	1




Shoulder pack

Item no.	DYJPESHPSM3	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1
Material	ADV - SMMMS	Beach chair shoulder drape, 9x10 cm elastic fen. w/ adh. strip	404x268 cm	1
PU	6	Hand towel	34x36 cm	4
SU	0	U drape, adh. split 15x61 cm	137x193 cm	1
		Back table cover,	152x190 cm	1



Shoulder pack

Item no.	DYJPESHPSM1	Description	Dimensions	Quantity
eClass-Nr.		Shoulder U drape, blue PE part w/ 7x20 cm pear shaped split	254x168 cm	1
Material	ADV - SMMMS	Body U drape for shoulder, adh. split 25x63 cm	196x290 cm	1
PU	6	Back table cover,	152x190 cm	1
SU	0	Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1




Drapes & Packs

Orthopedic Surgery

► Shoulder Drapes and Packs


Shoulder pack

Item no.	DYJPESHPSM2	Description	Dimensions	Quantity
eClass-Nr.		Shoulder drape, shoulder shape adh. fen. w/ split 28x35 cm	254x168 cm	1
Material	ADV - SMMMS	Body U drape for shoulder, adh. split 25x63 cm	196x290 cm	1
PU	6	Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1
SU	0	Back table cover,	152x190 cm	1
		Hand towel	34x36 cm	4




Shoulder pack

Item no.	DYJPESHPEC2	Description	Dimensions	Quantity
eClass-Nr.		Shoulder drape, shoulder shape adh. fen. w/ split and incise film, pouch	254x168 cm	1
Material	Eclipse - Spunl	Body U drape for shoulder, adh. split 25x63 cm	196x290 cm	1
PU	6	Back table cover,	152x190 cm	1
SU	0	Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1



Shoulder Pack with Pouch

Item no.	CP29122CE	Description	Dimensions	Quantity
eClass-Nr.		Drape, w/ adhesive	150x196 cm	1
Material	Ultimate - 3 la	Instrument pouch, PE, adhesive	30,5x29 cm	1
PU	4	Back table cover, folded, 50µm	140x190 cm	1
SU	0	Velcro strip	2.5x13 cm	2
		Drape with pouch, reinforced, pear-shaped fenestration, adhesive, velcro	262x411 cm	1
		Back table cover	140x190 cm	1
		CSR wrap	60x60 cm	1



Drapes & Packs

Orthopedic Surgery

► Shoulder Drapes and Packs


Shoulder pack,

Item no.	9183CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	OR gown, Eclipse, L	120 cm	1
Material	ECL - Spunlac	MAYO Stand cover, S, 50µ, absorbent reinf. 46x61 cm	59x141 cm	1
PU	3	OR gown, Ultimate, XLL, impervious, breathable, raglan sleeves	150 cm	1
SU	0	Suture bag		1
		Shoulder Split Drape	196x295 cm	1
		Shoulder Drape w/ Pouch	160x254 cm	1
		Back table cover, reinforced, 50µ	125x230 cm	1
		CSR wrap	60x60 cm	1



Shoulder pack, beach chair

Item no.	TB29371CE	Description	Dimensions	Quantity
eClass-Nr.		Side sheet, adhesive	75x90 cm	3
Material	Ultimate - 3 la	Impervious stockinette	20x97 cm	1
PU	4	Drape with pouch, reinforced, pear-shaped fenestration, adhesive, velcro	262x411 cm	1
SU	0	Back table cover, 90µm	160x240 cm	1
		Mayo stand cover, ultra absorbent reinforcement 75x100 cm	80x145 cm	1
		Adhesive tape	9x50 cm	2



Shoulder split drape

Item no.	DYNJPE8405	Description	Dimensions	Quantity
eClass-Nr.		Wrap	80x80 cm	1
Material	ECL - Spunlac	Shoulder U drape, pear shape adh. fen. 20x23 cm, pouch	254x178 cm	1
PU	8			
SU	0			




Drapes & Packs

Orthopedic Surgery

► Shoulder Drapes and Packs

Shoulder/Arthroscopy pack

Item no.	TB29368CE	Description	Dimensions	Quantity
eClass-Nr.		Adhesive tape	9x50 cm	2
Material	Ultimate - 3 la	Shoulder drape with pouch, pear-shaped split with incise	160x254 cm	1
PU	6	Shoulder split drape, pe-reinforced, split 25x64cm, adhesive	196 x 295 cm	1
SU	555	SMS wrap		1
		Back table cover, reinforced, 50µm	140x190 cm	1
		Absorbent Hand towels	20x30 cm	4




Drapes & Packs

Orthopedic Surgery

► Split Drapes and Packs


Double split sheet

Item no.	DYJPEDSDSM	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	ADV - SMMMS	Double split drape, 2 adh. split 6x90 cm	240x286 cm	1
PU	8			
SU	0			



Double split sheet

Item no.	DYJPEDSDEC	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	Eclipse - Spunt	Double split drape, 2 adh. split 6x90 cm	244x287 cm	1
PU	8			
SU	0			



Double split sheet

Item no.	TB29038CE	Description	Dimensions	Quantity
eClass-Nr.		Double split sheet, 11x94 cm adh. splits, 76x90 cm reinforc., 35x7..	240x290 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	10			
SU	0			




Drapes & Packs

Orthopedic Surgery

► Split Drapes and Packs

Extremity pack,


Item no. TB29116CE
eClass-Nr. 34-35-01-02
Material Ultimate - 3 la
PU 4
SU 0



Description	Dimensions	Quantity
Absorbent Towel	39x58 cm	1
Adhesive tape	9x50 cm	2
Drape	98x150 cm	1
U-drape PE, split 20x76cm, adhesive	152x213 cm	1
Split sheet, reinforced, split 6x94cm, adhesive	196x274 cm	1
Head sheet, reinforced, adhesive	183x254 cm	1
Back table cover	140x190 cm	1
CSR wrap	86x86 cm	1

Extremity pack,


Item no. TB29117CE
eClass-Nr. 34-35-01-02
Material Ultimate - 3 la
PU 4
SU 0



Description	Dimensions	Quantity
Back table cover, reinforced, 50µm	140x190 cm	1
Adhesive tape	9x50 cm	2
Drape	98x150 cm	1
U-drape PE, split 20x76cm, adhesive	152x213 cm	1
Split sheet, reinforced, split 6x94cm, adhesive	196x274 cm	1
Head sheet, reinforced, adhesive	183x254 cm	1
Back table cover	140x190 cm	1
CSR wrap	86x86 cm	1

Hip Pack

Item no. CP29031CE
eClass-Nr. Ultimate - 3 la
Material Ultimate - 3 la
PU 4
SU 0



Description	Dimensions	Quantity
Split drape, adhesive 110cm	240x304cm	1
Instrument pouch in clear plastic with adhesive	33x41cm	1
Three quarter drape	145x193 cm	1
Mayo stand cover, ultra absorbent reinforcement 75x100 cm	80x145 cm	1
drape sheet	75x90cm	1
Hand towels	30x39 cm	4
Adhesive tape	9x50 cm	3
Drape with adhesive	300x221 cm	1
Back table cover, 50µm	120x140 cm	1
Back table cover, 75µm	140x190 cm	1


Drapes & Packs

Orthopedic Surgery

► Split Drapes and Packs

Hip pack

Item no.	TB29509CE	Description	Dimensions	Quantity
eClass-Nr.		Impervious stockinette	30x122 cm	1
Material	Ultimate - 3 la	Split drape reinforced, adhesive split 15x116cm	196x305 cm	1
PU	6	Head drape, reinforcement, adhesive, 1 cord holding tab	254x183 cm	1
SU	0	U-drape PE, split 20x76cm, adhesive	152x213 cm	1
		Back table cover	140x190 cm	1




Orthopaedic pack

Item no.	TB29181CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Impervious stockinette	30x122 cm	1
Material	Ultimate - 3 la	Impervious surgical gown poly reinforced, XL, 130 cm	130 cm	2
PU	4	Split drape, adhesive	196x254 cm	1
SU	0	Head sheet, reinforced, adhesive	251x135 cm	1
		Back table cover, reinforced, 50µm	125x230 cm	1
		Hand towels	30x39 cm	4
		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
		Suture bag	25x33 cm	1
		U-drape PE, split 20x76cm, adhesive	152x213 cm	1



Orthopaedic pack

Item no.	DYJPEOPPSM	Description	Dimensions	Quantity
eClass-Nr.		Head drape, adh.	300x183 cm	1
Material	ADV - SMMMS	Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1
PU	5	Drape sheet	135x196 cm	1
SU	0	Impervious stockinette	31x122 cm	1
		Hand towel	34x36 cm	4
		U drape, adh. split 8x107 cm	220x300 cm	1
		U drape, adh. split 20x75 cm	152x213 cm	1
		Back table cover,	152x190 cm	1
		OP tape	10x50 cm	2



Drapes & Packs

Orthopedic Surgery

► Split Drapes and Packs

Orthopaedic pack

Item no.	DYJPEOPPEC2G	Description	Dimensions	Quantity
eClass-Nr.		Split drape, adh. split 8x107 cm	231x300 cm	1
Material	ADV - SMMMS	Impervious stockinette	31x122 cm	1
PU	4	Drape sheet	135x196 cm	1
SU	0	Surgical gown, XL		1
		Surgical gown, L		1
		OP tape	10x50 cm	2
		U drape, adh. split 20x75 cm	152x213 cm	1
		Hand towel	34x36 cm	4
		Head drape, adh.	300x178 cm	1
		Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1




Orthopaedic pack

Item no.	DYJPEOPPSM2G	Description	Dimensions	Quantity
eClass-Nr.		Impervious stockinette	31x122 cm	1
Material	ADV - SMMMS	Head drape, adh.	300x183 cm	1
PU	4	Drape sheet	135x196 cm	1
SU	0	Surgical gown, XL, w/ 2 hand towels 34x36 cm		1
		Surgical gown, L, w/ 2 hand towels 34x36 cm		1
		U drape, adh. split 8x107 cm	220x300 cm	1
		U drape, adh. split 20x75 cm	152x213 cm	1
		OP tape	10x50 cm	2
		Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1
		Back table cover,	152x190 cm	1



Orthopaedic pack

Item no.	DYJPEOPPSM3	Description	Dimensions	Quantity
eClass-Nr.		Hip drape, 99 cm adh. slit, pear shape adh. fen. 17x32 cm	196x285 cm	1
Material	ADV - SMMMS	Head drape, adh.	300x183 cm	1
PU	6	Drape sheet	135x196 cm	1
SU	0	Impervious stockinette	23x91 cm	1
		Back table cover,	152x190 cm	1
		OP tape	10x50 cm	1
		Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1
		Hand towel	34x36 cm	4




Drapes & Packs

Orthopedic Surgery

► Split Drapes and Packs


Orthopaedic pack

Item no.	DYJPEOPPSM4	Description	Dimensions	Quantity
eClass-Nr.		Head drape, adh.	300x183 cm	1
Material	ADV - SMMMS	Drape sheet	135x196 cm	1
PU	6	Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1
SU	0	U drape, adh. split 6x94 cm	196x274 cm	1
		Back table cover,	152x190 cm	1




Orthopaedic pack

Item no.	DYJPEOPPSM6	Description	Dimensions	Quantity
eClass-Nr.		Hand towel	34x36 cm	4
Material	ADV - SMMMS	Adhesive drape	229x290 cm	1
PU	8	Drape sheet	89x102 cm	1
SU	0	U drape, adh. split 15x102 cm	221x274 cm	1
		OP tape	10x50 cm	2
		Back table cover,	152x190 cm	1



Orthopaedic pack

Item no.	DYJPEOPPEC4	Description	Dimensions	Quantity
eClass-Nr.		Split drape, adh. split 6x94 cm	196x274 cm	1
Material	Eclipse - Spunl	Drape sheet	135x196 cm	1
PU	6	Head drape, adh.	300x178 cm	1
SU	0	Back table cover,	152x190 cm	1
		Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1




Drapes & Packs

Orthopedic Surgery

► Split Drapes and Packs


Orthopaedic pack

Item no.	DYJPEOPPEC	Description	Dimensions	Quantity
eClass-Nr.		Impervious stockinette	31x122 cm	1
Material	Eclipse - Spunl	Head drape, adh.	300x178 cm	1
PU	5	Drape sheet	135x196 cm	1
SU	0	Split drape, adh. split 8x107 cm	231x300 cm	1
		Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1
		OP tape	10x50 cm	2
		U drape, adh. split 20x75 cm	152x213 cm	1
		Hand towel	34x36 cm	2
		Back table cover,	152x190 cm	1




Orthopaedic pack

Item no.	TB9650CE	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1
Material	Ultimate - 3 la	Back table cover, folded, 50µm	140x190 cm	1
PU	6	Absorbent Hand towels	20x30 cm	1
SU	0	Head sheet, reinforced, adhesive	254x183 cm	1
		Split sheet, reinforced, split 6x94cm, adhesive	196 x 275 cm	1
		Back table cover	140x190 cm	1



Orthopaedic pack

Item no.	ES29630CE	Description	Dimensions	Quantity
eClass-Nr.		Plastic split sheet, adhesive split	152 x 213 cm	1
Material	ESS - 2 layers	bi-laminate head drape	260x170 cm	1
PU	6	Impervious stockinette	30x122 cm	1
SU	0	drape sheet	100x100 cm	1
		Split drape, adhesive 110cm	240x304cm	1
		Back table cover, reinforced, 50µm	140x190 cm	1
		SMS wrap		1
		Hand towels	30x39 cm	2
		Adhesive tape	9x50 cm	2
		Mayo stand cover, ultra absorbent reinforcement 75x100 cm	80x145 cm	1




Drapes & Packs

Orthopedic Surgery

► Split Drapes and Packs


Orthopaedic Pack

Item no.	ES29641CE	Description	Dimensions	Quantity
eClass-Nr.		Split drape reinforced, adhesive split 15x117cm	230x305 cm	1
Material	ESS - 2 layers	Drape with adhesive	300x221 cm	1
PU	6	Stockinette	30x120 cm	1
SU	555	Foot sheet, adhesive	190x190 cm	1
		Back table cover, reinforced, 50µm	140x190 cm	1
		Mayo stand cover, ultra absorbent reinforcement 75x100 cm	80x145 cm	1
		Adhesive tape, PE	8x40 cm	5




Orthopaedic pack 2

Item no.	TB29640CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Impervious stockinette	30x122 cm	1
Material	Ultimate - 3 la	Absorbent Hand towels	20x30 cm	2
PU	5	Mayo stand cover, ultra absorbent reinforcement 75x100 cm	80x145 cm	1
SU	0	Back table cover, folded, 50µm	140x190 cm	1
		Utility drape	98x98 cm	1
		Adhesive tapes	9x50 cm	2
		Split drape reinforced, adhesive split 15x116cm	196x305 cm	1
		Head drape	254x183 cm	1
		Back table cover	140x190 cm	1



Orthopaedic pack I

Item no.	TB29630CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Absorbent Hand towels	20x30 cm	2
Material	Ultimate - 3 la	Impervious stockinette	30x122 cm	1
PU	4	Mayo stand cover, ultra absorbent reinforcement	80x145 cm	1
SU	0	U-drape PE, split 20x76cm, adhesive	152x213 cm	1
		Back table cover, reinforced, 50µm	140x190 cm	1
		Utility drape	98x98 cm	1
		Adhesive tapes	9x50 cm	2
		Split drape, adhesive 110cm	240x304cm	1
		Head drape	254x183 cm	1




Drapes & Packs

Orthopedic Surgery

► Split Drapes and Packs


Orthopaedic pack with gowns

Item no.	DG29118C	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Impervious stockinette	30x122 cm	1
Material	ECL - Spunlac	Head sheet, reinforced, adhesive	183x254 cm	1
PU	3	Split Sheet, reinforced, split 6 x 94 cm adh	196x274 cm	1
SU	0	drape sheet	75x90cm	1
		Drape, adhesive	49x65 cm	4
		Hand towels	30x39 cm	2
		OR gown, Eclipse, XLL, fabric reinforced	150 cm	1
		Hand towels	30x20 cm	1
		Adhesive tape, PE	8x40 cm	2
		Mayo stand cover, reinforced	80x142 cm	1
		U-drape PE, split 20x76cm, adhesive	152x213 cm	1
		Back table cover	140x190 cm	1
		CSR wrap	95x95 cm	1




Orthopaedic pack,

Item no.	29214FRB	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Instrument pouch, PE, adhesive	30,5x29 cm	1
Material	ECL - Spunlac	Split drape, reinforced, adhesive	241x304 cm	1
PU	4	Drape	98x98 cm	1
SU	0	Anaesthetic drape, fenestrated	196x300 cm	1
		Back table cover, reinforced, 50µm	140x190 cm	1
		Hand towels	30x39 cm	2
		Adhesive tape, PE	8x40 cm	6
		U-drape PE, split 20x76cm, adhesive	152x213 cm	1
		Back table cover, reinforced, 50µm	140x190 cm	1
		CSR wrap	95x95 cm	1



Orthopaedic pack,

Item no.	TB29215CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Absorbent Hand towels	20x30 cm	2
Material	Ultimate - 3 la	Adhesive tape	9x50 cm	4
PU	5	Drape,	98x98 cm	1
SU	0	U-drape PE, split 20x76cm, adhesive	152x213 cm	1
		Split sheet, reinforced, split 7x110cm, adhesive 110cm	240 x 304 cm	1
		Head sheet, reinforced, adhesive	254x183 cm	1
		Back table cover	140x190 cm	1
		CSR wrap	112x112 cm	1




Drapes & Packs

Orthopedic Surgery

► Split Drapes and Packs


Orthopaedic split drape

Item no.	DYNJPE8301	Description	Dimensions	Quantity
eClass-Nr.		Wrap	76x76 cm	1
Material	ECL - Spunlac	Orthopaedic U drape, adh. split 15x102 cm	196x274 cm	1
PU	10			
SU	0			




Orthopaedic Split Pack, reinforced

Item no.	ES29116CE	Description	Dimensions	Quantity
eClass-Nr.		Head drape, adhesive	260x182 cm	1
Material	ESS - 2 layers	Adhesive tape	9x50 cm	2
PU	4	Drape sheet	150x150 cm	1
SU	555	Split drape, reinforced, 8.5x107cm adhesive split	220x300 cm	1
		SMS wrap		1
		U-drape PE, split 20x76cm, adhesive	152x213 cm	1
		Absorbent Hand towels	20x30 cm	2
		Back table cover, reinforced, 50µm	140x190 cm	1



Orthopaedic Split Pack, reinforced

Item no.	ES29117CE	Description	Dimensions	Quantity
eClass-Nr.		Head drape, adhesive	260x182 cm	1
Material	ESS - 2 layers	Adhesive tape	9x50 cm	2
PU	4	Drape sheet	150x150 cm	1
SU	555	Split drape, reinforced, 8.5x107cm adhesive split	220x300 cm	1
		Absorbent Hand towels	20x30 cm	2
		SMS wrap		1
		U-drape PE, split 20x76cm, adhesive	152x213 cm	1
		Back table cover, reinforced, 50µm	140x190 cm	1
		Back table cover, folded, 50µm	140x190 cm	1




Drapes & Packs

Orthopedic Surgery

► Split Drapes and Packs


Split drape

Item no.	DYJPESPSSM1	Description	Dimensions	Quantity
eClass-Nr.		Wrap	76x76 cm	1
Material	ADV - SMMMS	Split drape, adh. split 6x60 cm (adh. width 2.5cm)	135x178 cm	1
PU	20			
SU	0			




Split drape

Item no.	DYJPESPSSM3	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	ADV - SMMMS	U drape, adh. split 8x107 cm	220x300 cm	1
PU	6			
SU	0			



Split drape

Item no.	DYJPESPSEC3	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	ECL - Spunlac	Split drape, adh. split 8x107 cm	231x300 cm	1
PU	6			
SU	0			




Drapes & Packs

Orthopedic Surgery

► Split Drapes and Packs


Split drape

Item no.	DYNJPE2498	Description	Dimensions	Quantity
eClass-Nr.		SMS wrap		1
Material	ECL - Spunlac	U drape, adh. split 6x94 cm	196x274 cm	1
PU	10			
SU	0			




Split pack

Item no.	DYJPEUSPSM5	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
Material	ADV - SMMMS	Head drape, adh. split 61x48 cm	269x182 cm	1
PU	6	Suture bag w/ adh.	16x31 cm	1
SU	0	U drape, adh. split 8x107 cm	220x300 cm	1
		Back table cover	152x229 cm	1
		Hand towel	34x36 cm	4



Split pack

Item no.	DYJPEUSPEC3G	Description	Dimensions	Quantity
eClass-Nr.		Head drape, adh. split 61x48 cm	269x150 cm	1
Material	ECL - Spunlac	Surgical gown, L		1
PU	4	Surgical gown, XL		1
SU	0	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
		U drape, adh. split 23x76 cm, instrument pad 36x28 cm	196x295 cm	1
		Back table cover,	152x190 cm	1
		Hand towel	34x36 cm	4




Drapes & Packs

Orthopedic Surgery

► Split Drapes and Packs


Split pack

Item no.	DYJPEUSPEC4	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, SMS reinf. 55x88 cm	80x142 cm	1
Material	Eclipse - Spunl	Top drape, adh.	259x135 cm	1
PU	5	U drape, adh. split 6x94 cm	196x274 cm	1
SU	0	Back table cover,	152x190 cm	1




Split pack

Item no.	ES29620CE	Description	Dimensions	Quantity
eClass-Nr.		Head drape, adhesive	160x260 cm	1
Material	ESS - 2 layers	bi-laminate side drape	110x95 cm	1
PU	8	Split drape reinforced, adhesive split 15x117cm	196x305 cm	1
SU	0	Adhesive tape	9x50 cm	2
		Back table cover, reinforced, 50µm	140x190 cm	1
		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
		Absorbent Hand towels	20x30 cm	4



Split pack large

Item no.	TB29690CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Impervious stockinette	30x122 cm	1
Material	Ultimate - 3 la	Head drape, Xlarge	340x196cm	1
PU	4	Under buttocks drape with cuff 24cm, reinforced	96x126 cm	1
SU	0	drape sheet	75x90cm	1
		Split drape, adhesive 110cm	240x304cm	1
		Back table cover, 90µm	160x240 cm	1
		U-drape PE, split 20x76cm, adhesive	152x213 cm	1
		Mayo stand cover, ultra absorbent reinforcement	80x145 cm	1
		Adhesive tapes	9x50 cm	1
		Absorbent Hand towels	20x30 cm	2




Drapes & Packs

Orthopedic Surgery

► Split Drapes and Packs


Split sheet

Item no.	DYJPESPSSM2	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	ADV - SMMMS	U drape, adh. split 6x94 cm	196x274 cm	1
PU	10			
SU	0			




Split sheet

Item no.	ES29425CE	Description	Dimensions	Quantity
eClass-Nr.		SMS wrap		1
Material	ESS - 2 layers	U-drape, split 65x7cm, reinforced, adhesive	150 x 200 cm	1
PU	42			
SU	555			



Split sheet

Item no.	ES10438CE	Description	Dimensions	Quantity
eClass-Nr.		Split sheet, adh. split 20x102 cm	200 x 260 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	28			
SU	555			




Drapes & Packs

Orthopedic Surgery

► Split Drapes and Packs

Split sheet

Item no.	ES9424CE	Description	Dimensions	Quantity
eClass-Nr.		Split drape, reinforced, 8.5x107cm adhesive split	220x300 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	14			
SU	0			



Split sheet drape,

Item no.	TB29423CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Split sheet, reinforced, split 6x94cm, adhesive	196 x 275 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	12			
SU	0			



Split sheet drape,

Item no.	TB29434CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Split sheet, reinforced, split 3x72cm, small adhesive	147x196 cm	1
Material	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU	24			
SU	0			



Drapes & Packs

Orthopedic Surgery

► Split Drapes and Packs

Split sheet drape,

Item no.	TB29436CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Split sheet, reinforced, split 6x94cm, small adhesive	196 x 274 cm	1
Material	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU	12			
SU	0			




Split sheet drape,

Item no.	TB9437CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Split sheet, reinforced, split 7x110cm, adhesive 110cm	240 x 304 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	10			
SU	0			



Thoracic pack

Item no.	CP29007CE	Description	Dimensions	Quantity
eClass-Nr.	Eclipse - Spunt	Split drape, reinforced, adhesive	241x304 cm	1
Material	4	Instrument pouch, PE, adhesive	30,5x29 cm	1
PU	0	Drape, reinforced, adhesive	240x196 cm	1
SU	0	Drape adhesive	75x75 cm	6
		Adhesive tape, PE	8x40 cm	2
		Back table cover, 50µm	120x140 cm	1
		CSR wrap	112x112 cm	1
		Back table cover, folded, 50µm	140x190 cm	1
		Back table cover, 90µm	160x240 cm	1




Drapes & Packs

Orthopedic Surgery

► Split Drapes and Packs


U-Drape

Item no.	8476CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	CSR wrap	112x112 cm	1
Material	PE	U-drape PE, split 15x53cm, adhesive	152x183 cm	1
PU	40			
SU	0			




U-Drape

Item no.	8474CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	U-drape PE, split 15x53cm, adhesive	152x152 cm	1
Material	PE	CSR wrap	75x75 cm	1
PU	36			
SU	0			



U-Drape

Item no.	8475CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	U-drape PE, split 20x76cm, adhesive	152x213 cm	1
Material	PE	CSR wrap	60x60 cm	1
PU	28			
SU	0			




Drapes & Packs

Orthopedic Surgery

► Split Drapes and Packs


Universal split pack

Item no.	ES10132CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Head sheet, adhesive	240x150 cm	1
Material	ESS - 2 layers	Split sheet 200x250cm, with adh split 7x95cm	200x250cm	1
PU	8	Absorbent Hand towels	20x30 cm	4
SU	0	Back table cover	140x190 cm	1
		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1




Universal split pack

Item no.	DYNJPE1060	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Adhesive drape</td> <td>38x66 cm</td> <td>4</td>		Adhesive drape	38x66 cm	4
Material	Eclipse - Spunl	Top drape, adh.	259x135 cm	1
PU	4	Suture bag, adh.	16x31 cm	1
SU	0	U drape, adh. split 6x94 cm	196x274 cm	1
		Back table cover	125x229 cm	1
		Mayo stand cover, SMS reinf. 46x61 cm	61x135 cm	1



Universal Split Pack

Item no.	ES29113CE	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Head drape, adhesive</td> <td>260x182 cm</td> <td>1</td>		Head drape, adhesive	260x182 cm	1
Material	ESS - 2 layers	Split drape, reinforced, 8.5x107cm adhesive split	220x300 cm	1
PU	4	Drape	75x90cm	1
SU	555	Adhesive tape	9x50 cm	2
		Back table cover,	152x190 cm	2
		Mayo Stand Cover, absorbent reinf. 55x88 cm	80x142 cm	1
		Absorbent Hand towels	20x30 cm	2




Drapes & Packs

Orthopedic Surgery

► Split Drapes and Packs


Universal split pack VI

Item no.	TB29113CE	Description	Dimensions	Quantity
eClass-Nr.	Ultimate - 3 la	Absorbent Towel	39x58 cm	1
Material		Instrument pouch in clear plastic with adhesive	33x41cm	1
PU	6	Split drape reinforced, adhesive split 15x116cm	196x305 cm	1
SU	0	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
		CSR wrap	86x86 cm	1
		Back table cover, 75µm	160x240 cm	1
		Bar drape, reinforcement, 1 cord holding tab, adhesive	269x157 cm	1




Universal split pack VII

Item no.	TB29112CE	Description	Dimensions	Quantity
eClass-Nr.	Ultimate - 3 la	Absorbent Towel	39x58 cm	1
Material		Instrument pouch in clear plastic with adhesive	33x41cm	1
PU	6	Split drape, absorbent reinforcement, 2 cord holding tabs	196x305 cm	1
SU	0	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
		CSR wrap	112x112 cm	1
		Back table cover, 75µm	160x240 cm	1
		Bar drape, reinforcement, 1 cord holding tab, adhesive	269x157 cm	1



Universal split pack,

Item no.	TB29131CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Head sheet, reinforced, adhesive	251x135 cm	1
Material	Ultimate - 3 la	Split drape reinforced, adhesive split 15x116cm	196x305 cm	1
PU	6	Suture bag, adh.	16x31 cm	1
SU	0	MAYO Stand cover, S, 50µ, absorbent reinf. 46x61 cm	59x141 cm	1
		Back table cover, reinforced, 50µm	125x230 cm	1
		Hand towels	30x20 cm	4




Drapes & Packs

Orthopedic Surgery

► Split Drapes and Packs


Universal split pack,

Item no.	TB29132CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Absorbent Hand towels	20x30 cm	4
Material	Ultimate - 3 la	Mayo stand cover, reinforced	80x142 cm	1
PU	5	Head sheet, reinforced, adhesive	254x183 cm	1
SU	0	Split sheet, reinforced, split 15x102cm, adhesive	196x254 cm	1
		Back table cover	140x190 cm	1



XL split drape

Item no.	DYNJPE2455	Description	Dimensions	Quantity
eClass-Nr.		Wrap	76x76 cm	1
Material	Eclipse - Spunl	U drape, adh. split 15x102 cm	229x274 cm	1
PU	10			
SU	0			



Drapes & Packs

Paediatric Surgery

► Other Drapes and Packs


Neonatal drape, 112x152cm with 8x10cm fenestration, 1, 9cm width tape

Item no.	TB29494CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-10	Neonatal drape, reinforced, fenestration 8x10cm, 1, 9cm wide adhe..	112x152 cm	1
Material	Ultimate - 3 la	CSR wrap	60x60 cm	1
PU	38			
SU	0			




Paediatric Drape

Item no.	TB29493CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-10	Pediatric drape, reinforced, adhesive	196/254x286cm	1
Material	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU	15			
SU	0			



Paediatric drape

Item no.	TB9777CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-10	Directional incision drape, 10cm adhesive fenestration	76x56 cm	1
Material	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU	80			
SU	0			



Drapes & Packs

Paediatric Surgery

► Other Drapes and Packs


Paediatric drape,

Item no.	TB29492CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-10	Paediatric drape, reinforced, fenestration 15x15cm, adhesive	196x274 cm	1
Material	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU	18			
SU	0			




Paediatric laparotomy drape

Item no.	DYNJPE3007	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	ECL - Spunlac	Paediatric laparotomy drape, adh. fen. 15x15 cm	196x274 cm	1
PU	10			
SU	0			



Paediatric pack

Item no.	DYJPEPDSM	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, reinforced	80x142 cm	1
Material	ADV - SMMMS	Minor procedure drape, adh. fen. 15x15 cm	198/249x307 cm	1
PU	6	Back table cover,	152x190 cm	1
SU	0	Hand towel	34x36 cm	2



Drapes & Packs

Paediatric Surgery

► Other Drapes and Packs

Pediatric Incise drape

Item no.	DYNJPE3107	Description	Dimensions	Quantity
eClass-Nr.		Pediatric Incise drape	195x274 cm	1
Material	ECL - Spunlac	SMS wrap		1
PU	24			
SU	0			




Drapes & Packs

Paediatric Surgery

► Split Packs

Paediatric pack

Item no.	DYJPEPDSPPSM	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-10	Mayo stand cover, SMS reinf. 55x88 cm	80x142 cm	1
Material	ADV - SMMMS	Head drape, adh.	250x157 cm	1
PU	5	U drape, adh. split 6x35 cm	200x178 cm	1
SU	0	Back table cover,	152x190 cm	1




Drapes & Packs

Paediatric Surgery

► Universal Drapes and Packs


Paediatric pack

Item no.	DYJPEPDUPSM	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-10	Head drape, adh.	250x157 cm	1
Material	ADV - SMMMS	Side drape, adh.	75x102 cm	2
PU	8	Foot drape, adh.	190x183 cm	1
SU	0	Back table cover,	152x190 cm	1
		Hand towel	34x36 cm	1
		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1




Paediatric pack

Item no.	DYJPEPDUPEC	Description	Dimensions	Quantity
eClass-Nr.		Side drape, adh.	75x102 cm	2
Material	Eclipse - Spunl	Foot drape, adh.	190x178 cm	1
PU	8	Head drape, adh.	250x157 cm	1
SU	0	Hand towel	34x36 cm	4
		Mayo stand cover, SMS reinf. 55x88 cm	80x142 cm	1



Pediatric universal pack

Item no.	CP29174CE	Description	Dimensions	Quantity
eClass-Nr.		Hand towels	30x39 cm	1
Material	ECL - Spunlac	Adhesive tape, PE	8x40 cm	2
PU	10	Drape, 30x15 cm reinforcement, adh.	123x75 cm	2
SU	0	Drape, 30x15 cm reinforcement, adh.	112x183 cm	2
		Back table cover	140x190 cm	1
		CSR wrap	76x76 cm	1




Drapes & Packs

Paediatric Surgery

► Universal Drapes and Packs


Pediatric universal pack

Item no.	SM29174CE	Description	Dimensions	Quantity
eClass-Nr.		Head drape, adh.	157x250 cm	1
Material	ADV - SMMMS	Side drape, adh.	100x98 cm	2
PU	8	Foot drape, adh.	190x178 cm	1
SU	0	SMS wrap		1
		OP tape	10x50 cm	1
		Hand towel	34x36 cm	2
		Back table cover	140x190 cm	1




Pediatric universal pack

Item no.	ES29175CE	Description	Dimensions	Quantity
eClass-Nr.		Head drape, adhesive	160x260 cm	1
Material	ESS - 2 layers	foot drape	175x190 cm	1
PU	10	Drape adhesive	75x75 cm	2
SU	555	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
		Absorbent Hand towels	20x30 cm	4
		Adhesive tape	9x50 cm	1
		Back table cover	140x190 cm	1



Universal pack

Item no.	ES29103CE	Description	Dimensions	Quantity
eClass-Nr.		Hook and loop tube holder	14x2.5 cm	1
Material	ESS - 2 layers	foot drape	175x190 cm	1
PU	8	bi-laminate side drape	90x75 cm	2
SU	555	Adhesive tape	9x50 cm	1
		bi-laminate head drape	240x150 cm	1
		Absorbent Hand towels	20x30 cm	2
		SMS wrap		1
		Back table cover, reinforced, 50µm	140x190 cm	1




Drapes & Packs

Supplementary Program

► Adhesive drape sheets

Adhesive drape

Item no.	DG3037A	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Drape Sheet, 98x110cm, w/adh MED 5x95cm 2 sides	98x110 cm	1
Material	Ultimate - 3 la	CSR wrap	60x60 cm	1
PU	40			
SU	0			




Adhesive drape

Item no.	TB30056CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Drape, w/ adhesive	150x196 cm	1
Material	Ultimate - 3 la	CSR wrap	60x60 cm	1
PU	24			
SU	0			



Adhesive drape

Item no.	DYJPEADSSM3	Description	Dimensions	Quantity
eClass-Nr.	ADV - SMMMS	Wrap	76x76 cm	1
Material		Adhesive drape	115x115 cm	1
PU	20			
SU	0			




Drapes & Packs

Supplementary Program

► Adhesive drape sheets


Adhesive drape

Item no.	DYJPEADSSM4	Description	Dimensions	Quantity
eClass-Nr.		Wrap	76x76 cm	1
Material	ADV - SMMMS	Adhesive drape	180x190 cm	1
PU	20			
SU	0			




Adhesive drape

Item no.	DYJPEADSEC1	Description	Dimensions	Quantity
eClass-Nr.		Wrap	76x76 cm	1
Material	ECL - Spunlac	Adhesive drape	178x150 cm	1
PU	20			
SU	0			



Adhesive drape

Item no.	DYNJPE2405	Description	Dimensions	Quantity
eClass-Nr.		Adhesive drape	66x38 cm	2
Material	ECL - Spunlac			
PU	100			
SU	0			




Drapes & Packs

Supplementary Program

► Adhesive drape sheets

Adhesive drape

Item no.	DYNJPE2408	Description	Dimensions	Quantity
eClass-Nr.		Adhesive drape	51x76 cm	4
Material	Eclipse - Spunl			
PU	50			
SU	0			



Adhesive drape

Item no.	DYNJPE2406	Description	Dimensions	Quantity
eClass-Nr.		Adhesive drape	38x66 cm	4
Material	ECL - Spunlac			
PU	50			
SU	0			



Adhesive drape

Item no.	ES30056CE	Description	Dimensions	Quantity
eClass-Nr.		Drape, w/ adhesive	150x196 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	48			
SU	555			



Drapes & Packs

Supplementary Program

► Adhesive drape sheets


Adhesive tape

Item no.	DYJPEADSSM1	Description	Dimensions	Quantity
eClass-Nr.		Wrap	76x76 cm	1
Material	ADV - SMMMS	Adhesive drape	178x152 cm	1
PU	20			
SU	0			




Bottom drape

Item no.	DYJPEBD020	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Wrap</td> <td>91x91 cm</td> <td>1</td>		Wrap	91x91 cm	1
Material	ECL - Spunlac	Bottom drape, adh.	178x178 cm	1
PU	15			
SU	0			



Drape sheet

Item no.	DYNJPE2412TS	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Wrap</td> <td>61x61 cm</td> <td>1</td>		Wrap	61x61 cm	1
Material	ADV - SMMMS	Drape sheet w/ adh. 31x2.5 cm	102x178 cm	1
PU	50			
SU	0			



Drapes & Packs

Supplementary Program

► Adhesive drape sheets


Drape sheet with adhesive

Item no.	DYJPEADSBL1	Description	Dimensions	Quantity
eClass-Nr.		Adhesive drape, 150x178cm	150x178 cm	1
Material	ESS - 2 layers	CSR wrap	86x86 cm	1
PU	20			
SU	0			



Drape with adhesive

Item no.	TB30065CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Adhesive drape	75x91 cm	1
Material	Ultimate - 3 la	CSR wrap	60x60 cm	1
PU	36			
SU	0			



Drape with adhesive

Item no.	TB933CEA	Description	Dimensions	Quantity
eClass-Nr.		Drape with adhesive	300x221 cm	1
Material	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU	16			
SU	0			




Drapes & Packs

Supplementary Program

► Adhesive drape sheets


Drape, adhesive

Item no.	TB9471CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Champ latéral adhésif	98x98 cm	1
Material	Ultimate - 3 la	CSR wrap	112x112 cm	1
PU	68			
SU	0			



Drape, adhesive

Item no.	ES15211CE	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Drape</td> <td>75x90cm</td> <td>1</td>		Drape	75x90cm	1
Material	ESS - 2 layers			
PU	160			
SU	0			



Drape, adhesive

Item no.	ES15213CE	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Small drape, adhesive</td> <td>150x90 cm</td> <td>1</td>		Small drape, adhesive	150x90 cm	1
Material	ESS - 2 layers			
PU	54			
SU	0			



Drapes & Packs

Supplementary Program

► Adhesive drape sheets

Drape, adhesive

Item no.	ES15215CE	Description	Dimensions	Quantity
eClass-Nr.		Drape, adhesive	45x50 cm	1
Material	ESS - 2 layers			
PU	400			
SU	0			




Drape, adhesive

Item no.	ES15212CE	Description	Dimensions	Quantity
eClass-Nr.		Drape adhesive	75x75 cm	1
Material	ESS - 2 layers			
PU	160			
SU	0			



Drape, adhesive

Item no.	ES15216CE	Description	Dimensions	Quantity
eClass-Nr.		Small drape, adhesive	100x100 cm	1
Material	ESS - 2 layers			
PU	45			
SU	0			




Drapes & Packs

Supplementary Program

► Adhesive drape sheets


Drape, w/adh

Item no.	9471CEB	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Drape, 98 x 98 cm w/adh	98x98 cm	1
Material	ECL - Spunlac	Sterile Wrap	60 x 60 cm	1
PU	44			
SU	0			




Drapes, adhesive

Item no.	TB307554CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-03-01	Drape, adhesive	49x65 cm	2
Material	Ultimate - 3 la	CSR wrap	60x60 cm	1
PU	108			
SU	0			



Foot drape

Item no.	DYJPEFTDSM1	Description	Dimensions	Quantity
eClass-Nr. <td>ADV - SMMMS</td> <td>Wrap</td> <td>76x76 cm</td> <td>1</td>	ADV - SMMMS	Wrap	76x76 cm	1
Material		Foot drape, adh.	183x175 cm	1
PU	10			
SU	0			




Drapes & Packs

Supplementary Program

► Adhesive drape sheets


Foot drape with adhesive

Item no.	TB29469CE	Description	Dimensions	Quantity
eClass-Nr.		Foot drape, reinf., 2 cord holding tabs, adh.	183x191 cm	1
Material	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU	20			
SU	0			




Foot drape with adhesive

Item no.	TB31600CE	Description	Dimensions	Quantity
eClass-Nr. <td><td>Foot sheet, adhesive</td><td>196x200 cm</td><td>1</td></td>	<td>Foot sheet, adhesive</td> <td>196x200 cm</td> <td>1</td>	Foot sheet, adhesive	196x200 cm	1
Material	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU	26			
SU	0			



Foot drape with adhesive

Item no.	ES11602CE	Description	Dimensions	Quantity
eClass-Nr. <td><td>Foot drape, adhesive</td><td>190x200 cm</td><td>1</td></td>	<td>Foot drape, adhesive</td> <td>190x200 cm</td> <td>1</td>	Foot drape, adhesive	190x200 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	36			
SU	555			



Drapes & Packs

Supplementary Program

► Adhesive drape sheets

Foot drape with adhesive

Item no.	ES29470CE	Description	Dimensions	Quantity
eClass-Nr.		foot drape	175x190 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	42			
SU	555			




Foot drape with adhesive

Item no.	ES29471CE	Description	Dimensions	Quantity
eClass-Nr.		foot drape	200x200 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	36			
SU	555			



Foot drape with adhesive

Item no.	AGBDV010A	Description	Dimensions	Quantity
eClass-Nr.		CSR wrap	76x76 cm	1
Material	Ultimate - 3 la	Foot Drape, reinf, adhesive	196x200 cm	1
PU	18			
SU	0			




Drapes & Packs

Supplementary Program

► Adhesive drape sheets


Head drape

Item no.	DYJPEHEDSM1	Description	Dimensions	Quantity
eClass-Nr.		Wrap	100x100 cm	1
Material	ADV - SMMMS	Head drape, adh.	300x183 cm	1
PU	8			
SU	0			



Head drape

Item no.	DYJPEHEDSM2	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Wrap</td> <td>76x76 cm</td> <td>1</td>		Wrap	76x76 cm	1
Material	ADV - SMMMS	Head drape, adh.	152x250 cm	1
PU	8			
SU	0			



Head drape

Item no.	TB31900CE	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Head drape, adhesive</td> <td>260x160 cm</td> <td>1</td>		Head drape, adhesive	260x160 cm	1
Material	Ultimate - 3 Ia			
PU	26			
SU	0			



Drapes & Packs

Supplementary Program

► Adhesive drape sheets


Head drape

Item no.	TB29479CE	Description	Dimensions	Quantity
eClass-Nr.		Head drape, reinforcement, adhesive, 1 cord holding tab	254x183 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	15			
SU	0			



Head drape

Item no.	DYJPEHEDSM6	Description	Dimensions	Quantity
eClass-Nr. <td><td>Head drape, adh.</td><td>340x200 cm</td><td>1</td></td>	<td>Head drape, adh.</td> <td>340x200 cm</td> <td>1</td>	Head drape, adh.	340x200 cm	1
Material	ADV - SMMM5	SMS wrap		1
PU	20			
SU	0			



Head drape

Item no.	ES10483CE	Description	Dimensions	Quantity
eClass-Nr. <td><td>Head drape, adhesive</td><td>160x260 cm</td><td>1</td></td>	<td>Head drape, adhesive</td> <td>160x260 cm</td> <td>1</td>	Head drape, adhesive	160x260 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	34			
SU	555			




Drapes & Packs

Supplementary Program

► Adhesive drape sheets

Head drape

Item no.	ES29482CE	Description	Dimensions	Quantity
eClass-Nr.		bi-laminate head drape	300x170 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	26			
SU	555			




Head drape

Item no.	ES29483CE	Description	Dimensions	Quantity
eClass-Nr.		bi-laminate head drape	240x150 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	38			
SU	555			



Side drape

Item no.	DYJPESIDSM1	Description	Dimensions	Quantity
eClass-Nr.		Wrap	66x66 cm	1
Material	ADV - SMMMS	Side drape, adh.	112x102 cm	1
PU	20			
SU	0			




Drapes & Packs

Supplementary Program

► Adhesive drape sheets


Side drape

Item no.	DYJPESD121	Description	Dimensions	Quantity
eClass-Nr.		Wrap	76x76 cm	1
Material	Polyrayon	Side drape, adh.	102x102 cm	1
PU	20			
SU	0			




Side drape

Item no.	DYJPESIDEC1	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Wrap</td> <td>76x76 cm</td> <td>1</td>		Wrap	76x76 cm	1
Material	Eclipse - Spunl	Side drape, adh.	89x112 cm	1
PU	20			
SU	0			



Side drape with adhesive

Item no.	TB31300CE	Description	Dimensions	Quantity
eClass-Nr. <td>34-35-01-01</td> <td>Side drape</td> <td>75x98 cm</td> <td>1</td>	34-35-01-01	Side drape	75x98 cm	1
Material	Ultimate - 3 la	CSR wrap	60x60 cm	1
PU	70			
SU	0			




Drapes & Packs

Supplementary Program

► Adhesive drape sheets


Side drape with adhesive

Item no.	TB29439CEA	Description	Dimensions	Quantity
eClass-Nr.		Side drape, reinforced, 98 x 112 cm	98 x 112 cm	1
Material	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU	36			
SU	0			




Side drape with adhesive

Item no.	ES15210CE	Description	Dimensions	Quantity
eClass-Nr.		Side sheet, adhesive	98x75 cm	1
Material	ESS - 2 layers			
PU	60			
SU	555			



Side drape with adhesive

Item no.	ES29443CE	Description	Dimensions	Quantity
eClass-Nr.		bi-laminate side drape	90x75 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	120			
SU	555			




Drapes & Packs

Supplementary Program

► Adhesive drape sheets

Side drape with adhesive

Item no.	ES29444CE	Description	Dimensions	Quantity
eClass-Nr.		bi-laminate side drape	110x95 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	90			
SU	555			




Side drapes with adhesive

Item no.	TB31301CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Side drape	75x98 cm	2
Material	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU	36			
SU	0			



Side sheet, reinforced, adhesive

Item no.	2941CEA	Description	Dimensions	Quantity
eClass-Nr.		Side drape, reinforced, adhesive	110x98 cm	1
Material	Eclipse - Spunl	CSR wrap	60x60 cm	1
PU	42			
SU	0			




Drapes & Packs

Supplementary Program

► Adhesive drape sheets

Two U drapes, adhesive

Item no.	ES15214CE	Description	Dimensions	Quantity
eClass-Nr.		U drape, adhesive	45x75 cm	2
Material	ESS - 2 layers			
PU	160			
SU	0			




Drapes & Packs

Supplementary Program

► Adhesive drape sheets

Drape, reinforced, adhesive
241cm



Item no.	29340CE
eClass-Nr.	34-35-01-01
Material	ECL - Spunlace
Dimensions	196 x 240 cm
PU	16
SU	0


Drapes & Packs

Supplementary Program

► Cotton Towel

Absorbent towels pack

Item no.	DYJPEAT012	Description	Dimensions	Quantity
eClass-Nr.		Wrap	45x45 cm	1
Material	Cellulose	Absorbent towel, unfolded	51x38 cm	2
PU	74			
SU	0			



Absorbent towels, cotton, blue

Item no.	28700002CEA	Description	Dimensions	Quantity
eClass-Nr.	34-21-92-02	absorbant towels, blue	43 x 66 cm	2
Material	Cotton			
PU	40			
SU	0			



Absorbent towels, cotton, blue

Item no.	28700004CEA	Description	Dimensions	Quantity
eClass-Nr.	34-21-92-02	absorbant towels, blue	43 x 66 cm	4
Material	Cotton			
PU	20			
SU	0			



Drapes & Packs

Supplementary Program

► Cotton Towel

Absorbent towels, cotton, XR, blue

Item no.	28700002XR	Description	Dimensions	Quantity
eClass-Nr.	34-21-92-02	Absorbent towels, cotton, XR, blue	43x66 cm	2
Material	Cotton			
PU	40			
SU	0			



Absorbent towels, cotton, XR, blue

Item no.	28700004XR	Description	Dimensions	Quantity
eClass-Nr.	34-21-92-02	Absorbent towels, cotton, XR, blue	43x66 cm	4
Material	Cotton			
PU	20			
SU	0			



Drapes & Packs

Supplementary Program

► Cotton Towel

Absorbent towels, cotton,
green



Item no.	28200004CEA
eClass-Nr.	34-21-92-02
Material	Cotton
Dimensions	43 x 66 cm
PU	20
SU	0

Drapes & Packs

Supplementary Program

► Drape sheets

Drape

Item no.	TB9355CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Drape	112x180 cm	1
Material	Ultimate - 3 la	CSR wrap	112x112 cm	1
PU	20			
SU	0			



Drape

Item no.	TB9358CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Half drape	112x145 cm	1
Material	Ultimate - 3 la	CSR wrap	112x112 cm	1
PU	20			
SU	0			



Drape

Item no.	ES15222CE	Description	Dimensions	Quantity
eClass-Nr.	ESS - 2 layers	Small drape, 75x75cm, without adhesive	75x75 cm	1
Material	160			
PU	0			
SU	0			



Drapes & Packs

Supplementary Program

► Drape sheets

Drape

Item no.	9355CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Drape	112x180 cm	1
Material	ECL - Spunlac	CSR wrap	112x112 cm	1
PU	28			
SU	0			




Drape

Item no.	9444CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Drape sheet	196x249 cm	1
Material	Eclipse - Spunl	CSR wrap	76x76 cm	1
PU	18			
SU	0			



Drape

Item no.	9461CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Drape	98x98 cm	1
Material	ECL - Spunlac	CSR wrap	60x60 cm	1
PU	48			
SU	0			




Drapes & Packs

Supplementary Program

► Drape sheets

Drape

Item no.	ES15221CE	Description	Dimensions	Quantity
eClass-Nr.		drape sheet	75x90cm	1
Material	ESS - 2 layers			
PU	160			
SU	0			



Drape

Item no.	ES15223CE	Description	Dimensions	Quantity
eClass-Nr.		Drape	90x150 cm	1
Material	ESS - 2 layers			
PU	56			
SU	0			



Drape

Item no.	ES15225CE	Description	Dimensions	Quantity
eClass-Nr.		Small drape	45x50 cm	1
Material	ESS - 2 layers			
PU	400			
SU	0			



Drapes & Packs

Supplementary Program

► Drape sheets


Drape

Item no.	ES15226CE	Description	Dimensions	Quantity
eClass-Nr.		Drape sheet	45x75cm	1
Material	ESS - 2 layers			
PU	320			
SU	0			




Drape

Item no.	ES15235CE	Description	Dimensions	Quantity
eClass-Nr.		drape sheet	100x100 cm	1
Material	ESS - 2 layers			
PU	60			
SU	0			



Drape

Item no.	ES15236CE	Description	Dimensions	Quantity
eClass-Nr.		drape sheet	100x120 cm	1
Material	ESS - 2 layers			
PU	51			
SU	0			



Drapes & Packs

Supplementary Program

► Drape sheets


Drape sheet

Item no.	DYNJPE2414SM	Description	Dimensions	Quantity
eClass-Nr.		Wrap	61x61 cm	1
Material	ADV - SMMMS	Drape sheet	135x196 cm	1
PU	1			
SU	0			




Drape sheet

Item no.	DYJPEDSHSM10	Description	Dimensions	Quantity
eClass-Nr.		Wrap	61x61 cm	1
Material	ADV - SMMMS	Drape sheet	89x102 cm	1
PU	20			
SU	0			



Drape sheet

Item no.	DYNJPE2410SM	Description	Dimensions	Quantity
eClass-Nr.		Wrap	61x61 cm	1
Material	ADV - SMMMS	Drape sheet	102x147 cm	1
PU	20			
SU	0			



Drapes & Packs

Supplementary Program

► Drape sheets

Drape sheet

Item no.	DYNJPE2417SM	Description	Dimensions	Quantity
eClass-Nr.		Wrap	61x61 cm	1
Material	ADV - SMMMS	Drape sheet	178x216 cm	1
PU	20			
SU	0			




Drape sheet

Item no.	DYNJPE2419SM	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Wrap</td> <td>61x61 cm</td> <td>1</td>		Wrap	61x61 cm	1
Material	ADV - SMMMS	Drape sheet	178x254 cm	1
PU	12			
SU	0			



Drape sheet

Item no.	DYNJPE2411	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Wrap</td> <td>61x61 cm</td> <td>1</td>		Wrap	61x61 cm	1
Material	wetlaid	Drape sheet	102x147 cm	1
PU	20			
SU	0			




Drapes & Packs

Supplementary Program

► Drape sheets


Drape sheet

Item no.	DYNJPE2415	Description	Dimensions	Quantity
eClass-Nr.		Wrap	61x61 cm	1
Material	wetlaid	Drape sheet	135x196 cm	1
PU	20			
SU	0			



Drape Sheet, three quarter

Item no.	TB9349CEB	Description	Dimensions	Quantity
eClass-Nr.	34-35-03-01	Three quarter drape	145x193 cm	1
Material	Ultimate - 3 la	CSR wrap	112x112 cm	1
PU	16			
SU	0			



Drape,

Item no.	TB9444CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Drape sheet	196x249 cm	1
Material	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU	16			
SU	0			



Drapes & Packs

Supplementary Program

► Drape sheets

Drape,

Item no.	TB9446CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Drape sheet	216x183 cm	1
Material	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU	18			
SU	0			



Drape,

Item no.	TB9461CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Utility drape	98x98 cm	1
Material	Ultimate - 3 la	CSR wrap	112x112 cm	1
PU	68			
SU	0			



Drape, 112 x 145 cm

Item no.	9358CEC	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Drape, 112 x 145 cm	112 x 145 cm	1
Material	ECL - Spunlac	Sterile Wrap	60 x 60 cm	1
PU	28			
SU	0			




Drapes & Packs

Supplementary Program

► Drape sheets

Drape, 145 x 193 cm

Item no.	9349CEB	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Drape, 145 x 193 cm	145x193 cm	1
Material	ECL - Spunlac	CSR wrap	112x112 cm	1
PU	16			
SU	0			



Drape, Sheet, Fan-Folded Large

Item no.	9446CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Drape sheet, Large	183x216 cm	1
Material	Eclipse - Spunl	CSR wrap	76x76 cm	1
PU	22			
SU	0			




Drapes & Packs

Supplementary Program

► Fenestrated Drapes


Drape utility, FEN 10x15cm CEN w/o adh., NS

Item no.	TB9462CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-03-01	CSR wrap	60x60 cm	1
Material	Ultimate - 3 la	Utility Drape, 43x56 cm reinforcement, 10x15 cm fenestration	152x112 cm	1
PU	32			
SU	0			




Drape, fenestration diameter 10cm

Item no.	ES15209CE	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Small drape, 90x75 cm, with 10cm fenestration</td> <td>90x75 cm</td> <td>1</td>		Small drape, 90x75 cm, with 10cm fenestration	90x75 cm	1
Material	ESS - 2 layers			
PU	180			
SU	555			



Drape, fenestration diameter 5cm

Item no.	ES15227CE	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Small drape, 60x50 cm, with 5cm fenestration</td> <td>60x50 cm</td> <td>1</td>		Small drape, 60x50 cm, with 5cm fenestration	60x50 cm	1
Material	ESS - 2 layers			
PU	320			
SU	0			



Drapes & Packs

Supplementary Program

► Fenestrated Drapes


Drape, fenestration diameter 7cm

Item no.	ES15230CE	Description	Dimensions	Quantity
eClass-Nr.		bi-laminate small drape 45x75 cm, 7 cm diameter circular fenestration	45x75cm	1
Material	ESS - 2 layers			
PU	180			
SU	0			




Drape, fenestration diameter 8cm

Item no.	ES152275CE	Description	Dimensions	Quantity
eClass-Nr.		bi-laminate small drape 90x75 cm, 8 cm diameter circular fenestration	90x75cm	1
Material	ESS - 2 layers			
PU	90			
SU	0			



Emergency drape

Item no.	S951CE	Description	Dimensions	Quantity
eClass-Nr.		Emergency drape, 12 cm diam. fen.	65x53 cm	1
Material	ADV - SMMMS			
PU	200			
SU	0			




Drapes & Packs

Supplementary Program

► Fenestrated Drapes w. Adh.


Adherent window drape

Item no.	DYNJPE2493	Description	Dimensions	Quantity
eClass-Nr.		Wrap	61x61 cm	1
Material	Eclipse - Spunl	Adherent window drape, oval adh. fen. w/ full incise film 9x19 cm	114x102 cm	2
PU	34			
SU	0			



Anaesthetic drape, fenestration

Item no.	9144FRC	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-02	Anaesthetic drape, fenestrated	196x300 cm	1
Material	Eclipse - Spunl	CSR wrap	85x85 cm	1
PU	10			
SU	0			



Drape 120x150 cm, 6x8 cm adhesive fenestration

Item no.	ES152183CE	Description	Dimensions	Quantity
eClass-Nr.		Drape, with a 6x8 cm oblong adhesive fenestration, centered	120x150 cm	1
Material	ESS - 2 layers			
PU	60			
SU	0			




Drapes & Packs

Supplementary Program

► Fenestrated Drapes w. Adh.


Drape sheet with 6x4cm adhesive fenestration

Item no.	DYNJPE7013	Description	Dimensions	Quantity
eClass-Nr.		Wrap	61x61 cm	1
Material	ECL - Spunlac	Basic eye drape, oval fen. 6x4 cm w/ adh.	102x142 cm	1
PU	40			
SU	0			



Drape with adhesive fenestration

Item no.	ES15207CEA	Description	Dimensions	Quantity
eClass-Nr.		Small drape, with adhesive fenestration 10cm	75x90cm	1
Material	ESS - 2 layers			
PU	180			
SU	0			



Drape with fenestration

Item no.	TB30668CE	Description	Dimensions	Quantity
eClass-Nr.		Drape with 10cm dia fenestration, with soft adhesive tape	150x98 cm	1
Material	Ultimate - 3 la	CSR wrap	60x60 cm	1
PU	42			
SU	0			




Drapes & Packs

Supplementary Program

► Fenestrated Drapes w. Adh.


Drape, adhesive fenestration 8x6cm

Item no.	ES15228CE	Description	Dimensions	Quantity
eClass-Nr.		Drape, 6x8 cm adhesive fenestration	75x90 cm	1
Material	ESS - 2 layers			
PU	80			
SU	0			




Drape, adhesive fenestration diam. 10cm

Item no.	ES15207CE	Description	Dimensions	Quantity
eClass-Nr.		Small drape, with adhesive fenestration 10cm	75x90cm	1
Material	ESS - 2 layers			
PU	320			
SU	555			



Drape, adhesive fenestration diam. 5cm

Item no.	ES15217CE	Description	Dimensions	Quantity
eClass-Nr.		Small drape, with adhesive fenestration 5cm	50x60cm	1
Material	ESS - 2 layers			
PU	320			
SU	0			



Drapes & Packs

Supplementary Program

► Fenestrated Drapes w. Adh.


Drape, adhesive fenestration diam. 7cm

Item no.	ES15231CE	Description	Dimensions	Quantity
eClass-Nr.		Small drape, with adhesive fenestration 7cm	45x75 cm	1
Material	ESS - 2 layers			
PU	150			
SU	0			



Drape, adhesive fenestration diam. 9x12cm

Item no.	ES15219CE	Description	Dimensions	Quantity
eClass-Nr.		Drape fenestrated, fenestration 12x9cm	100x100 cm	1
Material	ESS - 2 layers			
PU	45			
SU	0			



Drape, fenestration diameter 10cm, adhesive

Item no.	TB30669CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-01	Drape sheet w/ adhesive fen. diam.10cm	98x150 cm	1
Material	Ultimate - 3 la	CSR wrap	60x60 cm	1
PU	42			
SU	0			



Drapes & Packs

Supplementary Program

► Fenestrated Drapes w. Adh.


Drape, oval fenestration 8x5cm, adhesive

Item no.	ES15218CE	Description	Dimensions	Quantity
eClass-Nr.		Drape, 5x8 cm adhesive fenestration	120x150 cm	1
Material	ESS - 2 layers			
PU	30			
SU	0			




EENT drape

Item no.	DYNJPE7012SM	Description	Dimensions	Quantity
eClass-Nr. <td><td>Wrap</td><td>80x80 cm</td><td>1</td></td>	<td>Wrap</td> <td>80x80 cm</td> <td>1</td>	Wrap	80x80 cm	1
Material	ADV - SMMMS	EENT drape, instrument pad 41x41 cm w/ adh. fen. 5x6 cm	203x229 cm	1
PU	10			
SU	0			



ENT drape

Item no.	TB30456CE	Description	Dimensions	Quantity
eClass-Nr. <td>34-35-01-05</td> <td>Fenestrated drape, fenestration diameter 12cm, adhesive</td> <td>220x240 cm</td> <td>1</td>	34-35-01-05	Fenestrated drape, fenestration diameter 12cm, adhesive	220x240 cm	1
Material	Ultimate - 3 la	CSR wrap	76x76 cm	1
PU	18			
SU	0			




Drapes & Packs

Supplementary Program

► Fenestrated Drapes w. Adh.


Epidural Drape 60x80 cm

Item no.	ES15003CE	Description	Dimensions	Quantity
eClass-Nr.		Epidural drape, w/ 8x8 cm diamond adh. fen.	60x80 cm	1
Material	ESS - 2 layers			
PU	180			
SU	0			



Epidural drape, oval fenestration 18x8cm

Item no.	TB30009CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-09	Drape Epidural, fen 8 x 18 cm oval cen w/ adh	80x120 cm	1
Material	Ultimate - 3 la	CSR wrap	60x60 cm	1
PU	32			
SU	0			



Epidural drape, oval fenestration 18x8cm

Item no.	ES30009CE	Description	Dimensions	Quantity
eClass-Nr.		Bilaminare epidural drape 120x80 cm, 18x8 cm adhesive fenestrat..	80x120 cm	1
Material	ESS - 2 layers	SMS wrap		1
PU	100			
SU	555			




Drapes & Packs

Supplementary Program

► Fenestrated Drapes w. Adh.


Fenestrated drape

Item no.	DYJPEFDSEC	Description	Dimensions	Quantity
eClass-Nr.		Wrap	60x60 cm	1
Material	Eclipse - Spunl	Fenestrated drape, adh. fen. 10x15 cm	102x157 cm	1
PU	20			
SU	0			



Minor procedure drape pack

Item no.	ES15312CE	Description	Dimensions	Quantity
eClass-Nr.		drape sheet	75x90cm	1
Material	ESS - 2 layers	Drape fenestrated, fenestration 12x9cm	75x90 cm	1
PU	50			
SU	0			



rape 120x150 cm, 10x7 cm adhesive fenestration

Item no.	ES152182CE	Description	Dimensions	Quantity
eClass-Nr.		Drape, 10x7 cm adhesive fenestration	120x150 cm	1
Material	ESS - 2 layers			
PU	60			
SU	0			



Drapes & Packs

Supplementary Program

► Leggings


Leggings

Item no.	DYNJPE2460SM	Description	Dimensions	Quantity
eClass-Nr.		Wrap	61x61 cm	1
Material	ADV - SMMMS	Leggings with cuff 15 cm	79x122 cm	2
PU	20			
SU	0			




Leggings

Item no.	DYNJPE2460	Description	Dimensions	Quantity
eClass-Nr.		Wrap	61x61 cm	1
Material	wetlaid	Leggings w/ cuff 15 cm	74x109 cm	2
PU	20			
SU	0			



Leggings, SMS

Item no.	8420CEB	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-10	Leggings	79x140 cm	1
Material	ADV - SMMMS	SMS wrap	60x60 cm	1
PU	26			
SU	0			




Drapes & Packs

Supplementary Program

► Leggings

Leggings, clear PE



Item no.	8430CEA
eClass-Nr.	34-35-04-10
Material	PE
Dimensions	80x143 cm
PU	24
SU	0

Drapes & Packs

Supplementary Program

► Perineal Cover

Pubic drape

Item no.	TB3007CE	Description	Dimensions	Quantity
eClass-Nr.		Perineal sheet, adhesive	16x37 cm	1
Material	Ultimate - 3 la			
PU	300			
SU	0			




Drapes & Packs

Supplementary Program

► Pouches

Fluid collection pouch

Item no.	4638CE	Description	Dimensions	Quantity
eClass-Nr.		Fluid collection pouch 51x61 cm, with malleable bar, filter and exit port	51 x 61 cm	1
Material	PE			
PU	40			
SU	0			




Fluid collection pouch

Item no.	IDE1016	Description	Dimensions	Quantity
eClass-Nr.		Fluid collection pouch	51x36 cm	1
Material	Matte PE film			
PU	40			
SU	0			



Instrument pouch

Item no.	IDE1018	Description	Dimensions	Quantity
eClass-Nr.		Instrument pouch	17x29 cm	1
Material	Matte PE film			
PU	40			
SU	0			




Drapes & Packs

Supplementary Program

► Pouches

Instrument pouch

Item no.	IDE1018L	Description	Dimensions	Quantity
eClass-Nr.		Instrument pouch	25x46 cm	1
Material	Matte PE film			
PU	40			
SU	0			




Instrument pouch, adhesive

Item no.	96860CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-14	Instrument pouch, PE, adhesive	30,5x29 cm	1
Material	PE			
PU	150			
SU	0			



Instrument pouch, adhesive

Item no.	9686CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-14	Instrument pouch, PE, adhesive	30,5x29 cm	1
Material	ECL - Spunlac			
PU	150			
SU	0			




Drapes & Packs

Supplementary Program

► Pouches


Pouch, Endoscopy, 2 compartments

Item no.	96890CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-14	Diathermy pouch	49x43 cm	1
Material	PE	Peel Pouch	22 x 35 cm	1
PU	250			
SU	0			




Pouch, Endoscopy, 2 compartments

Item no.	9689CE	Description	Dimensions	Quantity
eClass-Nr.		Instrument pouch with adhesive, 2 compartments	43,5x49 cm	1
Material	ECL - Spunlac	Peel Pouch	22 x 35 cm	1
PU	150			
SU	0			



Universal fluid collection pouch

Item no.	DYJPEFP020	Description	Dimensions	Quantity
eClass-Nr.		Wrap	45x45 cm	1
Material	PE	Fluid collection pouch, adh., filter, exit port, malleable bar	51x56 cm	1
PU	50			
SU	0			



Drapes & Packs

Supplementary Program


► Pouches

Instrument pouch w/ cord,
clear

No picture
available

Item no.	586560
eClass-Nr.	34-35-04-14
Material	PE
Dimensions	65x65 cm
PU	40
SU	0

Instrument pouch,



Item no.	9687CE
eClass-Nr.	34-35-04-14
Material	ECL - Spunlace
Dimensions	37 x 86 cm
PU	120
SU	0

Instrument pouch,
polyethylene, one pocket

No picture
available

Item no.	P80114
eClass-Nr.	34-35-04-14
Material	PE
Dimensions	
PU	240
SU	0


Drapes & Packs

Supplementary Program

► Protective Sleeves


Sleeve with knit cuff

Item no.	ES597CE	Description	Dimensions	Quantity
eClass-Nr.		Sleeve with knit cuff	14x62 cm	1
Material	ESS - 2 layers			
PU	225			
SU	555			



Sleeves with knit cuffs

Item no.	598CE	Description	Dimensions	Quantity
eClass-Nr.		Sleeve with knit cuff	62x14 cm	2
Material	Prevention Plu			
PU	80			
SU	0			



Spunlace sleeve with elastic ends, 13x60 cm

Item no.	DYNJPE2000	Description	Dimensions	Quantity
eClass-Nr.				
Material	ECL - Spunlac			
PU	60			
SU	0			




Drapes & Packs

Supplementary Program

► Protective Sleeves

Sterile sleeve

Item no.	599CEB	Description	Dimensions	Quantity
eClass-Nr.				
Material	SMMMS			
PU	150			
SU	0			




Drapes & Packs

Supplementary Program

► Ring Cover


Wound protector, diam. 12, 5cm

Item no.	D37106210	Description	Dimensions	Quantity
eClass-Nr.		Wound protector, fen. 12 cm diameter	90x90 cm	1
Material	Clear PE film			
PU	40			
SU	0			




Wound protector, diam. 18cm

Item no.	D37104910	Description	Dimensions	Quantity
eClass-Nr.		Wound protector, fen. 18 cm diameter	90x90 cm	1
Material	Clear PE film			
PU	40			
SU	0			



Wound protector, diam. 23cm

Item no.	D37106310	Description	Dimensions	Quantity
eClass-Nr.		Wound protector, fen. 23 cm diameter	90x90 cm	1
Material	Clear PE film			
PU	40			
SU	0			




Drapes & Packs

Supplementary Program

► Ring Cover


Wound protector, diam. 28cm

Item no.	D37106410	Description	Dimensions	Quantity
eClass-Nr.		Wound protector, fen. 28 cm diameter	90x90 cm	1
Material	Clear PE film			
PU	40			
SU	0			



Wound protector, diam. 7, 5cm

Item no.	D37612210	Description	Dimensions	Quantity
eClass-Nr.		Wound protector, fen. 7 cm diameter	90x90 cm	1
Material	Clear PE film			
PU	40			
SU	0			




Drapes & Packs

Supplementary Program

► Stockinettes

Impervious stockinette

Item no.	NONE22530	Description	Dimensions	Quantity
eClass-Nr.		Impervious stockinette	23x91 cm	1
Material	PP/PE			
PU	15			
SU	0			



Impervious stockinettes, Large

Item no.	NONE22534	Description	Dimensions	Quantity
eClass-Nr.		Impervious stockinette	23x122 cm	1
Material	PP/PE			
PU	12			
SU	0			



Impervious stockinettes, Medium

Item no.	NONE22525	Description	Dimensions	Quantity
eClass-Nr.		Impervious stockinette	20x97 cm	1
Material	PP/PE			
PU	15			
SU	0			




Drapes & Packs

Supplementary Program

► Stockinettes


Impervious stockinettes, X-Large

Item no.	NONE22535	Description	Dimensions	Quantity
eClass-Nr.		Impervious stockinette	30x122 cm	1
Material	PP/PE			
PU	10			
SU	0			




Impervious stockinettes, X-Large

Item no.	NONE22540	Description	Dimensions	Quantity
eClass-Nr.		Impervious stockinette	36x122 cm	1
Material	PP/PE			
PU	10			
SU	0			



Stockinette

Item no.	ES27009CE	Description	Dimensions	Quantity
eClass-Nr.		Stockinette	36.5x72 cm	1
Material	ESS - 2 layers	Adh extendable, turquoise	9x50 cm	1
PU	140			
SU	0			




Drapes & Packs

Supplementary Program

► Stockinettes

Impervious stockinettes,
Small



Item no.	NONE22520
eClass-Nr.	
Material	PP/PE
Dimensions	15x76 cm
PU	24
SU	0


Drapes & Packs

Supplementary Program

► Tapes/Velcro´s


Adhesive tape, non-woven, single packed

Item no.	3525CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-11			
Material	ECL - Spunlac			
PU	480			
SU	0			



Hook and loop tube holder

Item no.	DYNJPE2216	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Hook and loop tube holder</td> <td>14x2.5 cm</td> <td>1</td>		Hook and loop tube holder	14x2.5 cm	1
Material	PE			
PU	50			
SU	0			



Hook and loop tube holder

Item no.	7023WCE	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td></td> <td></td> <td></td>				
Material	PE			
PU	400			
SU	0			




Drapes & Packs

Supplementary Program

► Tapes/Velcro´s


Hook and loop tube holder, with XR

Item no.	SPT-7014WP	Description	Dimensions	Quantity
eClass-Nr.		Hook and loop tube holder	14x2.5 cm	1
Material	Polyester			
PU	2000			
SU	0			




OP tape

Item no.	DYJPETAPSM	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>OP tape</td> <td>10x50 cm</td> <td>1</td>		OP tape	10x50 cm	1
Material	ADV - SMMMS			
PU	120			
SU	0			



OP tape

Item no.	DYJPETAPSM2	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>OP tape</td> <td>10x50 cm</td> <td>2</td>		OP tape	10x50 cm	2
Material	ADV - SMMMS			
PU	100			
SU	0			




Drapes & Packs

Supplementary Program


► Tapes/Velcro´s

Adhesive tape, non-woven,
P2



Item no.	3550CEA
eClass-Nr.	34-35-04-11
Material	ECL - Spunlace
Dimensions	9x50 cm
PU	390
SU	0

Adhesive tape, polyethylene,
P2



Item no.	3532CEA
eClass-Nr.	34-35-04-11
Material	PE
Dimensions	8x40 cm
PU	390
SU	0


Drapes & Packs

Supplementary Program

► Tourniquet cover

Tourniquet cover, 38x45 cm with 5x8 cm fenestration

Item no.	SPT-2701	Description	Dimensions	Quantity
eClass-Nr.				
Material	Kraton			
PU	100			
SU	0			




Drapes & Packs

Supplementary Program

► Towels


Absorbent towels

Item no.	DYNJPE2220A	Description	Dimensions	Quantity
eClass-Nr.		Hand towel	37x69 cm	2
Material	Cellulose			
PU	100			
SU	0			



Towel, White, Absorbent

Item no.	7550CEA	Description	Dimensions	Quantity
eClass-Nr.		Absorbent Towel	39x58 cm	2
Material	Cellulose			
PU	150			
SU	0			




Drapes & Packs

Supplementary Program

► Tube pocket

Tube pocket

Item no.	8341CE	Description	Dimensions	Quantity
eClass-Nr.		Tube pocket	90x92 cm	1
Material	PP/PE			
PU	32			
SU	0			




Drapes & Packs

Supplementary Program

► Under-Buttocks Drape


Under buttocks drape

Item no.	8483CEA	Description	Dimensions	Quantity
eClass-Nr.		Under buttocks drape with cuff 24cm, reinforced	96x126 cm	1
Material	PE			
PU	56			
SU	0			




Under buttocks drape

Item no.	DYNJPE6001	Description	Dimensions	Quantity
eClass-Nr.		Wrap	60x60 cm	1
Material	PE	Blue PE under buttocks drape, wetlaid reinf.	89x119 cm	1
PU	40			
SU	0			



Under buttocks drape

Item no.	DYNJPE6002	Description	Dimensions	Quantity
eClass-Nr.		Wrap	61x61 cm	1
Material	PE	Under buttocks drape, pouch 74x81 cm w/ filter and exit port	85x36 cm	1
PU	20			
SU	0			




Drapes & Packs

Supplementary Program

► Under-Buttocks Drape


Under buttocks drape

Item no.	DYJPEFP016	Description	Dimensions	Quantity
eClass-Nr.		Wrap	61x61 cm	1
Material	PE	Underbuttocks drape, pouch w/ filter and exit port, adh.	89x117 cm	1
PU	30			
SU	0			




Under buttocks drape with pouch

Item no.	DYNJPE6004	Description	Dimensions	Quantity
eClass-Nr.		Wrap	61x61 cm	1
Material	PE	Under buttocks drape, 1 fluid collection pouch	102x112 cm	1
PU	20			
SU	0			



Under Buttocks Drape with Pouch

Item no.	8484CEA	Description	Dimensions	Quantity
eClass-Nr.		Under bottock drape with pouch	96,5x118 cm	1
Material	PE	CSR wrap	112x112 cm	1
PU	36			
SU	0			




Drapes & Packs

Supplementary Program

► Under-Buttocks Drape

Under buttocks drape



Item no.	8482CEA
eClass-Nr.	34-35-01-01
Material	PE
Dimensions	96,5 x 92cm
PU	36
SU	0


Drapes & Packs

Urological Surgery

► T.U.R. Drapes and Packs


T.U.R. pack

Item no.	TB29081CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-08	Absorbent Towel	39x58 cm	2
Material	Ultimate - 3 la	Adhesive tape	9x50 cm	1
PU	7	Drape TUR, w/ adh Pouch	150/251x174 cm	1
SU	0	Back table cover	140x190 cm	1




T.U.R. pack,

Item no.	TB29361CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-08	drape attached leggings, pouch with mesh, fen 8x8cm adh, finger cot	150/260x300 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	10			
SU	0			



TUR drape

Item no.	DYJPETUDEC4	Description	Dimensions	Quantity
eClass-Nr.		Wrap	91x91 cm	1
Material	Eclipse - Spunl	TUR drape, 6 cm diam. adh. suprapubic fen., 6 cm diam. lower fen., pouch	256/160x214 cm	1
PU	8			
SU	0			




Drapes & Packs

Urological Surgery

► T.U.R. Drapes and Packs


TUR pack

Item no.	ES10381CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-08	Absorbent Hand towels	20x30 cm	2
Material	ESS - 2 layers	CSR wrap	60x60 cm	1
PU	7	Adhesive tapes	9x50 cm	2
SU	0	TUR drape with attached leggings	250/280x280 cm	1
		Back table cover	140x190 cm	1




TUR pack

Item no.	DYJPETUPSM3	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Hand towel</td> <td>34x36 cm</td> <td>2</td>		Hand towel	34x36 cm	2
Material	ADV - SMMMS	Back table cover,	152x190 cm	1
PU	10	TUR drape, fen. 3.5cm and 6cm., w/ pouch, finger cot	160/236x191 cm	1
SU	0			



TUR pack

Item no.	DYJPETUPSM	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Hand towel</td> <td>34x36 cm</td> <td>1</td>		Hand towel	34x36 cm	1
Material	ADV - SMMMS	Back table cover,	152x190 cm	1
PU	7	TUR drape, 8x8 cm adh. fen., 10x7 cm upper adh. fen., 1 pouch	79/145x294 cm	1
SU	0			




Drapes & Packs

Urological Surgery

► T.U.R. Drapes and Packs


TUR Pack

Item no.	ES29081CE	Description	Dimensions	Quantity
eClass-Nr.		Absorbent Hand towels	20x30 cm	2
Material	ESS - 2 layers	Adhesive tape	9x50 cm	1
PU	7	Back table cover	140x190 cm	1
SU	555	SMS wrap		1
		TUR drape with attached leggings	175/242x290 cm	1



TUR pack

Item no.	TB27002CE	Description	Dimensions	Quantity
eClass-Nr.		Absorbent Hand towels	20x30 cm	1
Material	Ultimate - 3 la	Back table cover	140x190 cm	1
PU	6	TUR drape with pouch/tubing	110/240x210 cm	1
SU	0			




Drapes & Packs

Urological Surgery

► Urology packs


Cystoscopy Pack

Item no.	TB29380CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-08	Absorbent Towel	39x58 cm	1
Material	Ultimate - 3 la	Leggings	140x79cm	2
PU	10	Cysto-T-sheet	98/163x232 cm	1
SU	0	Back table cover	140x190 cm	1




Cystoscopy Pack

Item no.	TB29471CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-08	Absorbent Towel	39x58 cm	1
Material	Ultimate - 3 la	OR gown, Ultimate, XLL, impervious, breathable, raglan sleeves	150 cm	1
PU	6	Leggings	140x79cm	2
SU	0	Suture bag		1
		Cysto-T-sheet	98/163x232 cm	1
		Back table cover	140x190 cm	0



Cystoscopy pack

Item no.	DYJPE5010SM	Description	Dimensions	Quantity
eClass-Nr.	ADV - SMMMS	Cystoscopy drape, fen. 8x18 cm, mesh screen 33x91 cm	71/267x203 cm	1
Material	ADV - SMMMS	Hand towel	51x43 cm	1
PU	10	Back table cover	110x190 cm	1
SU	0			




Drapes & Packs

Urological Surgery

► Urology packs

Cystoscopy pack

Item no.	ES15311CE	Description	Dimensions	Quantity
eClass-Nr.		Drape, 6x12cm diameter ovale fenestration and 2x16cm cuts	75 x 120 cm	1
Material	ESS - 2 layers			
PU	45			
SU	0			




Cystoscopy Pack

Item no.	TB277044CE	Description	Dimensions	Quantity
eClass-Nr.		Cystoscopy drape, 6 cm diam. suprapubic adh. fen., 6 cm diam. low..	150/310x250 cm	1
Material	Ultimate - 3 la	Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
PU	8	Back table cover	140x190 cm	1
SU	0			



Cystoscopy pack I

Item no.	DYNJPE5000	Description	Dimensions	Quantity
eClass-Nr.		Leggings w/ cuff 15 cm	74x109 cm	2
Material	ECL - Spunlac	Cystoscopy drape, fen. 8x18 cm w/ reinf. and mesh screen 25x25 cm	89/152X221 cm	1
PU	14	Back table cover	110x190 cm	1
SU	0	Hand towel	51x43 cm	1




Drapes & Packs

Urological Surgery

► Urology packs


Cystoscopy pack,

Item no.	TB29360CE	Description	Dimensions	Quantity
eClass-Nr.	34-35-01-08	Hand towels	30x20 cm	1
Material	Ultimate - 3 la	drape attached leggings, reinf, fenestration 8x15cm, mesh screen 3..	89/262x220 cm	1
PU	11	Back table cover	140x190 cm	1
SU	0			



Nephroscopy drape

Item no.	DYNJPE1221	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Wrap</td> <td>91x91 cm</td> <td>1</td>		Wrap	91x91 cm	1
Material	Eclipse - Spunl	Nephroscopy drape, 4 flaps 11x51 cm, oval incise film 30x20 cm w..	300x183 cm	1
PU	8			
SU	0			



Phimosis drape

Item no.	CPU30121CE	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Phimosis drape, 4 cm diam. adh. fen.</td> <td>145x200 cm</td> <td>1</td>		Phimosis drape, 4 cm diam. adh. fen.	145x200 cm	1
Material	Ultimate - 3 la	CSR wrap	86x86 cm	1
PU	28			
SU	0			




Drapes & Packs

Urological Surgery

► Urology packs

Rectum pack

Item no.	DYJPEREPEC	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, SMS reinf. 55x88 cm	80x142 cm	1
Material	Eclipse - Spunl	Rectum drape, adh. abdominal fen. 29x29 cm, adh. rectum fen. 16x7 cm	258/250x343 cm	1
PU	6	Back table cover,	152x190 cm	1
SU	0	Hand towel	34x36 cm	4




Equipment Covers

Back Table Cover


Back table cover

Item no.	DYJPETC060	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-01	Wrap	61x61 cm	1
Material	PE	Back table cover,	152x190 cm	1
PU	22			
SU	0			




Back table cover

Item no.	DYJPETC065	Description	Dimensions	Quantity
eClass-Nr.		Wrap	61x61 cm	1
Material	PE	Back table cover, reinf. 76x152 cm	152x152 cm	1
PU	22			
SU	0			



Back table cover

Item no.	DYNJPE2309	Description	Dimensions	Quantity
eClass-Nr.		Wrap	61x61 cm	1
Material	PE	Back table cover, reinf. 60x193 cm	112x193 cm	1
PU	22			
SU	0			




Equipment Covers

Back Table Cover


Back table cover

Item no.	DYJPETC001	Description	Dimensions	Quantity
eClass-Nr.		Wrap	60x60 cm	1
Material	PE	Back table cover, reinf. 61x112 cm	112x112 cm	1
PU	20			
SU	0			




Back table cover

Item no.	DYNJPE2320	Description	Dimensions	Quantity
eClass-Nr.		Wrap	61x61 cm	1
Material	PE	Back table cover, reinf. 76x229 cm	152x229 cm	1
PU	22			
SU	0			



Back table cover

Item no.	DYNJPE2312	Description	Dimensions	Quantity
eClass-Nr.		Wrap	61x61 cm	1
Material	PE	Back table cover, reinf. 102x229 cm	196x229 cm	1
PU	22			
SU	0			




Equipment Covers

Back Table Cover


Back table cover

Item no.	DYNJPE2314	Description	Dimensions	Quantity
eClass-Nr.		Wrap	66x66 cm	1
Material	PE	Back table cover, reinf. 76x229 cm	152x229 cm	1
PU	12			
SU	0			



Back table cover

Item no.	DYNJPE2317	Description	Dimensions	Quantity
eClass-Nr.		Wrap	80x80 cm	1
Material	PE	Back table cover, reinf. 102x279 cm	201x279 cm	1
PU	12			
SU	0			



Back table cover, reinforced, 50μ, 125x230cm

Item no.	8378CEB	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-01	Back table cover, reinforced, 50μ	125x230 cm	1
Material	PE	SMS wrap	60x60 cm	1
PU	18			
SU	0			



Equipment Covers

Back Table Cover

Back table cover, reinforced, 50 μ , 140x120cm

Item no.	8368CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-01	Back table cover, folded, 50 μ	140x120 cm	1
Material	PE	SMS wrap	60x60 cm	1
PU	30			
SU	30			



Back table cover, reinforced, 50 μ , 140x190cm

Item no.	8371CEB	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-01	Back table cover, folded, 50 μ	140x190 cm	1
Material	PE	SMS wrap	60x60 cm	1
PU	22			
SU	0			



Back table cover, reinforced, 75 μ

Item no.	1900CEB	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-01	Back table cover, 75 μ , folded	140x190cm	1
Material	PE	SMS wrap	60x60 cm	1
PU	22			
SU	0			




Equipment Covers

Back Table Cover

Back table cover, reinforced, 75µ

Item no.	2100CEB	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-01	Back table cover, folded, 75µm	140x120 cm	1
Material	PE	SMS wrap	60x60 cm	1
PU	26			
SU	0			




Back table cover, reinforced, 75µ

Item no.	8380CEB	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-01	Back table cover, folded, 75µm	203x240 cm	1
Material	PE	SMS wrap		1
PU	14			
SU	0			



Back table cover, reinforced, 90µ

Item no.	8186CEB	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-01	Back table cover, folded, 90 µm	160x240 cm	1
Material	PE	SMS wrap	60x60 cm	1
PU	14			
SU	0			




Equipment Covers

Banded Bags


Banded bag, blue

Item no.	ICE60028S	Description	Dimensions	Quantity
eClass-Nr.		Banded bag with sewn elastic, blue	A=Ø45 cm, B=35	1
Material	PE			
PU	25			
SU	0			




Banded bag, bordeaux

Item no.	ICE60030S	Description	Dimensions	Quantity
eClass-Nr.		Banded bag with sewn elastic, bordeaux	A=Ø56 cm, B=40	1
Material	PE			
PU	25			
SU	0			



Banded bag, brown (rect)

Item no.	ICEBBG095	Description	Dimensions	Quantity
eClass-Nr.		Banded bag with sewn elastic, brown	95x85 cm	1
Material	PE			
PU	50			
SU	0			




Equipment Covers

Banded Bags


Banded bag, green (rect)

Item no.	ICE63654S	Description	Dimensions	Quantity
eClass-Nr.		Banded bag with sewn elastic, green	95x130 cm	1
Material	PE			
PU	25			
SU	0			




Banded bag, pink

Item no.	ICE60040S	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Banded bag with sewn elastic, pink</td> <td>A=Ø58 cm, B=50</td> <td>1</td>		Banded bag with sewn elastic, pink	A=Ø58 cm, B=50	1
Material	PE			
PU	25			
SU	0			



Banded bag, red

Item no.	ICE60020S	Description	Dimensions	Quantity
eClass-Nr. <td></td> <td>Banded bag with sewn elastic, red</td> <td>A=Ø30 cm, B=25</td> <td>1</td>		Banded bag with sewn elastic, red	A=Ø30 cm, B=25	1
Material	PE			
PU	25			
SU	0			




Equipment Covers

Banded Bags


Banded bag, white

Item no.	ICE60050S	Description	Dimensions	Quantity
eClass-Nr.		Banded bag with sewn elastic, white	A=Ø72 cm, B=60	1
Material	PE			
PU	25			
SU	0			




Banded bag, white (rect)

Item no.	ICE64848S	Description	Dimensions	Quantity
eClass-Nr.		Banded bag with sewn elastic, white	122x122 cm	1
Material	PE			
PU	25			
SU	0			



Banded bag, white (rect)

Item no.	ICE63636R	Description	Dimensions	Quantity
eClass-Nr.		Banded bag with sewn elastic, white	91x91 cm	1
Material	PE			
PU	25			
SU	0			




Equipment Covers

Banded Bags


Banded bag, white (rect)

Item no.	ICE63628R	Description	Dimensions	Quantity
eClass-Nr.		Banded bag with sewn elastic, white	91x71 cm	1
Material	PE			
PU	25			
SU	0			



Banded bag, white (rect)

Item no.	ICE63036R	Description	Dimensions	Quantity
eClass-Nr.		Banded bag with sewn elastic, white	76x91 cm	1
Material	PE			
PU	25			
SU	0			



Banded bag, white (rect)

Item no.	ICE65656S	Description	Dimensions	Quantity
eClass-Nr.		Banded bag with sewn elastic, white	142x142 cm	1
Material	PE			
PU	25			
SU	555			




Equipment Covers

Banded Bags

Banded bag, yellow

Item no.	ICEBBG070	Description	Dimensions	Quantity
eClass-Nr.		Banded bag with sewn elastic, yellow	A=Ø75 cm, B=65	1
Material	PE			
PU	50			
SU	0			



Equipment Covers

C-Arm Cover

C-arm cover set 2000 Siremobil

Item no.	DYNJE4547	Description	Dimensions	Quantity
eClass-Nr.		Banded bag, 70x130 cm, white	70x130 cm	1
Material	PE	Banded bag, 60x115 cm, green	60x115 cm	1
PU	20	C-Arm cover, 45x183 cm	45x183 cm	1
SU	0			



C-ARM COVER SET PHILIPS, 2 PARTS

Item no.	ICE8232	Description	Dimensions	Quantity
eClass-Nr.		Arm cover with 6 clips	8.5/10 x 205 cm	1
Material	PE	Banded bag, rectangular, sewn elastic band	75x75 cm	1
PU	20			
SU	0			



C-ARM COVER SET PHILIPS, 3 PARTS

Item no.	ICE8230	Description	Dimensions	Quantity
eClass-Nr.		Arm cover with 6 clips	8.5/10 x 205 cm	1
Material	PE	Banded bag, rectangular, sewn elastic band	75x75 cm	2
PU	20			
SU	0			



Equipment Covers

C-Arm Cover

C-Arm Cover Set Philips, BV 300

Item no.	DYNJE4527	Description	Dimensions	Quantity
eClass-Nr.		Banded bag, rectangular, 95x85 cm, brown	95x85 cm	1
Material	PE	Banded bag, 60x115 cm, green	60x115 cm	1
PU	20	C-Arm cover, 30x198 cm	30x198 cm	1
SU	0			



C-arm cover set Siemens Siremobil Iso-C, Iso-C3D/ Arcadis Orbic, Orbic 3D

Item no.	DYNJE4557	Description	Dimensions	Quantity
eClass-Nr.		Clear poly drape with clips, Top	51/82x130 cm	1
Material	PE			
PU	20			
SU	0			



C-ARM COVER SET SIEMENS, 2 PARTS

Item no.	ICE8132	Description	Dimensions	Quantity
eClass-Nr.		Arm cover with 6 clips	10/14 x 225 cm	1
Material	PE	Banded bag, rectangular, sewn elastic band	75x90 cm	1
PU	20			
SU	0			



Equipment Covers

C-Arm Cover

C-ARM COVER SET SIEMENS, 3 PARTS

Item no.	ICE8130	Description	Dimensions	Quantity
eClass-Nr.		Arm cover with 6 clips	10/14 x 225 cm	1
Material	PE	Banded bag, rectangular, sewn elastic band	75x90 cm	2
PU	20			
SU	0			



C-ARM COVER SET ZIEHM, 2 PARTS

Item no.	ICE8332	Description	Dimensions	Quantity
eClass-Nr.		Arm cover with 6 clips	8.5/10 x 180 cm	1
Material	PE	Banded bag, rectangular, sewn elastic band	75x75 cm	1
PU	20			
SU	0			



C-ARM COVER SET ZIEHM, 3 PARTS

Item no.	ICE8330	Description	Dimensions	Quantity
eClass-Nr.		Arm cover with 6 clips	8.5/10 x 180 cm	1
Material	PE	Banded bag, rectangular, sewn elastic band	75x75 cm	2
PU	20			
SU	0			



Equipment Covers

C-Arm Cover


C-arm cover set, Philips BV 25/26

Item no.	DYNJE4530	Description	Dimensions	Quantity
eClass-Nr.		Banded bag, 60x115 cm, green	60x115 cm	2
Material	PE	C-Arm cover, 30x198 cm	30x198 cm	1
PU	20			
SU	0			




C-arm cover set, Philips BV 29/300

Item no.	DYNJE4526	Description	Dimensions	Quantity
eClass-Nr.		Banded bag, 70x130 cm, white	70x130 cm	1
Material	PE	Banded bag, 60x115 cm, green	60x115 cm	1
PU	20	C-Arm cover, 30x198 cm	30x198 cm	1
SU	0			



C-arm cover set, Philips Pulsera/Endura 30, 5cm

Item no.	DYNJE4569	Description	Dimensions	Quantity
eClass-Nr.		Banded bag, 60x115 cm, yellow	60x115 cm	1
Material	PE	Banded bag, 70x140 cm, white	70x140 cm	1
PU	20	C-Arm cover, 30x198 cm	30x198 cm	1
SU	0			



Equipment Covers

C-Arm Cover


C-arm cover set, Siemens Iso-C, Iso-C3D/ Arcadis Orbic, Orbic 3D

Item no.	DYNJE4556	Description	Dimensions	Quantity
eClass-Nr.		Clear poly drape with clips, bottom	53/46x134 cm	1
Material	PE	Clear poly drape with clips, Top	51/82x130 cm	1
PU	20			
SU	0			




C-arm cover set, Siemens Siremobil

Item no.	DYNJE4532	Description	Dimensions	Quantity
eClass-Nr.		Banded bag, 60x115 cm, white	60x115 cm	1
Material	PE	Banded bag, 60x115 cm, green	60x115 cm	1
PU	20	C-Arm cover, 35x183 cm	35x183 cm	1
SU	0			



C-arm cover set, Ziehm

Item no.	DYNJE4531	Description	Dimensions	Quantity
eClass-Nr.		Banded bag, 70x130 cm, white	70x130 cm	1
Material	PE	Banded bag, 60x115 cm, green	60x115 cm	1
PU	20	C-Arm cover, 35x185 cm	35x185 cm	1
SU	0			




Equipment Covers

C-Arm Cover


C-Arm drape

Item no.	ICE4400	Description	Dimensions	Quantity
eClass-Nr.		C-arm drape, w/ adh. straps	104x188 cm	1
Material	Clear PE film			
PU	20			
SU	0			




C-Arm drape

Item no.	ICE4430	Description	Dimensions	Quantity
eClass-Nr.		C-arm drape, w/ adh. straps	152x267 cm	1
Material	Clear PE film			
PU	20			
SU	0			



C-Arm Drape with Rubberbands

Item no.	ICE4400R	Description	Dimensions	Quantity
eClass-Nr.		C-arm drape, w/ 3 adh. straps, 2 elastic bands	104x188 cm	1
Material	Clear PE film			
PU	20			
SU	0			




Equipment Covers

C-Arm Cover

Full Size C-Arm Kit

Item no.	ICE4500	Description	Dimensions	Quantity
eClass-Nr.		Intensifier drape	91x91 cm	1
Material	PE	Footswitch drape	46x41 cm	1
PU	20	C-arm drape with plastic clamps	10/13x218 cm	1
SU	0	Banded bag	71x31 cm	1




Mini C-arm drape

Item no.	ICE4410	Description	Dimensions	Quantity
eClass-Nr.		Mini C-arm drape, elastic band on the back, 3 adh. straps	137x160 cm	1
Material	Clear PE film			
PU	20			
SU	0			



Mini C-arm drape


Item no.	ICE4420	Description	Dimensions	Quantity
eClass-Nr.				
Material	Clear PE film			
PU	20			
SU	0			



Equipment Covers


Grip Cover

Light handle cover L



Item no.	16558
eClass-Nr.	34-35-04-90
Material	PE
Dimensions	L
PU	200
SU	200

Light handle cover S



Item no.	16557
eClass-Nr.	34-35-04-90
Material	PE
Dimensions	S
PU	200
SU	200

Light Handle Cover Sterile (1Ea), only sold in Germany



Item no.	26-011
eClass-Nr.	34-35-04-90
Material	PE
Dimensions	n/a
PU	100
SU	100

Light Handle Cover Sterile (2Ea), only sold in Germany



Item no.	26-012
eClass-Nr.	34-35-04-90
Material	PE
Dimensions	n/a
PU	60
SU	60

Equipment Covers

Mayo Stand Cover


Heavy duty Mayo Stand Cover

Item no.	8338CE	Description	Dimensions	Quantity
eClass-Nr.		Mayo Stand Cover, absorbent reinf. 55x88 cm	80x142 cm	1
Material	PP/PE	CSR wrap	60x60 cm	1
PU	44			
SU	0			




Mayo stand cover

Item no.	DYJPEMAYOM	Description	Dimensions	Quantity
eClass-Nr.		Wrap	76x76 cm	1
Material	PE	Mayo stand cover, absorbent reinf. 61x85 cm	80x145 cm	1
PU	30			
SU	0			



Mayo Stand Cover

Item no.	DYNJPE2510	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
Material	PE	SMS wrap		1
PU	22			
SU	0			



Equipment Covers

Mayo Stand Cover

Mayo Stand Cover

Item no.	8339DCE	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, absorbent reinf. 55x88 cm	80x142 cm	1
Material	PE	CSR wrap	60x60 cm	1
PU	48			
SU	0			



Mayo stand cover, Large

Item no.	8339CEC	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-01	Mayo stand cover, reinforced	80x142 cm	1
Material	PE	CSR wrap	60x60 cm	1
PU	50			
SU	50			



Mayo stand cover, reinforced, Small

Item no.	8337CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-01	CSR wrap	60x60 cm	1
Material	PE	MAYO Stand cover, S, 50μ, absorbent reinf. 46x61 cm	59x141 cm	1
PU	56			
SU	56			




Equipment Covers

Mayo Stand Cover


Orthopaedic Mayo Stand Cover

Item no.	8340DCE	Description	Dimensions	Quantity
eClass-Nr.		Mayo stand cover, ultra absorbent reinforcement 75x100 cm	80x145 cm	1
Material	PE	CSR wrap	60x60 cm	1
PU	40			
SU	0			



Orthopedic mayo stand cover, absorbent reinforcement

Item no.	8340CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-01	Mayo stand cover, ultra absorbent reinforcement 75x100 cm	80x145 cm	1
Material	PE	CSR wrap	60x60 cm	1
PU	28			
SU	0			




Equipment Covers

Microscope Cover


Microscope Cove R, Wild, Leitz, 68 mm

Item no.	16516M	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-02			
Material	PE			
PU	10			
SU	10			



Microscope cover, Zeiss opmi

Item no.	16500M	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-02	Microscope cover, Zeiss Opmi	104 x 163cm, 48m	1
Material	PE			
PU	10			
SU	10			



Microscope cover, Zeiss opmi md

Item no.	16502M	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-02	Microscope Cover, Zeiss, Opmi MD	117 x 163cm, 65m	1
Material	PE			
PU	10			
SU	10			




Equipment Covers

Microscope Cover

Microscope cover, Zeiss opmi md

Item no.	16508M	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-02			
Material	PE			
PU	10			
SU	10			



Microscope drape, 48mm diam lens

Item no.	DYNJE5010	Description	Dimensions	Quantity
eClass-Nr.		Microscope drape, 48 mm diam. lens, 3 adh. strips, 1 pair of oculars	51x163 cm	1
Material	Clear PE film			
PU	20			
SU	0			



Microscope drape, 48mm diam lens

Item no.	DYNJE5040	Description	Dimensions	Quantity
eClass-Nr.		Microscope drape, 48 mm diam. lens, 3 adh. strips, 1 pair of oculars	51x267 cm	1
Material	Clear PE film			
PU	40			
SU	0			




Equipment Covers

Microscope Cover

Microscope drape, 65mm diam lens

Item no.	DYNJE5210	Description	Dimensions	Quantity
eClass-Nr.		Microscope drape, 65 m diam. lens, 3 adh. strips, 3 pairs of oculars	117x267 cm	1
Material	Clear PE film			
PU	10			
SU	0			




Microscope drape, 65mm diam lens

Item no.	DYNJE5220	Description	Dimensions	Quantity
eClass-Nr.		Microscope drape, 65 mm diam. lens, 5 adh. strips, 3 pairs of oculars	137x381 cm	1
Material	Clear PE film			
PU	10			
SU	0			



Microscope drape, 65mm diam lens

Item no.	DYNJE5030	Description	Dimensions	Quantity
eClass-Nr.		Microscope drape, 68 mm diam. lens, 4 adh. strips, 3 pairs of oculars	117x305 cm	1
Material	Clear PE film			
PU	10			
SU	0			




Equipment Covers

Microscope Cover


Microscope drape, 65mm diam lens

Item no.	DYNJE5200	Description	Dimensions	Quantity
eClass-Nr.		Microscope drape, 65 mm diam. lens, 3 adh. strips, 3 pairs of oculars	117x203 cm	1
Material	Clear PE film			
PU	10			
SU	0			



Microscope drape, 68mm diam lens

Item no.	DYNJE5300	Description	Dimensions	Quantity
eClass-Nr.				
Material	Clear PE film			
PU	10			
SU	0			



Microscope drape, 68mm diam lens

Item no.	DYNJE5320	Description	Dimensions	Quantity
eClass-Nr.		Microscope drape, 68 mm diam. lens, 3 adh. strips, 3 pairs of oculars	137x267 cm	1
Material	Clear PE film			
PU	10			
SU	0			



Equipment Covers

Microscope Cover

Microscope drape, 68mm diam lens

Item no.	DYNJE5330	Description	Dimensions	Quantity
eClass-Nr.		Microscope drape, 68 mm diam. lens, 5 adh. strips, 3 pairs of oculars	137x381 cm	1
Material	Clear PE film			
PU	10			
SU	0			




Equipment Covers

Other Covers


Armboard covers

Item no.	9460CE	Description	Dimensions	Quantity
eClass-Nr.		Armboard cover	76x36 cm	2
Material	ECL - Spunlac			
PU	34			
SU	0			




Cassette cover

Item no.	ICE4000	Description	Dimensions	Quantity
eClass-Nr.		X-ray cassette cover	61x102 cm	1
Material	Clear PE film			
PU	50			
SU	0			



Cassette cover

Item no.	ICE4050	Description	Dimensions	Quantity
eClass-Nr.		X-ray cassette cover	61x62 cm	1
Material	Clear PE film			
PU	50			
SU	0			




Equipment Covers

Other Covers


Cassette cover with liner

Item no.	ICE4010	Description	Dimensions	Quantity
eClass-Nr.		X-ray cassette cover	61x64 cm	1
Material	Clear PE film			
PU	48			
SU	0			




Cover chair

Item no.	DYJPECHCEC	Description	Dimensions	Quantity
eClass-Nr. <td><td>Wrap</td><td>45x45 cm</td><td>1</td></td>	<td>Wrap</td> <td>45x45 cm</td> <td>1</td>	Wrap	45x45 cm	1
Material	Eclipse - Spunl	Cover chair	14x38 cm	1
PU	30			
SU	0			



Cover chair

Item no.	ES304070CE	Description	Dimensions	Quantity
eClass-Nr. <td><td>Cover chair</td><td>75x75 cm</td><td>1</td></td>	<td>Cover chair</td> <td>75x75 cm</td> <td>1</td>	Cover chair	75x75 cm	1
Material	ESS - 2 layers			
PU	80			
SU	0			




Equipment Covers

Other Covers


Cover equipment XR cassette

Item no.	8720CEA	Description	Dimensions	Quantity
eClass-Nr.	34-35-04-08			
Material	PE			
PU	30			
SU	0			



Fluoroscope cover

Item no.	DYNJSD1012	Description	Dimensions	Quantity
eClass-Nr.				
Material	PE			
PU	10			
SU	0			



Remote control cover

Item no.	DYNJE4635	Description	Dimensions	Quantity
eClass-Nr.		Remote control cover	10x28 cm	1
Material	PE			
PU	200			
SU	0			



Equipment Covers

Other Covers


Remote Control Cover

Item no.	DYNJE4641	Description	Dimensions	Quantity
eClass-Nr.		Remote control cover	18x35 cm	1
Material	PE			
PU	200			
SU	0			



Setup cover

Item no.	ICE2008	Description	Dimensions	Quantity
eClass-Nr.		Setup cover	137x229 cm	1
Material	PE			
PU	20			
SU	0			



Setup cover

Item no.	ICE35975	Description	Dimensions	Quantity
eClass-Nr.		Clear setup cover, 150x190 cm, 3 adhesive tapes	150x190 cm	1
Material	PE	SMS wrap		1
PU	40			
SU	555			




Equipment Covers

Other Covers

Thermo cover

Item no.	DYNJE7000	Description	Dimensions	Quantity
eClass-Nr.		Thermo cover	68 cm	1
Material	Clear PE film			
PU	40			
SU	0			



Equipment Covers

Other Covers

Birthing ball cover



Item no.	SPTe-ESBB01
eClass-Nr.	
Material	ESS - 2 layers
Dimensions	
PU	96
SU	0

Equipment Covers

Probe Covers

Probe cover kit

Item no.	ICE5900	Description	Dimensions	Quantity
eClass-Nr.		Probe cover	10x122 cm	1
Material	clear Polyuret	rubberbands		2
PU	20	SMS wrap		1
SU	0	Ultrasound gel		1



Probe cover kit

Item no.	ICE5910	Description	Dimensions	Quantity
eClass-Nr.		Probe cover	61x15 cm	1
Material	clear Polyuret	rubberbands		2
PU	20	SMS wrap		1
SU	0	Ultrasound gel		1



Probe cover kit

Item no.	ICE5920	Description	Dimensions	Quantity
eClass-Nr.		Probe cover	122x15 cm	1
Material	clear Polyuret	rubberbands		1
PU	20	SMS wrap		1
SU	0	Ultrasound gel		1



Equipment Covers

Probe Covers

Probe cover kit

Item no.	ICE5930	Description	Dimensions	Quantity
eClass-Nr.		Probe cover	244x15 cm	1
Material	clear Polyuret	rubberbands		1
PU	20	SMS wrap		1
SU	0	Ultrasound gel		1



Equipment Covers

Video Camera Cover

Cable cover, 8x200 cm, cardboard tab for application and extra adhes..

Item no.	ICE4100	Description	Dimensions	Quantity
eClass-Nr.		Cable cover, with cardboard tab for application and extra adhesive ..	8x200 cm	1
Material	Clear PE film			
PU	25			
SU	0			



Video camera cover with pre-cut tip, 23x244 cm, with stretchable adh..

Item no.	ICE4300	Description	Dimensions	Quantity
eClass-Nr.		Camera / arm drape, 2 adh. strips	23x244 cm	1
Material	Clear PE film			
PU	25			
SU	0			



Video camera cover, 15x298 cm, with cardboard tab for application, w..

Item no.	ICE4270	Description	Dimensions	Quantity
eClass-Nr.		Video Camera cover, concertina folding	15x298 cm	1
Material	Clear PE film			
PU	25			
SU	0			



Equipment Covers

Video Camera Cover

Video camera cover, 18x244 cm, with tip elastoflex, stretchable adhe..

Item no.	ICE4200	Description	Dimensions	Quantity
eClass-Nr.		Camera / Laser arm drape, perforated elastomeric tip, adh. straps	18x244 cm	1
Material	Clear PE film			
PU	25			
SU	0			



Video camera cover, with pre-cut tip, 13x240 cm, with cardboard tab ..

Item no.	ICE4240	Description	Dimensions	Quantity
eClass-Nr.		Video camera cover with pre-cut tip, with cardboard tab for applicat..	13x240 cm	1
Material	Clear PE film			
PU	25			
SU	0			



Video camera cover, with pre-cut tip, 13x244 cm, 2 extra adhesive ta..

Item no.	ICE4250	Description	Dimensions	Quantity
eClass-Nr.		Camera / arm drape, 2 adh. strips	13x244 cm	1
Material	Clear PE film			
PU	25			
SU	0			



Equipment Covers

Video Camera Cover

Video camera cover, with pre-cut tip, 15x240 cm, 1 extra adhesive ta..

Item no.	ICE4260	Description	Dimensions	Quantity
eClass-Nr.		Video camera cover with ring	15x240 cm	1
Material	Clear PE film			
PU	25			
SU	0			



Invisishield

Isolation bag

Item no.	ICE1003	Description	Dimensions	Quantity
eClass-Nr.		Isolation bag with closure string, 51x51 cm	51x51 cm	1
Material	Clear PE film			
PU	40			
SU	0			



Invisishield

Clear plastic drapes

Aperture pouch

Item no.	IDE1072	Description	Dimensions	Quantity
eClass-Nr.		Aperture pouch	47x23 cm	1
Material	Matte PE film			
PU	40			
SU	0			




Clear drape with adhesive

Item no.	5810CE	Description	Dimensions	Quantity
eClass-Nr.		Clear drape with adhesive	150x300 cm	1
Material	Clear PE film	SMS wrap		1
PU	17			
SU	0			



Drape with adhesive fenestration 5x7 cm

Item no.	IDE1030	Description	Dimensions	Quantity
eClass-Nr.		Drape with adhesive fenestration 5x7 cm	82.5x129.5 cm	1
Material	Matte PE film			
PU	40			
SU	0			




Invisishield

Clear plastic drapes

Drape with adhesive fenestration diam.6 cm

Item no.	IDE1020	Description	Dimensions	Quantity
eClass-Nr.		Drape with adhesive fenestration diam.6 cm	41x41 cm	1
Material	Matte PE film			
PU	40			
SU	0			




Instrument drape

Item no.	ICE2158	Description	Dimensions	Quantity
eClass-Nr.		Instrument drape	61x46 cm	1
Material	Matte PE film			
PU	40			
SU	0			



Ophthalmic aperture pouch

Item no.	IDE1024	Description	Dimensions	Quantity
eClass-Nr.		Ophthalmic aperture pouch	48x41 cm	1
Material	Matte PE film			
PU	40			
SU	0			




Invisishield

Clear plastic drapes

Ophthalmic incise drape with fluid pouch

Item no.	IDE1061	Description	Dimensions	Quantity
eClass-Nr.		Ophthalmic incise drape with fluid pouch	122x129 cm	1
Material	Matte PE film			
PU	40			
SU	0			



Ophthalmic incise drape with pouch

Item no.	DYNJSD1023	Description	Dimensions	Quantity
eClass-Nr.		Ophthalmic drape with incise film, 12x10 cm, 2 adhesives tabs, 5x..	48x41 cm	1
Material	Matte PE film			
PU	40			
SU	0			



Patient isolation drape

Item no.	DYNJSD1014	Description	Dimensions	Quantity
eClass-Nr.		Patient isolation drape with incise film, 49x25 cm	317x213 cm	1
Material	Clear PE film			
PU	15			
SU	0			



Invisishield

Clear plastic drapes

Patient isolation drape with pouch

Item no.	DYNJSD1017	Description	Dimensions	Quantity
eClass-Nr.		Patient isolation drape with pouch, incise film, 49x25 cm, adhesive..	317x213 cm	1
Material	Clear PE film			
PU	15			
SU	0			




Patient isolation drape with pouch

Item no.	DYNJSD1017W	Description	Dimensions	Quantity
eClass-Nr.		Patient isolation drape with pouch, incise film, 28x85 cm, adhesive..	317x254 cm	1
Material	Clear PE film			
PU	15			
SU	0			



PE drape with 6x10cm adhesive fenestration

Item no.	8465CE	Description	Dimensions	Quantity
eClass-Nr.		SMS wrap		1
Material	Matte PE film	PE fenestrated drape with adhesive		1
PU	30			
SU	0			



Invisishield

Clear plastic drapes


PE drape with a 9cm diam. adhesive fenestration

Item no.	8461CE	Description	Dimensions	Quantity
eClass-Nr.		PE fenestrated drape with adhesive	60x100cm	1
Material	Matte PE film			
PU	150			
SU	0			



Poly U drape, 15x48 cm adhesive split

Item no.	IDE1015	Description	Dimensions	Quantity
eClass-Nr.		Poly U drape, 15x48 cm adhesive split	122x132 cm	1
Material	Matte PE film			
PU	40			
SU	0			



Poly U drape, 15x48 cm adhesive split

Item no.	IDE1019	Description	Dimensions	Quantity
eClass-Nr.		Poly U drape, 15x48 cm adhesive split	122x180 cm	1
Material	Matte PE film			
PU	40			
SU	0			




Invisishield

Clear plastic drapes

Towel adhesive drape

Item no.	IDE1000	Description	Dimensions	Quantity
eClass-Nr.		Towel drape	46x30 cm	1
Material	Matte PE film			
PU	40			
SU	0			



Towel adhesive drape

Item no.	IDE1010	Description	Dimensions	Quantity
eClass-Nr.		Towel adhesive drape	46x61 cm	1
Material	Matte PE film			
PU	40			
SU	0			



Urological drape

Item no.	DYNJSD1071	Description	Dimensions	Quantity
eClass-Nr.		Urological drape with adhesive, 60x5 cm, finger cot	60x91 cm	1
Material	Matte PE film			
PU	40			
SU	0			



Invisishield

Incise drapes

Clear Incise Drape

Item no.	33525	Description	Dimensions	Quantity
eClass-Nr.		Adhesive Clear Incise Drape	35x25cm / 30x25c	1
Material	Polyurethane/			
PU	120			
SU	0			



Clear Incise Drape

Item no.	34040	Description	Dimensions	Quantity
eClass-Nr.		Adhesive Clear Incise Drape	40x40cm / 34x40c	1
Material	Polyurethane/			
PU	90			
SU	0			



Clear Incise Drape

Item no.	35050	Description	Dimensions	Quantity
eClass-Nr.		Adhesive Clear Incise Drape	50x48cm / 44x48c	1
Material	Polyurethane/			
PU	90			
SU	0			



Invisishield

Incise drapes

Clear Incise Drape

Item no.	37045	Description	Dimensions	Quantity
eClass-Nr.		Adhesive Clear Incise Drape	70x45cm / 60x45c	1
Material	Polyurethane/			
PU	90			
SU	0			



Clear Incise Drape

Item no.	38060	Description	Dimensions	Quantity
eClass-Nr.		Adhesive Clear Incise Drape	80x60cm / 60x60c	1
Material	Polyurethane/			
PU	60			
SU	0			



Clear Incise Drape

Item no.	38080	Description	Dimensions	Quantity
eClass-Nr.		Adhesive Clear Incise Drape	80x80cm / 70x80c	1
Material	Polyurethane/			
PU	60			
SU	0			



Invisishield

Incise drapes

Clear Incise Drapes

Item no.	32020	Description	Dimensions	Quantity
eClass-Nr.		Adhesive Clear Incise Drape	20x20cm / 15x20c	1
Material	Polyurethane/			
PU	120			
SU	0			



Impregnated Incise Drape

Item no.	42020	Description	Dimensions	Quantity
eClass-Nr.		Adhesive Incise Drape impregnated	20x20cm / 15x20c	1
Material	Polyurethane/			
PU	120			
SU	0			



Impregnated Incise Drape

Item no.	44040	Description	Dimensions	Quantity
eClass-Nr.		Adhesive Incise Drape impregnated	40x40cm / 34x40c	1
Material	Polyurethane/			
PU	90			
SU	0			



Invisishield

Incise drapes

Impregnated Incise Drape

Item no.	47045	Description	Dimensions	Quantity
eClass-Nr.		Adhesive Incise Drape impregnated	70x45cm / 60x45c	1
Material	Polyurethane/			
PU	90			
SU	0			




Impregnated Incise Drape

Item no.	48060	Description	Dimensions	Quantity
eClass-Nr.		Adhesive Incise Drape impregnated	80x60cm / 60x60c	1
Material	Polyurethane/			
PU	60			
SU	0			



Impregnated Incise Drape

Item no.	48080	Description	Dimensions	Quantity
eClass-Nr.		Adhesive Incise Drape impregnated	80x80cm / 70x80c	1
Material	Polyurethane/			
PU	60			
SU	0			



Invisishield

Instrument and fluid collection pouches


Cesarean pouch

Item no.	DYNJSD2057	Description	Dimensions	Quantity
eClass-Nr.		360 degree cesarean pouch with incise film, 38x36 cm	91x76 cm	1
Material	Matte PE film			
PU	10			
SU	0			



Fluid pouch

Item no.	IDE2114	Description	Dimensions	Quantity
eClass-Nr.		Fluid pouch with adhesive, 30x2.5 cm	30x28 cm	1
Material	Clear PE film			
PU	40			
SU	0			



Instrument pouch

Item no.	DYJPEISTPOM1	Description	Dimensions	Quantity
eClass-Nr.		Instrument pouch, 2 compartments, with adhesive, 41x5 cm	33x41 cm	1
Material	Clear PE film			
PU	100			
SU	0			



Invisishield

Instrument and fluid collection pouches

Instrument pouch

Item no.	DYJPEISTPOM3	Description	Dimensions	Quantity
eClass-Nr.		Instrument pouch, 2 compartments, with adhesive, 43x5 cm	43x48 cm	1
Material	Clear PE film			
PU	60			
SU	0			




Instrument pouch

Item no.	DYJPEISTPOM2	Description	Dimensions	Quantity
eClass-Nr.		Instrument pouch with adhesive, 41x5 cm	33x41 cm	1
Material	Clear PE film			
PU	100			
SU	0			



Instrument pouch, adhesive

Item no.	04A040CE	Description	Dimensions	Quantity
eClass-Nr.		Instrument pouch in clear plastic with adhesive	33x51 cm	1
Material	PE			
PU	180			
SU	0			




Invisishield

Instrument and fluid collection pouches


Shoulder pouch

Item no.	DYNJSD1067	Description	Dimensions	Quantity
eClass-Nr.		Shoulder pouch with oval aperture, 19x25 cm	91x73 cm	1
Material	Matte PE film			
PU	24			
SU	0			



U fluid collection pouch

Item no.	DYJPEFP090	Description	Dimensions	Quantity
eClass-Nr.		U fluid collection pouch with adhesive aperture, 15x46 cm	76x112 cm	1
Material	Clear PE film			
PU	25			
SU	0			





Återförsäljare:

MEDI PLAST

Medioplast AB

Box 1004, 212 10 Malmö

T 020-78 80 35

mediplast.info@mediplast.com

www.mediplast.com

Copyright 2018
Medline Inc. or one of its branches
All rights reserved
LITECCen774000025

Your contact

Linnéa Nilsson
International product manager
Medioplast
+46 (0)73-076 88 13
linnea.nilsson@mediplast.com

www.medline.com/eu