

Extensions

- Different lengths between 30 - 300 cm
- Flowstop with hydrophobic filter which simplifies priming and reduces the risk of leakage and contamination
- PVC-free or PVC DEHP-free material
- Tube 2,9 x 4,1 mm or 1 x 2 mm (ID x OD)
- Micro boring with a low internal volume
- Available with needle-free connector and/or back-check valve/anti-siphon valve
- Self-ejecting rotating luer-lock
- Specific items with UV-light resistance in the range of 200-600 nm
- MRI compatible

Extension line for infusion

2,9 x 4,1 mm (ID x OD)

| Art. no. | Description | | Ward/Trsp.box |
|--------------|-------------|----------|---------------|
| 60SRL80162P | 30 cm | PVC-free | 10/100* |
| 60SRL80161P | 50 cm | PVC-free | 10/100* |
| 60SRL80170P | 75 cm | PVC-free | 10/100* |
| 60MRO4106310 | 100 cm | PVC-free | 10/100* |
| 60SRL80166P | 120 cm | PVC-free | 10/100* |
| 60SRL80164P | 150 cm | PVC-free | 10/100* |
| 60SRL80167P | 200 cm | PVC-free | 10/100* |



Extension line for syringe pump

1 x 2 mm (ID x OD)

| Art. no. | Description | | Ward/Trsp.box |
|----------------|---|---------------|---------------|
| 60SRL04100103 | 30 cm | PVC-free | 10/100* |
| 60SRL80070P | 75 cm | PVC-free | 10/100* |
| 60SRL04100115 | 150 cm | PVC-free | 10/100* |
| 60MP1501 | 150 cm | DEHP-free PVC | 10/100* |
| 60SRL04100120 | 200 cm | PVC-free | 10/100* |
| 60MP2001 | 200 cm | DEHP-free PVC | 10/100* |
| 60SRL04100130 | 300 cm | PVC-free | 10/100* |
| | | | |
| 60SRL04100115P | 150 cm, UV protection | PVC-free | 10/100* |
| 60SRL04100120P | 200 cm, UV protection | PVC-free | 10/100* |
| | | | |
| 60SRL04100215 | 150 cm, back-check valve | PVC-free | 10/100* |
| 60SRL04100220 | 200 cm, back-check valve | PVC-free | 10/100* |
| 60SRL04100230 | 300 cm, back-check valve | PVC-free | 10/100* |
| | | | |
| 60SRL04100315 | 150 cm, anti-siphon valve | PVC-free | 10/100* |
| 60SRL04100320 | 200 cm, anti-siphon valve | PVC-free | 10/100* |
| 60SRL04100330 | 300 cm, anti-siphon valve | PVC-free | 10/100* |
| | | | |
| 60SRL80150P | 150 cm, needle-free connector | PVC-free | 10/100* |
| 60SRL80200P | 200 cm, needle-free connector | PVC-free | 10/100* |
| | | | |
| 60SRL80150PP | 150 cm, back-check valve, needle-free connector | PVC-free | 10/100* |
| 60SRL80200PP | 200 cm, back-check valve, needle-free connector | PVC-free | 10/100* |

| | | |
|---|---|--|
| Material: Tube: TPE/ PUR/ PVC DEHP-free Luer-Lock: ABS Priming cap/Flowstop: ABS/PE Pinch-clamp: PP | Protective cap: PP Back-check valve: SAN/MABS/ Silicone Anti-siphon valve: PC/Silicone Needle-free connector: PC/ Silicone | Sterilization: EO Single use, sterile 1 pcs in peel pack * = minimum ordering quantity |
|---|---|--|