

micro G-JET[®]

Low Profile Pediatric G-J Enteral Feeding Tube



Designed for Pediatric Patients*

14F gastric segment transitions to an 8F jejunal segment

Anti-Kink Protection

Spans the entire 8F jejunal segment

Low Profile Bolster

With clearly labeled G + J ports

Exclusive AMT Balloon


The balloon you trust, adapted for small stomachs



Innovating. Educating. Changing Lives.™

The AMT micro G-JET®


THE ONLY LOW PROFILE GASTRIC-JEJUNAL FEEDING
DEVICE WITH AN 8F JEJUNAL SEGMENT



micro G-JET®

GASTRIC 14F	JEJUNAL 8F	LEGACY	ENFit®
LENGTH		ORDER NUMBERS	
0.8 cm	10 cm	MGJ-1408-10	MGJ-1408-10-I
	15 cm	MGJ-1408-15	MGJ-1408-15-I
	22 cm	MGJ-1408-22	MGJ-1408-22-I
1.0 cm	10 cm	MGJ-1410-10	MGJ-1410-10-I
	15 cm	MGJ-1410-15	MGJ-1410-15-I
	22 cm	MGJ-1410-22	MGJ-1410-22-I
1.2 cm	10 cm	MGJ-1412-10	MGJ-1412-10-I
	15 cm	MGJ-1412-15	MGJ-1412-15-I
	22 cm	MGJ-1412-22	MGJ-1412-22-I
1.5 cm	30 cm	MGJ-1412-30	MGJ-1412-30-I
	10 cm	MGJ-1415-10	MGJ-1415-10-I
	15 cm	MGJ-1415-15	MGJ-1415-15-I
1.7 cm	22 cm	MGJ-1415-22	MGJ-1415-22-I
	30 cm	MGJ-1415-30	MGJ-1415-30-I
	10 cm	MGJ-1417-10	MGJ-1417-10-I
2.0 cm	15 cm	MGJ-1417-15	MGJ-1417-15-I
	22 cm	MGJ-1417-22	MGJ-1417-22-I
	30 cm	MGJ-1417-30	MGJ-1417-30-I
2.3 cm	22 cm	MGJ-1420-22	MGJ-1420-22-I
	30 cm	MGJ-1420-30	MGJ-1420-30-I
	22 cm	MGJ-1423-22	MGJ-1423-22-I
2.5 cm	30 cm	MGJ-1423-30	MGJ-1423-30-I
	22 cm	MGJ-1425-22	MGJ-1425-22-I
	30 cm	MGJ-1425-30	MGJ-1425-30-I

The micro G-JET® comes packaged as 1/box (sterile). | NOTE: Do not cut the micro G-JET® to a custom length.
The micro G-JET® is compatible with guidewires up to and including 0.038" in diameter.



Äterförsäljare:


MEDI PLAST

Medioplast AB
Box 1004, 212 10 Malmö
T 020-78 80 35
info.se@mediplast.com
www.mediplast.com

Extension + Feed Sets Included in Kit

For a complete list of compatible extension and feed Sets visit: AppliedMedical.net/Enteral/Accessories

*AMT jejunal feed sets are compatible with Mic-Key® devices.



For a complete list of products visit us at www.AppliedMedical.net
800 869 7382 | CS@AppliedMedical.net | ICS@AppliedMedical.net



AMT micro G-JET®



Applied Medical Technology, Inc.
8006 Katherine Boulevard
Brecksville, OH 44141
www.AppliedMedical.net
An ISO 13485 Registered Company
MR Conditional
Made In USA



© 2021 Applied Medical Technology, Inc. C-4672-E

Patents: <http://www.AppliedMedical.net/Patents/> | Trademarks: <http://www.AppliedMedical.net/Trademarks/>

*The micro G-JET® is cleared for patients > 6 kg. | Mic-Key® is a registered trademark of Avent, Inc. | ENFit® is a registered trademark of GEDSA.