



Metal Suction Tips with transparent handle type Fukushima

- For rhinological surgeries and examinations
- Transparent handle with Fukushima suction control
- Single use, sterile, packed in peel-pack
- Sterilization: EO
- Medical grade stainless steel

Art.No.	Description	Ward/Trsp.box
66520150	2,0 x 150 mm	25/375
66530150	3,0 x 150 mm	25/375
66540150	4,0 x 150 mm	25/375
66540200	4,0 x 200 mm	25/375
66550150	5,0 x 150 mm	25/375
66550180	5,0 x 180 mm	25/375
66550220	5,0 x 220 mm	25/375

Ward box = minimum order qty

Manufacturer: Mediplast AB / www.mediplast.com